

IMPORTANT SAFETY INSTRUCTIONS

WARNING: Use carefully. May cause serious burns. **Do Not** use over sensitive skin areas or in the presence of poor circulation. The unattended use of the Heated Eye Pad by children or incapacitated persons may be dangerous. To reduce the risk of burns, electrical shock, fire, and accident, this product must be used in accordance with the following instructions:

1. **Read all instructions carefully**
2. **WARNING:** Device may cause skin irritation or burning sensation. Do not use on sensitive skin areas or in the presence of poor circulation. This product not intended for use by incapacitated individuals.
3. **DO NOT** use pad on an infants.
4. This pad is not to be used on or by an invalid, a paralyzed person, a sleeping or unconscious person, a person with diabetes, or a person with poor blood circulation. Do not use eye pad on areas of sensitive skin.
5. Burns may occur. Check skin under pad frequently to avoid burning and blistering.
6. **Place pad on top of closed eyelids. Never on eyes, or eyeball.**
7. **DO NOT** use pins or other metallic means to fasten heater eye pad in place.
8. **DO NOT** fold, bend, crush, lie, or sit on heater eye pad to avoid damage to device.
9. **Never** pull this pad by the supply cord and do not use the cord as a handle.
10. **Carefully** examine before each use. Discard the pad if it shows any sign of deterioration (such as blistering or cracking).
11. **Only Use** the power supply provided by Digital Heat.
12. **Unplug** pad when not in use.
13. **DO NOT** tamper or modify the heater eye pad materials/configuration. There are no user serviceable parts. If for any reason this pad does not function satisfactorily, contact/return to Digital Heat.
14. **DO NOT** use this pad while wearing contact lenses or with other eye care products.
15. **DO NOT** use this pad with any liniments, salve, ointments, liquids, or any other materials in associations with the specified usage instructions of this device.
16. **DO NOT** use this pad while taking sensory dulling medication.
17. **DO NOT** wraps cord tightly or around eye pad heater to avoid damage to components. Loop cord lightly for storage to avoid any damage.
18. **Save** these instructions.

Contact Information:

Web:

www.heatedeyepad.com

Email:

Info@digitalheat911.com

Technical Support:

(512) 517-6649