

# M-Thumb


## Instructions for use



### M-THUMB

1. If using a silicone socket, apply a small amount of water to limb to dampen skin.
2. Push limb in to prosthesis.
3. Adjust thumb to desired position.
4. If more friction is required for flexion/extension, insert a phillips screw driver in to side of thumb to tighten screw.
5. If more rotation resistance is required, remove phillips screw, remove thumb to access allen head screw. Tighten screw to desired tension. Install thumb back in to position and re-install screw and tighten to desired tension.
6. New users should start with a wearing schedule and ease in to a full day of use.
7. To remove prosthesis, hold prosthesis with other hand and pull your limb out of prosthesis or consult with you prosthetist or therapist on proper technique.
8. Clean daily with fragrant free anti-bacterial wipes. Consult with your doctor or therapist on preferred cleansing wipes.

M-THUMB



Shown connected to prosthetic socket

M-THUMB



Shown connected to prosthetic socket

M-THUMB



PARTIAL HAND SOLUTIONS  
86 Ladyslipper Lane  
Southington, CT, USA 06489  
+1 (860) 538-5532



EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

