

RNA Isolation
Qiu Lab

Updated 6/2016

RNA Isolation

BEFORE STARTING

1. The first time you use the kit, you will need to make several solutions. All should be prepared to RNase/DNase free water
 - a. 70% EtOH
 - i. 35 mL Absolute (200 proof) ethanol
 - ii. 15 mL water
 - b. 80% EtOH
 - i. 40 mL Absolute (200 proof) ethanol
 - ii. 10 mL water
2. You will need to add 4 volumes EtOH (200 proof) to buffer RPE.
3. DNase I needs to be prepared and aliquoted.
 - a. Add 550 μ L of RNase/DNase free water to the lyophilized enzyme.
 - b. Mix via inversion. Do NOT vortex or sonicate.
 - c. Aliquot 90 μ L into 1.5 mL tubes. Store at 2°C for 6-8 weeks –20°C for 9 months.
4. Buffer RLT (provided with kit) is not used in this protocol; thus, you do not need to add β -ME to it.

NOTE: Steps 12 and forward are taken directly from the Qiagen manual (tissues section).

PROCEDURE

1. Collect necessary materials.
 - a. Qiagen RNAeasy Micro Kit (74004)
 - b. Samples (usu. Brain punches)
 - c. Dry Ice
 - d. 70% Ethanol**
 - e. 80% Ethanol**
 - f. Trizol (Invitrogen; 100 mL; 15596-026)
 - g. Chloroform (Sigma 500 mL; C-2432)
 - h. RNase free Pellet Pestle (Kimball-Kontes; 100/pk; 749521-1590)
 - i. Mini-Pestle Mixer/Vibrator
 - j. DNase aliquot (in being done).
2. Place samples on dry ice. Because of the speed necessary for some steps, doing more than 8 samples at 1 time is NOT recommended.
3. Add 200 μ L Trizol a sample and let thaw at RT for 30 seconds.
4. Homogenize tissue using mini-pestle and vibrator. Avoid generating excessive bubbles and change mini-pestle between groups.
5. Repeat 3/4 until all samples are homogenized.
6. Add 600 μ L Trizol to each sample and let stand **4 - 5 MINUTES**.
7. Add 160 μ L of chloroform to each tube. Shake vigorously for **15 seconds**.
8. Incubate at RT for **3 MINUTES**.
9. Spin samples at 12000g for **15 MINUTES** at **4°C**.
10. Label new tubes.
11. Remove the top (clear) layer from the samples and place in new tubes. Note the maximum volume would be 480 μ L, though 300-350 μ L is more reasonable. Because you will have to add a 1 volume of 70% ethanol in the next step and the column max is 700 μ L, using 350 μ L of the homogenized lysate is highly recommended.
12. Add 1 volume of 70% EtOH to the tube. Mix with the pipettor. Some precipitate may form.
13. Apply the sample, including any precipitate that may have formed, to an RNeasy MinElute Spin Column in a 2 mL collection tube. Close and centrifuge **15 SECONDS** at $\geq 8000g$. Discard flow-through.
14. Add 350 μ L of buffer RW1 to the column. Close the tube and centrifuge **15 SECONDS** at $\geq 8000g$. Discard the flow-through. Note: if skipping DNase treatment, use 700 μ L and then move to step 18.
15. Remove 1 aliquot of DNase (90 μ L). Add 630 μ L of buffer RDD. Mix via inversion.
16. Pipet 80 μ L of the DNase onto the column. Incubate **15 MINUTES**.
17. Pipet 350 μ L of RW1 into the column. Close and centrifuge **15 SECONDS** at $\geq 8000g$. Discard flow-through. Also discard the collection tube.
18. Place the column in a new 2 mL collection tube. Add 500 μ L of RPE buffer to each column. Close and centrifuge **15 SECONDS** at $\geq 8000g$. Discard flow-through.
19. Add 500 μ L of 80% EtOH to column. Close and centrifuge **2 MINUTES** at $\geq 8000g$ to dry the membrane. Discard the flow-through and the collection tube.
20. Transfer to column to a new 2 mL collection tube. Open the cap and spin at full speed for **5 MINUTES** to remove trace EtOH.
21. To elute, add 14 μ L of RNase/DNase free water. Place directly on the silica membrane. Spin at max speed for **1 MINUTE** to elute.
22. Read 230, 260, 280, ratios, and concentrations on NanoDrop.