




Pick&Place™ Holder Product Data Sheet


EXpressLO LLC
5483 Lee St Unit 12
Lehigh Acres, FL 33971 USA

+1-321-663-3806
www.EXpressLO.com
info@EXpressLO.com


The **Pick&Place™** holder kit facilitates *ex situ* lift out using the patented **EXpressLO™** grids and method for S/TEM or other site specific materials characterization [1-5]. The specially designed holder accepts up to a 25 mm diameter sample stub for the specimen **Pick**, and a dual pin grid holder for the **Place**. The screw drive adjusts the grid **Place** height relative to the specimen **Pick** height for optimized tip positioning and manipulation during the **EXpressLO™** lift out process. The **Pick&Place™** grid holder features dual pins at 90° with respect to each other; the first pin for positioning the grid flat for the **Place** and the second pin for positioning the grid upright for further FIB processing without additional grid handling. The grid holder can be used with pre-tilted holders or mounting stubs.



The **Pick&Place™** dual pin grid holder.




The patented **EXpressLO™** grid and method was used to **Pick** a specimen, **Place** it across an open slot in a backside orientation, and transfer back to the FIB for final thinning without additional handling using the dual pin **Pick&Place™** holder.




Level grid for easy lift out

The **Pick&Place™** holder kit.



Includes a grid stub holder for thin film coated 3 mm grids.



The patented **EXpressLO™** Cu half grid [5] with 15 specimen positions: 7 small numbered slotted openings, 8 large lettered slotted openings.

References:

- [1] Introduction to FIB, eds. L.A. Giannuzzi & F.A. Stevie Springer (2005).
- [2] L.A. Giannuzzi, Microsc. Microanal. 18, 2012, 632.
- [3] L.A. Giannuzzi, ISTFA 2012, ASM Int. 388.
- [4] L.A. Giannuzzi Microsc. Microanal. 19, 2013, 906.
- [5] www.YouTube.com/LAGiannuzzi.