



Rhombix

S E R I E S



Presented by: Cheryl Austin



Laboratory Systems Integration

The Rhombix Series of Products

Custom Software and Database Development

- Oracle (8.3 or higher, 9i preferred)
- MS-SQL 2000 (SP3 or higher)

Laboratory Automation and Informatics

- Scientific process definition
- Robotics to support protein crystallography
- Integration of selected third-party products

Imaging

- Image acquisition
- Image processing
- Image viewing and annotation



The Rhombix OPUS

A Scalable Complete Solution For Protein Crystallization




R h o m b i x O P U S


Developed in Collaboration with GlaxoSmithKline



RHOMBIX SERIES COMPONENTS



Experiment Optimization/Recursion



The goal is to aid users in *quickly* and *accurately* defining refined screen definitions based on individual initial screen well results without waste or duplication.

Benefits of a Data-centric Approach

- Frees scientists to focus on the science, not worry about the equipment
- Improves use of data and generated knowledge, presents data so that it is meaningful
- Provides a conduit for data mining
- Adapts to process changes
- Creates more data with increased value
- Allows for *flexible* automation
- Enables efficient scheduling of tasks
- Process output used as feedback

Automation Benefits

- Speed (saves time, increases throughput, minimal downtime)
- Quality
 - less contamination
 - improved accuracy: Coefficient of Variation $\leq 3\%$ (CV \cong 8-9% for non-automated preps)*
 - eliminate pipetting errors and other “human” factors
- Reproducibility and Repeatability
- Smaller sample sizes
 - Makes better use of reagents
 - Reduces costs in the long run
- Scalability
- Unattended operation
- Flexibility
 - More robust experimentation procedure, more complex preps possible
 - Rapid optimization
- Frees scientists from mundane tasks

* “Automated Sample Prep Grows Up”, in Laboratory Automation, Spring 2002, p. 5



The Rhombix Centrix Engine

- **PC-based control engine, written in C++, multi-tasking, multi-threaded**
 - ✓ Decouples process flow from lower-level control code
- **Ideal controller for high-throughput lab systems**
 - ✓ Collect the data from all sources and associate it with the original experiment definition
 - ✓ Connect third-party equipment (liquid handlers, sealer, shaker)
- **Control multiple pieces of equipment simultaneously and handle:**
 - ✓ Process scheduling / process flow
 - ✓ Equipment-level task queuing
 - ✓ Resource/equipment allocation and utilization
 - ✓ Plate data tracking
 - ✓ Error notification




Rhombix Screen Highlights

- **Manage reagents:** buffers, macromolecules, chemicals, additives, etc.
- **Create screens:** various scope and distribution patterns supported
- **Supports various experiment types:** sitting drop, hanging drop, microbatch, reagent prep (“daughter” plates)
- **Define experiment:**
 - Select screen, plate type, well volume
 - Select protein, shelf drop volume, additive, ratios of each
 - Select temperature, image schedule
 - Select process, modify process steps
- **Deck management** (inventory reagents in Q1-2004)
- **Store all information in relational database**
- **Uses this data to drive instrumentation:** dynamic task management
- **Multi-user, network-based, .NET application**



Rhombix OPUS





Reagent Manager - Chemicals

Reagent Information

Chemical Name:

Description:

☐ Custom ☐ Volatile

Categories

☐ Acid
☒ Additive
☐ Base
☒ Buffer
☐ Detergent
☐ Heavy Atom
☐ Ligand
☐ Oil
☐ Precipitant
☒ Salt
☐ Screening Solution
☐ Solvent

Chemical Version

Concentration: M pH:

Viscosity:

Formulations Inventory

Formulation (2 of 2)

☐ Primary ☐ Inactive

Components:

Name	Category	Final Conc	
Sodium Acetate (4M PH 4.6)	Additive	1.0 M	

Recipe:

User-defined values

Flexible formulation-based reagent definitions



Reagent Manager - Macromolecules

Reagent Information

Macromolecule Name:

Kinase-BTK

Description:

Classes

- ☐ GPCR
- ☐ Ion Channel
- ☒ Kinase
- ☐ Nuclear Receptor
- ☐ Protease
- ☐ Transporter
- ☐ Viral Protein

Macromolecule
versioning allows full
tracking of all samples

Target Details

Target ID:

BR-549

Sample ID:

12-3771A

Sample Description:

Macromolecule Version Detail

Inventory

Version: 30mg/ml at pH 4.6

☐ Inactive

Buffer: Sodium Acetate

Final Buffer Conc: 25mM at pH 4.6

Final Protein Conc: 30 mg/ml

Components:


Name	Description
Magnesium Formate	.02M

Ligands:

Name	Description
Ligand 5	15%

Version Description:

Resuspend 30mg lyophilized protein in 1ml 25mM NaAc which has had MgFormate and Ligand 5 added.



Screen Designer

Name: # Wells:

Description:

Type
☐ Custom Screen ☐ Additive Screen

Status
☒ Active ☐ Published

Grid Detail
Grid Comment:

Parent Screen:

Parent Well:

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Add ▾

Grid Recipe							
Color	Category	Name	Scope	Pattern	Pattern Value 1	Pattern Value 2	
Red	Additive	Sodium Acetate (1.0M pH 4.6)	GRID	CONSTANT	0.25		
Yellow	Precipita	Ammonium Sulphate (3.5M)	PER ROW	START/END	0.0	2.0	
Purple	Precipita	Ethylene Glycol (100.0%)	PER ROW	START/STEP	25.0	-5.0	



Experiment Manager - Setup

Experiment Name:	Experiment Type:	Description:
<input type="text" value="Coarse Screen Experiment 1"/>	<input type="text" value="Sitting Drop"/>	<input type="text" value="Coarse screen using Hampton CS1 and 2 screens"/>
Barcode:	Status:	Date Submitted:
<input type="text" value="⟨Deferred⟩"/>	<input type="text" value="Design*"/>	<input type="text"/>

Setup | Process

All volumes are expressed in uL (microliters)

Screen:	<input type="text" value="Hampton CS 1 and 2 from Tubes"/>	<input type="button" value="..."/>
Plate Type:	<input type="text" value="GREINER 1-SQ"/>	
Fill Vol.:	<input type="text" value="100"/>	Max Vol.: <input type="text" value="140"/>
Image Sched:	<input type="text" value="Cheryl's Schedule 1A"/>	
Temp:	<input type="text" value="4"/>	<input type="text" value="C"/>

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Drop Site Details (* Required fields)

Site	* Macromolecule	Vol	Reagent Additive	Screen Additive	* Ratio-[Pro:Pre:Add]	* Disp Order	* Drop Vol	MaxVol
1	Lysozyme in NaAc (30mg/ml)	0.167	Inositol (100.0%)		1 : 1 : 1	Precip-Prot-Adtv	0.5	4

Macromolecule Reagent
Row Detail

Reagent Additive Row
Detail

Screen Additive Row
Detail

Show Protein Totals

Experiment Summary



Experiment Manager - Process

Experiment Name:	Experiment Type:	Description:
<input type="text" value="Coarse Screen Experiment 1"/>	<input type="text" value="Sitting Drop"/>	<input type="text" value="Coarse screen using Hampton CS1 and 2 screens"/>
Barcode:	Status:	Date Submitted:
<input type="text" value="<Deferred>"/>	<input type="text" value="Design*"/>	<input type="text"/>

Setup Process

Process:

Process Steps			
Step	Mode	Machine	Notify
Get Plate	Robotic	Screen Robot	None
Screen-Making	Robotic	Hamilton MicroLab STAR	None
Move to Next Station	Robotic	Screen Robot	None
Shaking	Skip	Plate Shaker	None
Move to Next Station	Robotic	Screen Robot	None
Protein Drop-Making	Robotic	Cartesian HoneyBee	None
Move to Next Station	Robotic	Screen Robot	None
Sealing	Manual	Plate Sealer	None
Move to Next Station	Robotic	Screen Robot	None
Imaging/Storage	Robotic	Rhombix Vision (4C)	After
			None
			Before
			After
			Before and After

To:

Subject:

Msg Text:



Experiment Launcher

Select Experiments

Selec	Name	Description	Temp.	Date Submitted	Type	Username	Order
<input checked="" type="checkbox"/>	C Test		20	08/19/03	GREINER 1-SQ	DCA	1
<input checked="" type="checkbox"/>	TestScreen 2		20	08/19/03	GREINER 1-SQ	DCA	2
<input type="checkbox"/>	Test HCS 1 & 2 Reagent Prep			08/04/03	GREINER MAST	DCA	
<input type="checkbox"/>	Greg Test 451			08/06/03	GREINER MAST	DCA	
<input type="checkbox"/>	Greg Test 452			08/06/03	LINBRO	DCA	
<input type="checkbox"/>	Cheryl's Test Exp1		20	08/11/03	GREINER 3-SQ	AUSTIN	
<input type="checkbox"/>	Cheryl's Test Exp2		4	08/11/03	GREINER 3-SQ	AUSTIN	
<input type="checkbox"/>	Copy of Test HCS 1+2		4	08/12/03	GREINER MAST	AUSTIN	
<input type="checkbox"/>	Linbro Drop-Finding Plate		20	08/18/03	LINBRO	DCA	
<input type="checkbox"/>	CA Test Screen 1		20	08/18/03	GREINER 1-SQ	DCA	
<input type="checkbox"/>	Small Volume Experiment		4	08/19/03	GREINER 1-SQ	DCA	
<input type="checkbox"/>	Cartesian Additives Test		4	08/19/03	GREINER 3-SQ	DCA	
<input type="checkbox"/>	Test Additive Screen 1		20	08/19/03	GREINER 1-SQ	DCA	
<input type="checkbox"/>	Test Additive Screen 2		20	08/19/03	GREINER 1-SQ	DCA	

System Prep


Device	Deck Layout	Find Best Match	Edit Deck	Status	# of Items Not on Deck
▶ Hamilton MicroLab STAR	Hamilton Deck #1	Match	Edit	Ready	0
Cartesian HoneyBee	Additives Deck	Match	Edit	Missing Items	2



Deck Manager

Deck Items		
Name	Labware	Position
Magnesium Sulphate (4.0M)	✗	✗
Potassium Chloride (2.0M)	✗	✗
Zinc Acetate (1.0M)	✗	✗
Jeffamine M-600 (60.0%)	✗	✗
Sodium Iodide (1.0M)	✗	✗
Hampton CS 1 and 2 from Plate:	CAR_DW_1_PLATE_5	
1000uL Disposable	TIP_BASE_1_1	
10uL Disposable	✗	
300uL Disposable	TIP_BASE_1_3	
300uL Needle	WASH_BASE_1_1	
300uL Disposable	TIP_BASE_2_3	
PEG 6000 (60.0%)	CAR_TUBE_50ML_5	1
Citric Acid (1.0M pH 5.5)	CAR_TUBE_15ML_4	2

Experiment Items		
Name	Volume (uL)	Build





Process Monitor

Experiments:

Experiment Name	Barcode	Location	Owner	Status
Open House 3	<no barcode>		DCA	Aborted
Open House 4	000000LMA8	Screen_Robot	DCA	Built
Open House 5	<no barcode>	Screen_Robot	DCA	Built
Open House 5a	<no barcode>		DCA	Aborted
Open House 5b	<no barcode>		DCA	Aborted
Open House 5c	000000LMAB		DCA	Built
Open House 5d	000000LMAC	Screen_Robot	DCA	Running
s32	000000LM9Y		DCA	Built

Process Steps:

Process Name	Machine	Mode	Start Time	End Time
Get Plate	Screen Robot	Robotic	9/23/2003 5:08:29 PM	9/23/2003 5:09:18 PM
Screen-Making	Hamilton MicroLab STAR	Robotic	9/23/2003 5:09:18 PM	9/23/2003 5:12:06 PM
Move to Next Station	Screen Robot	Robotic	9/23/2003 5:12:06 PM	9/23/2003 5:12:55 PM
Shaking	Plate Shaker	Robotic	9/23/2003 5:12:55 PM	9/23/2003 5:14:58 PM
Move to Next Station	Screen Robot	Robotic	9/23/2003 5:14:58 PM	9/23/2003 5:18:17 PM
Protein Drop-Making	Cartesian HoneyBee	Robotic	9/23/2003 5:18:17 PM	9/23/2003 5:23:21 PM
Move to Next Station	Screen Robot	Robotic	9/23/2003 5:23:22 PM	9/23/2003 5:24:24 PM
Sealing	Plate Sealer	Robotic	9/23/2003 5:24:24 PM	9/23/2003 5:25:00 PM
Move to Next Station	Screen Robot	Robotic	9/23/2003 5:25:00 PM	
Imaging/Storage	Rhombix Vision (4C)	Skip		



Rhombix Vision XL




- **Supports various plate formats and geometries**
- **Storage enclosure**
 - Temperature controlled (4° C to 40° C)
 - Robotic random access & retrieval
 - 200 to 2000+ plates
- **Image Scheduling**
 - User-based image event scheduling
 - Dynamic task management, intelligent conflict resolution
- **Customizable Configurations**
 - Vision MS
- **Complete plate and well audit trail**



Rhombix Vision



- Multi-mode* plate imager
 - Darkfield, Brightfield, and Polarized Illumination
- Custom Lighting Matched to Optics to Mitigate Inverse “Halo” Effect
- Image formats: bmp, png, tif, jpg
- True Color CCD Imaging
 - Image resolution: 1360 x 1036 pixels
 - Auto Focus, Auto Exposure, White Balance
- Image Processing:
 - Xfocus (extended focus): composite image created from multiple slices
 - Automatic Drop Finding: zoom in on drop to full field of view
 - Future: Clear drop discernment, crystal and precipitate detection
- Linbro & SBS Plate Formats
- Auto and Microscope Modes
- Data fusion
 - User-defined image annotation/scoring
 - Automated machine scoring



Clarix

Image Viewing & Annotation

- Windows-based desktop application freely distributable to users' desktop PCs
- Annotation parameters (descriptors and modifiers) are user-definable
- Images are accessed a plate at a time, via unique plate identifier (barcode)
- Images may be viewed as thumbnails (small, medium, large) and "full view" (one at a time)
- View progression of images for same well over time (different image events)
- Infinite zoom
- Measure objects
- Images and annotations may be exported
- Images may be deleted (annotations always kept)
- Will be re-written in .NET framework in Q1-2004



Schedule Template

Schedule

Template: 4 events 2-hrs. apart, brite only

Events

Event Type	Relative Start Time
Image	On Entry
Image	+ 2 hours
Image	+ 4 hours
Image	+ 6 hours
Removal	+ 8 hours

Insert

Delete

☐ Skip This Event

Notification

☐ Notify

To:

Subject:

Msg Text:

Image types and treatment

- | | | |
|-------------------------------------------------|-----------------------------------------|-----------------------------------------|
| <input checked="" type="checkbox"/> Brightfield | <input type="checkbox"/> Autofocus | <input type="checkbox"/> XFocus |
| <input type="checkbox"/> Darkfield | <input type="checkbox"/> Auto Exposure | <input type="checkbox"/> White Balance |
| <input type="checkbox"/> Polarized Extinction | <input type="checkbox"/> Auto Drop Find | <input type="checkbox"/> Semi-Automatic |

Image Type

☒ Color

☐ Grayscale

Format: JPEG Large



Image Scoring (Clarix)

Plate ID: LMA4 ☒ Sharpen Image

Image (Brightfield, 1 of 1): Plate Type: RIK - LINBRO 1-DROP

9/29/2003

	No Score	Some Scored	All Scored			
	01	02	03	04	05	06
A						
B						
C						
D						

Well: B01

Image Scoring

Annotation:

- ☐ (C) Clear
- ☐ (P) Precipitate
- ☒ (S) Crystal
- ☐ (O) Other

Modifier:

- ☐ Microcrystals
- ☐ Microcrystals?
- ☐ Crystals?
- ☐ Needles
- ☐ Plates > 50 x 50 micr
- ☐ Cubes > 50 x 50 micr

Comments:

Image: \\gemin\EDrv\Images\030929\LMA4\LMA4_B01-0_BRITE_030929 Image Taken: September 29, 2003 11:27:04

Buttons: Delete This Image, Delete Clear Images, Refresh, Export, Save

Flexible user-defined scoring descriptors and modifiers



Rhombix Manager - Automation Monitor

Equipment				
Name	Description	Device	Subsystem	Detail
Rhombix_Screen	Rhombix Screen			
Screen_Robot	Screen Robot			
Hamilton_MLSTA	Hamilton MicroLab STAR			
Plate_sealer	Plate Sealer			
Cartesian_HB	Cartesian HoneyBee			
Plate_shaker	Plate Shaker			
Rhombix_Vision	Rhombix Vision (4C)			
Vision_Robot	Vision Robot			
Matrix	Matrix Storage			
Manual_Stack	Manual Stack			
Hotel	Hotel			
Imager	Imager			
Rhombix_Vision	Rhombix Vision (20C)			
Vision_Robot	Vision Robot			
Matrix	Matrix Storage			
Manual_Stack	Manual Stack			
Hotel	Hotel			
Imager	Imager45			
Hotel	Sealer Spacer Block Hotel			



Rhombix Manager – Vision Monitor

Vision Unit: Rhombix Vision (20C) ▼

Experiments & Plates | Scheduled Image Events | Location Detail

Experiments / Plates

Experiment Name ▲▼	Barcode ▼	Location ▼	Owner ▼	Create Date ▼
CA-G3-from DW	000000LM94	Kendro Cytomat Hotel St	Cheryl Austin	11/18/2003 10:12:29 AM
DAUGHTER TEST M AM 1	000000LM91	Kendro Cytomat Hotel St	Joe Christopher	11/17/2003 10:04:35 AM
DAUGHTER TEST M PM 2	000000LM92	Imager Nest	dca dca	11/17/2003 4:19:39 PM

Scheduled Events

Event Type ▼	Scheduled ▲▼	Started ▼	Ended ▼	Image Setting ▼	Image Type ▼	Image Format ▼
Removal	11/18/2003 4:40:22 PM	11/18/2003 4:40:37 PM	11/18/2003 4:40:37 PM	Color	Brightfield	
Image	11/18/2003 4:42:32 PM	11/18/2003 4:43:04 PM	11/18/2003 4:52:59 PM	Color,Extended Fo	Brightfield	
Image	11/18/2003 5:42:32 PM		11/19/2003 9:25:41 AM	Color	Brightfield,Darkfield	
Image	11/18/2003 10:42:32 PM		11/19/2003 9:25:41 AM	Auto Focus,Color	Brightfield	
Image	11/19/2003 3:42:32 AM			Color,Extended Fo	Brightfield	



Rhombix Manager – Vision Monitor

Vision Unit: Rhombix Vision (20C) ▼

Experiments & Plates | Scheduled Image Events | Location Detail

☒ Show Occupied Location Only

Location Detail

Location				# Positions	# Plates	
Imager Nest				1		
Row	Column	Status	Experiment Name	Barcode	Event Type	Scheduled
1	1	Enabled	DAUGHTER TEST M PM 2	000000LM92	Image	11/18/2003 6:05:05 PM
Location				# Positions	# Plates	
Kendro Cytomat Hotel Storage Nest				504		
Row	Column	Status	Experiment Name	Barcode	Event Type	Scheduled
19	1	Enabled	CA-G3-from DW	000000LM94	Image	11/19/2003 3:42:32 AM
20	1	Enabled	DAUGHTER TEST M AM 1	000000LM91	Image	11/18/2003 6:01:30 PM
Location				# Positions	# Plates	
Kendro Cytomat Hotel Transfer Nest				1		
Linbro Matrix Storage Storage Nest				70		
Manual Input/Output Hotel Nest				8		
Rhombix Vision (20C) Automatic Input Transfer Nest				1		
Vision Robot Gripper				1		



Query Tool

Note: Wildcard character = %, Escape character = \

Macromolecule Name:

Macromolecule Version:

Macromolecule Sample:

Macromolecule Target:

Macromolecule Class: Or

Well Chemical Name: Or

Concentration in Well: to

pH in Well: to (1.0-14.0)

Well Chemical Cat: Or

Drop Additive Name:

Drop Additive Conc: to

Drop Additive pH: to (1.0-14.0)

Results Descriptor: Or

Results Modifier:

Experiment Name:


Experiment Submit Date: to

Temperature: C to C

☐ Include well compartment ID ☒ Current user only

☐ Include well and shelf info

Query



Query Results (16 rows)

RowID	Chemical	Hits
4	Aluminum Phosphate	4
5	PEG2000	4
6	Ammonium Chlorhydrate	5
7	Sodium Chlorhydrate	6
8	Sodium HEPES	7
9	Tris Hydrochloride	8
16	Sodium Cacodylate	11
15	Imidazole	14
14	PEG800	15
13	Sodium Formate	18