

# Product Specification Sheet



## Functions:

The RehabPro Closed Chain Trainer makes early weight shifting, early weight bearing, early squats and step ups/down safe by reducing the amount of body weight applied to the lower extremities. One can get an objective measurement of actual body weight and body weight distribution by placing two weight scales on the 24"x 24" foot board. The exercises are dosed by raising the sled on the ladder as the tolerance/function improves. The joint excursion at the hips, knees and ankles is controlled by "blocking" the sled from returning all the way to its lowest position by a "marine rope stop". Rehab tools of different density can be placed on the foot board to further enhance the functional quality of coordination.

## Specifications:

**Vendor Parts Number:** 13125 w/ ladder

**Construction:** 1"x 1.5" – 14-gauge, high quality, powder coated steel frame for sturdiness & safety. The frame supports patients up to 400 lbs. Partial weight bearing starts at 5% body weight and gradually increases to 80% body weight at the highest rung on the ladder. The sled rides on 4x-2" high end ball-bearing wheels, all guided to stay on track for safety purposes. The two aluminum hooks at the head end of the sled attaches to 1" rungs. The ¾" thick, birch plywood foot board has a non-slip rubber surface.

**Assembly:** No assembly required.

**Accessories:** Ladder included  
**Dims:** 84"x 18"x 2", 45 lbs.

**Product Dimensions:** 13125  
**Dims:** 87"x 22"x 31", 125 lbs.

**Shipping Dimensions:** 13125 w/ ladder  
**Dims:** 91"x 24"x 35", 265 lbs.

SHIPS VIA TRUCK 

## Warranty:

1 year bumper to bumper on all parts.  
5 years on all moving parts.  
Lifetime on frame.