

# KimanTech Technical Application: 1-Step vs 2-Step PCR for Pathogen DNA Detection from Whole Blood

## Background

Pathogen detection by PCR analysis has become an area of significant interest and product development. Commercial systems exist that can detect panels of pathogens using both one-step and two-step PCR approaches but utilize enclosed systems (a sealed consumable). While enclosed systems are ideal for point-of-care applications, they are typically expensive, inflexible, and do not utilize existing laboratory equipment. Reagent kits exist that are compatible with standard laboratory equipment, but all use single-step PCR to avoid the added steps and contamination risks associated with two-step PCR. Since the Alluvia system can enable the development of two-step PCR pathogen detection reagent kits for standard laboratory equipment, we investigated if the added PCR step will significantly improve PCR results.


## Experimental Design and Methods:

To compare two-step and one-step PCR approaches for pathogen detection, we have developed a two-step PCR panel to differentiate eight major pathogenic species of the yeast *Candida*, a causative agent of serious and often fatal blood infections. Four primers were designed to amplify the 26S ribosomal RNA gene from each strain of *Candida*, two outer primers for the first round of PCR and two nested inner primers for the second round of PCR. Each strain of *Candida* was cultured, and genomic DNA purified to be used as template. The assay was performed as follows.


- Applied Biosystems™ (ABI) Power SYBR™ Green PCR master mix (used for all PCR reactions) along with all eight outer primers pairs were added to a single PCR tube, then 1000 copies of genomic DNA from one of the *Candida* species added.
- This primary multiplex PCR reaction was performed for 40 cycles using the ABI StepOne™ instrument (used for all PCR reactions).
- The next round of qPCR was performed by adding master mix and a different pair of nested primers to each well of an 8-well strip.
- To compare one-step to two-step qPCR, either 1000 copies of genomic DNA from one of the *Candida* species or a 100X dilution of the first-round multiplex PCR products were used as templates.
- The second round of qPCR was also performed for 40 cycles and a melting analysis performed.
- PCR products were also visualized on an electrophoresis gel.

## Results:

### Example qPCR melting curve results (from *C. guilliermondii*)



### Example gel electrophoresis data



In this example, background signals exist with both one-step and two-step PCR methods. However, results with two-step PCR are easier to interpret.

### Two-step qPCR specificity summary (based on visual analysis)

Nested PCR Result	Target Template Input (gDNA - 1000 copies for 9-plex preamplification)								Non-Target (gDNA - 100000 copies)		
	<i>C. albicans</i>	<i>C. auris</i>	<i>C. kruei</i>	<i>C. glabrata</i>	<i>C. guilliermondii</i>	<i>C. lusitaniae</i>	<i>C. parapsilosis</i>	<i>C. tropicalis</i>	<i>C. dubliniensis</i>	<i>S. cerevisiae</i>	<i>S. pombe</i>
<i>C. albicans</i> PM	Positive	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.
<i>C. auris</i> PM	Neg.	Positive	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.
<i>C. krusei</i> PM	Neg.	Neg.	Positive	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.
<i>C. glabrata</i> PM	Neg.	Neg.	Neg.	Positive	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.
<i>C. guilliermondii</i> PM	Neg.	Neg.	Neg.	Neg.	Positive	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.
<i>C. lusitaniae</i> PM	Neg.	Neg.	Neg.	Neg.	Neg.	Positive	Neg.	Neg.	Neg.	Neg.	Neg.
<i>C. parapsilosis</i> PM	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.	Positive	Neg.	Neg.	Neg.	Neg.
<i>C. tropicalis</i> PM	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.	Neg.	Positive	Neg.	Neg.	Neg.

### One-step qPCR specificity summary (based on visual analysis)

Nested PCR Result	Target Template Input (gDNA - 1000 copies)								Non-Target (gDNA - 1000 copies)		
	<i>C. albicans</i>	<i>C. auris</i>	<i>C. kruei</i>	<i>C. glabrata</i>	<i>C. guilliermondii</i>	<i>C. lusitaniae</i>	<i>C. parapsilosis</i>	<i>C. tropicalis</i>	<i>C. dubliniensis</i>	<i>S. cerevisiae</i>	<i>S. pombe</i>
<i>C. albicans</i> PM	Positive	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk
<i>C. auris</i> PM	Cross-talk	Positive	Neg.	Neg.	Cross-talk	Neg.	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk
<i>C. krusei</i> PM	Cross-talk	Neg.	Positive	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk
<i>C. glabrata</i> PM	Cross-talk	Cross-talk	Cross-talk	Positive	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Neg.
<i>C. guilliermondii</i> PM	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Positive	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk
<i>C. lusitaniae</i> PM	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Positive	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk
<i>C. parapsilosis</i> PM	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Positive	Cross-talk	Cross-talk	Cross-talk	Cross-talk
<i>C. tropicalis</i> PM	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Cross-talk	Positive	Cross-talk	Cross-talk	Cross-talk

**Summary:** Two-step PCR significantly improves specificity allowing for more robust distinction of closely-related species of *Candida*.