



Research

General Information

ZIRC

Home ■ Genes / Markers / Clones ■ Expression ■ Antibodies ■ BLAST ■ Mutants / Tg ■ Anatomy ■ Maps ■ Publications

Site Search:

CHAPTER 10 - RECIPES

[ABC Solution \(Avidin-Biotin HRP Complex\)](#)

[Agar/Sucrose](#)

[Atabrine Stock Solution](#)

[Blocking Solution for Western Blots](#)

[BSA/PBS/DMSO](#)

[BT Fix](#)

[Chorion Removal](#)

[DAB \(Diaminobenzidine\)](#)

[DAB Heavy Metal Stain](#)

[DAB Presoak Solution](#)

[DAB Solution](#)

[DNA Extraction Buffer](#)

[EDTA](#)

[Egg Water](#)

[Embryo Extract](#)

[Embryo Medium](#)

[Epon](#)

[Epon/Araldite](#)

[Finquel \(MS222\)](#)

[Fish Water](#)

[Fix Buffer](#)

[2X Fix Buffer](#)

[Fixatives](#)

[Gelatin Embedding Medium](#)

[Genomic DNA Extraction Buffer](#)

[GHCI Buffer](#)

[Giemsa](#)

[Ginzburg Fish Ringers](#)

[GIT Buffer](#)

[Growth Medium](#)

[Hank's \(Final\)](#)

[Hank's \(Full Strength\)](#)

[Hank's Premix](#)

[Hank's Stock Solutions](#)

[Heat Shock \(Diploid eggs\)](#)

[In situ Hybridization Staining Buffer](#)

[Paramecia Seed Cultures](#)

[Paramecia for Baby Fish \(Traditional Method\)](#)

[Paramecia for Baby Fish \(Streamlined Procedure\)](#)

[Paramecia Medium Stock Solution I](#)

[Paramecia Medium Stock Solution II](#)

[Paramecia Storage Medium \(Long Term\)](#)

[PBS](#)

[2X PBS](#)

[PBS/BSA/DMSO](#)

[PCR Extraction Buffer](#)

[PMSE](#)

[PO4 buffer \(0.1 M, pH 7.3\)](#)

[Pronase](#)

[Protein Extraction Buffer](#)

[PTU](#)

[Ringer's Solutions](#)

[Salt Stock](#)

[SDS Sample Buffer](#)

[Schönefeld's Medium for Growing Paramecia](#)

[Sodium Bicarbonate](#)

[Sperm Freezing Medium](#)

[Sperm Thawing Medium](#)

[Stock salts](#)

[Subbing Solution for Slides](#)

[TBS](#)

[TES buffer](#)

[Tricaine](#)

[TTBS](#)

[Western Blot color development buffer](#)

[Western Blot color development solution](#)

[Whole-Mount in situ Hybridization Solutions](#)

CHAPTER 10 - RECIPES

ABC Solution (Avidin-Biotin HRP Complex):

3ml 1% DMSO, 0.1% Triton, PBS
30µl Vectastain solution A
30µl Vectastain solution B

Agar/Sucrose:

1.5% agar
5% M sucrose
Boil into solution and store as 1.5 ml aliquots at 4°C.

Atabrine Stock Solution:

10 mg/ml dH₂O. Store in a light tight bottle.

Blocking Solution for Western Blots:

3% dried milk in TBS

BSA/PBS/DMSO:

1% bovine serum albumin, 1% dimethyl sulfoxide in PBS.

BT Fix:

10 ml 10% para-formaldehyde
15 ml 1.25X fix buffer
Adjust pH to 7.3 if necessary. Add dH₂O to a final volume of 25.0 ml. Final concentration is 4% para-formaldehyde, 0.15 mM CaCl₂, 4% sucrose in 0.1 M PO₄ buffer.

Chorion Removal:

Drain eggs add 5 ml of 0.5 mg/ml pronase for 3.5 min.
Dilute eggs with 200 ml of 8x water. Rinse 3x more with 200 ml washes of 8x water.
8x water = 12 ml stock salts per liter dH₂O

Stock salts = 40 g Instant Ocean per liter dH₂O

DAB (Diaminobenzidine):

DAB is prepared by dissolving the contents of 1 g bottle of DAB in 25 ml dH₂O and filtering with a disposable unit attached to a disposable syringe. Aliquots containing 1 mg DAB per 25 µl of solution are stored in 1.5 ml microfuge tubes in a -20°C freezer. Thawed aliquots may be refrozen. To minimize risk in handling this carcinogen, use disposable plastic containers and pipettes wherever possible. Treat waste with bleach.

DAB Heavy Metal Stain:

Presoak solution

4 mg DAB
5 ml 0.6% Ni(NH₄)₂(SO₄)₂
5 ml 0.2 M Tris buffer, pH 7.4
50 µl DMSO

After presoak interval add

10 µl 3% H₂O₂

Final concentrations are 0.04% DAB, 0.1 M Tris buffer, 0.3% Ni(NH₄)₂(SO₄)₂ and 0.003% H₂O₂.

DAB Presoak Solution:

25 µl (=1 mg) DAB stock solution
1 ml 0.1 M PO₄ Buffer, pH7.3
1 ml dH₂O
20 µl DMSO
After presoak interval add
5 µl 3% H₂O₂
Final concentrations are 0.05% DAB, 0.05 M PO₄, 1% DMSO and 0.004% H₂O₂.

DAB Solution:

0.05% diaminobenzidine
1% dimethyl sulfoxide in 0.05 M PO₄ buffer, pH 7.3.

DNA Extraction Buffer:

10 mM Tris pH 8.2
10 mM EDTA
200 mM NaCl
0.5% SDS
200 µg/ml proteinase K

EDTA:

Stock - 10 mM EDTA, pH 7.0. Add 1 ml of stock/10 ml Ringer's (final conc. 1 mM EDTA)

Egg Water:

1.5 ml stock salts added to 1 L distilled water = 60µg/ml final concentration.

Embryo Extract:

1. Chill 200 3 d embryos after removing from chorions.
2. Rinse in 0.5% chilled bleach for 2 min and then in zero calcium Ringer for 2 min.
3. Transfer to a Dounce homogenizer with a minimum of liquid and homogenize well.
4. Resuspend in 1 ml L-15 supplemented with 0.3 mg/ml glutamine, 50 U/ml penicillin, 0.05 mg/ml streptomycin and 0.8 mM CaCl₂.
5. Store at -20°C.

Embryo Medium:

1.0 ml [Hank's Stock #1](#)
0.1 ml [Hank's Stock #2](#)
1.0 ml [Hank's Stock #4](#)
95.9 ml dd H₂O
1.0 ml [Hank's Stock #5](#)
1.0 ml fresh [Hank's Stock #6](#)
Use about 10 drops 1 M NaOH to pH 7.2

Epon:

12 g Epon 812
24.7 g Dodecenyl succinic anhydride
Weigh components into a 50 ml plastic beaker.
Cover with parafilm and stir on a magnetic stirrer for
15 min. Add 0.5 ml DMP-30

Continue stirring for another 15 min. Store in syringes in -20°C freezer.

Epon/Araldite:

24.7 g Epon 812
33.25 g Dodecenyl succinic anhydride
31.05 g Araldite 506

Mix well by stirring. Use disposable syringes to add:
2.3 ml Dibutylphthalate
2.5 ml DMP-30

Continue to mix well by stirring. Try to avoid incorporating excess air. If this occurs resin may be degassed by subjecting it to a mild vacuum. Store resin in 10 ml syringes at -20°C.

Finquel (MS222):

See [Tricaine](#). (Finquel is the trademark brand of Argent Chemical Laboratories, Inc.)

Fish Water:

60 mg "Instant Ocean" per liter dH₂O.

Fix Buffer:

Dilute 1.25X fix buffer 3:2 with dH₂O.

1.25X Fix Buffer:

1.0 g Sucrose
18.75 μ l 0.2 M CaCl₂
5 ml 0.5 M PO₄ buffer, pH 7.3
Check pH. Adjust if necessary to 7.3 with 1 M NaOH or HCl.
Add H₂O to a final volume of 15 ml.

Fixatives:

Fix buffer: 4% sucrose, 0.15 mM CaCl₂, 0.1 M PO₄ pH 7.3.
For general fixation: 1.5% glutaraldehyde, 0.5% paraformaldehyde in fix buffer.
For antibody staining: 4% paraformaldehyde in fix buffer.

Gelatin Embedding Medium:

17% gelatin in 10% Hank's saline.

Genomic DNA Extraction Buffer:

10 mM Tris pH 8
100 mM EDTA pH 8
0.5% SDS
200 μ g/ml Proteinase K

GHCI Buffer:

7.5 M guanidinium hydrochloride
0.025 M NaOAc pH 7.0
5 mM dithiothreitol
0.5% N-laurylsarcosinate.

Giemsa:

4 mls Giemsa Stock (Sigma Diagnostics)
4 mls 0.5 M Na Phosphate pH7
200 mls distilled water

Ginzburg Fish Ringers:

6.5 g NaCl
0.25 g KCl
0.3 g CaCl₂ (0.4 g CaCl₂•2H₂O)
Add ddH₂O to almost 1 liter
0.2 g NaHCO₃
Add ddH₂O to 1 liter
Note: The order of addition is important to prevent precipitation

GIT Buffer:

4 M guanidinium isothiocyanate
0.1 M Tris-HCl pH 7.5
1% β -mercaptoethanol.

Growth Medium:

L-15 (Sigma)
0.3 mg/ml glutamine
50 U/ml penicillin
0.05 mg/ml streptomycin
0.8 mM CaCl₂
10% embryo extract
3% fetal calf serum.

Hank's (Final):

9.9 ml Hank's Premix
0.1 ml Stock #6

Hank's (Full Strength):

0.137 M NaCl
5.4 mM KCl
0.25 mM Na₂HPO₄
0.44 mM KH₂PO₄
1.3 mM CaCl₂
1.0 mM Mg SO₄
4.2 mM NaHCO₃

Hank's Premix:

Combine the following in order:
10.0 ml Solution #1
1.0 ml Solution #2
1.0 ml Solution #4
86.0 ml ddH₂O
1.0 ml Solution #5

Store Hank's Premix in the refrigerator along with the Hank's solutions.

Hank's Stock Solutions:**Stock #1**

8.0 g NaCl
0.4 g KCl
in 100 ml dd H₂O

Stock #2

0.358 g Na₂HPO₄ Anhydrous
0.60 g KH₂PO₄
in 100 ml ddH₂O

Stock #4

0.72 g CaCl₂
in 50 ml ddH₂O

Stock #5

1.23 g MgSO₄·7H₂O
in 50 ml dd H₂O

Stock #6

0.35 g NaHCO₃
10.0 mls dd H₂O

Heat Shock (Diploid eggs):

T = 0 Fertilize eggs with UV-irradiated sperm
T = 5 min Transfer embryos to 28.5°C cylinder
T = 13 min Move cylinder from 28.5°C to 41.4°C
T = 15 min Move cylinder from 41.4°C to 28.5°C

In situ Hybridization Staining Buffer:

100 mM Tris pH 9.5
50 mM MgCl₂
100 mM NaCl
0.1% Tween-20
1 mM Levamisol (add fresh)

Paramecia Seed Cultures:

1. Add 10-15 grains of boiled wheat to 175 mls of dH₂O.
2. Inoculate with 20 ml from an excellent existing seed culture dish or with a sample from the commercial inoculant.
3. Grow for 7 to 12 days before using.

Paramecia for Baby Fish (Traditional Method):

1. Use glass finger bowls filled two-thirds full with system water.
2. To each bowl add 8-9 grains of of boiled wheat and 8 ml from a paramecia seed culture.
3. Stack the finger bowls 6 or 8 high, cover, and store at 28.5°C on well-lit shelves.
4. After 10-14 days, the paramecia are ready to feed to the fish larvae. Cultures remain useable and healthy for a month or more.

Paramecia for Baby Fish (Streamlined Procedure):

1. Fill plastic mouse cage or equivalent with 2 liters of system water.
2. Add a large pinch (about 40-50 grams) of boiled wheat, four 250 mg tablets of brewer's yeast, and 100 ml (or half the contents of a plastic seed culture dish) of paramecia culture to each mouse cage.

3. Cover and store in a warm, well-lit place. The covered mouse cages may be stacked three layers high with the highest layer closest to a warm light. It will be ready to feed first.
4. Brewer's yeast may also be added to the seed cultures (one half of a 250 mg tablet per 200 ml of seed culture). The cultures are then ready to use in four days and are depleted after two weeks.

Paramecia Medium Stock Solution I:

Component - g/l in dH₂O - Source

Calcium pantothenate - 1 - Sigma P-2250

Nicotinamide - 1 - Sigma N-3376

Riboflavin - 1 - Sigma R-9881

Pyridoxamine HCL - 1.16 - Sigma P-9158

Folic Acid - 0.5 - Sigma F-7876

Thiamine HCL - 3 - Sigma T-4625

Biotin - 0.00125 - Sigma B-4639

(make 100x stock) Store at -20°C the riboflavin is dispersed through-out the solution and will settle to the bottom. Shake well before using

Paramecia Medium Stock Solution II:

Component - g/l in ETOH - Source

Stigmasterol - 2 - Sigma S-6126

Lipoic Acid (DI-6,8-thioctic acid) - 0.02 - Sigma T-1395

Palmitic Acid - 3 - Sigma P-5917

Stearic Acid - 2 - Sigma S-4751

Oleic Acid - 0.4 - Sigma O-4379

Linoleic Acid - 0.2 - Sigma L-1376

Linolenic Acid - 0.06 - Sigma L-2376

(Note: Blow nitrogen gas over the solution and store at -20°C)

Paramecia Storage Medium (Long Term):

Component - Stock - Amount (ml)

NaH₂PO₄ - 0.01 M - 10

Na₂HPO₄ - 0.01 M - 10

Sodium citrate - 0.1 M - 20

CaCl₂ - 0.1 M - 15

dH₂O - 945

PBS:

0.8% NaCl

0.02% KCl

0.02 M PO₄, pH 7.3.

2X PBS:

8.0 g NaCl

0.2 g KCl

200 ml 0.1 M PO₄ Buffer, pH 7.3

300 ml dH₂O

When diluted 1:1 with dH₂O, final concentrations are 0.8% NaCl, 0.02% KCl and 0.02 M PO₄.

PBS/BSA/DMSO:

50 ml 2X PBS, pH 7.3
1 g BSA
1 ml DMSO

Check pH with paper and adjust if necessary. Add dH₂O to a total volume of 100 ml

PCR Extraction Buffer:

10 mM Tris pH 8
2 mM EDTA
0.2% Triton X-100,
200 µg/ml Proteinase K

PMSF:

Stock - 100 mM phenylmethylsulfonylfluoride in isopropanol. Immediately before use, add 30 µl of stock/10 ml Ringer's (final conc. 0.3 mM PMSF).

PO₄ buffer (0.1 M, pH 7.3):

77 ml 0.1 M NaH₂PO₄ (13.8 g NaH₂PO₄·xH₂O/liter dH₂O)
23 ml 0.1 M Na₂HPO₄ (14.2 g Na₂HPO₄/liter dH₂O)

Pronase:

5 mg/ml pronase diluted to 1 mg/ml in embryo medium

Protein Extraction Buffer:

10 mM tris, pH 7.4
2% Triton-X 100
1 mM PMSF
1 mM aprotinin
1 mM leupeptin
1 mM trypsin inhibitor

PTU:

0.003% 1-phenyl-2-thiourea in 10% Hank's saline.

Ringer's Solutions:

Normal

116 mM NaCl
2.9 mM KCl
1.8 mM CaCl₂
5 mM HEPES, pH 7.2.

High calcium

116 mM NaCl
2.9 mM KCl
10 mM CaCl₂
5 mM HEPES, pH 7.2.

Calcium free

116 mM NaCl
2.9 mM KCl
5 mM HEPES, pH 7.2.

Salt Stock:

20 tablespoons (280 g) Instant Ocean Sea Salts (Aquarium Systems, Inc.) dissolved in 2 liters distilled water

SDS Sample Buffer:

0.63 ml 1M Tris-HCl, pH 6.8
1.0 ml glycerol
0.5 ml β-mercaptoethanol
1.75 ml 20% SDS
6.12 ml H₂O
(10 ml total)

Store at -20°C in aliquots.

Schönefeld's Medium for Growing Paramecia:

Component - Concentration - Source
Powdered Skim Milk - 8.5 g/l - Any grocery store
Ribonucleic Acid - 1.0 g/l - Sigma R-6625
MgSO₄ - 0.5 g/l - Sigma M-2643
Phosphatidylcholine - 250 mg/l - Sigma P-5638

Stock Solution I: 5 ml/l
Stock Solution II: 2.5 ml/l

Sodium Bacarbonate

0.35 g NaHCO₃
10.0 ml ddH₂O

Sperm Freezing Medium:

9 ml Ginzburg Fish Ringers
1 ml Methanol
1.5 g Carnation Powdered Skim Milk
Note: This order is important to prevent precipitation of milk

Sperm Thawing Medium:

a. Measure the mm of sperm + freezing medium (sm) in the capillary.
b. Convert to volume:

$$10\mu\text{l (cap. vol.)} \times \text{mm sm} = \text{vol } (\mu\text{l})$$
$$90 \text{ mm (cap. length)} \quad 1$$

c. Multiply the calculated volume x 10 and use that much 10% Hank's saline for thawing.

Stock Salts

40 g "Instant Ocean" Sea Salts added to 1 L distilled water

Subbing Solution for Slides:

0.5% gelatin
0.05% $\text{CrK}(\text{SO}_4)\text{x}12\text{H}_2\text{O}$.

TBS:

20 mM Tris, pH 7.5
500 mM NaCl

TES buffer:

1% sodium dodecyl sulfate
5 mM EDTA
10 mM Tris-HCl pH 7.4.

Tricaine:

Tricaine (3-amino benzoic acid ethyl ester also called ethyl m-aminobenzoate) comes in a powdered form from Sigma (Cat.# A-5040). It is also available as Finquel (Part No. C-FINQ-UE) from Argent Chemical Laboratories, Inc. Make tricaine solution for anesthetizing fish by combining the following in a glass bottle with a screw cap:

400 mg tricaine powder
97.9 ml DD water
~2.1 ml 1 M Tris (pH 9).

Adjust pH to ~7. Store this solution in the freezer. (Buy the smallest amount possible because tricaine gets old.)
To use tricaine as an anesthetic combine the following in a 250 ml beaker:

4.2 ml tricaine solution
~100 ml clean tank water.

TTBS:

20 mM Tris, pH 7.5
500 mM NaCl
0.05% Tween-20

Western Blot color development buffer:

102 mg MgCl_2
4.2 g NaHCO_3
500 ml dH_2O

Western Blot color development solution:

66 μ l nitroblue tetrazolium (NBT) stock

33 μ l 5-bromo-4-chloro-3-indolyl galactopyranoside (BCIP) stock

10 ml color development buffer

Both NBT and BCIP stocks are 50 mg/ml in dimethylformamide. NBT stock is made by suspending 50 mg NBT in 700 ml dimethylformamide. Vortex. Add 300 μ l distilled water to dissolve. Store both stocks at 4°C in dark.

Whole-Mount in situ Hybridization Solutions:

PBST - PBS plus 0.1% Tween

SSCT - SSC plus 0.1% Tween

HYB* - 50% formamide, 5xSSC, 0.1% Tween-20

HYB+ - HYB* with 5mg/ml torula (yeast) RNA, 50 μ g/ml heparin

[Home](#) ■ [About ZFIN](#) ■ [Citing ZFIN](#) ■ [Glossary](#) ■ [Help and Tips](#) ■ [Contact ZFIN](#)

[Disclaimer, limitations, copyright](#) © University of Oregon, 1994-2008, Eugene, Oregon.
ZFIN logo design by Kari Pape, University of Oregon