# AS300RX RECEIVER ALARM & MONITORING SYSTEM

# Monitor and Protect

# Irreplaceable Samples / Valuable Assets

## SIMPLE INSTALLATION

It's easy to install the AS300 system. Fit the temperature probes into the equipment from the AS200TX transmitters and set required alarm values. Then program the settings into the AS300RX from a mobile phone. No phone line or line rental is required as the system uses a mobile network of your choice.

## **MULTIPLE SITES**

Multiple linked sites can be alarmed and monitored across the country as the system uses a GSM network. There's no limit to distance or location.





## QUERY UNIT FROM A MOBILE PHONE

Including

- Alarm Status of AS200TX Transmitters.
- Real Time Values of Input 1, 2 & 3 and Power.
- Contact Telephone Numbers and Alarm Log.
- Network Signal Strength.

## SETUP AND CHANGE SETTINGS FROM A MOBILE PHONE ON AS300RX RECEIVER, ALL PASSWORD PROTECTED.

## Including

- Contact Telephone Numbers, Three Telephone numbers allocated to each AS200TX Transmitter.
- Outgoing Message.

## ACKNOWLEDGE ALARM FROM A MOBILE PHONE

## Including

- Alarm text must be acknowledged from mobile phone, if no acknowledgement is received the AS300RX will repeatedly call every 60 minutes until acknowledged.
- Once acknowledged, recipient's who have received an alarm text will be sent a message giving telephone details of the acknowledger.

## INPUT DATA & ALARM LOG SENT TO WEB SERVER

- Including
- Data logged from the AS200TX transmitter can be sent to the web.
- Alarm log of input alarms and telephone calls made in the event of an alarm. (Note this feature is an option)

Make an informed decision on what action to take based on instant readings from the equipment monitored. This can be sent at any time to your mobile telephone. This is particularly useful if an alarm occurs when the person is on call or out of the workplace. With no need for a dedicated telephone line or line rental the AS300RX is exceptionally easy to install.





www.:www.ritemp.co.uk

## **Features**

The AS300 system is a new generation of monitoring and alarm systems using the Global System for Mobile Communications (GSM) and wireless technology. The monitoring and data logging system allows laboratory and hospital equipment to be protected. The system ensures the safety and continued effectiveness of medicine produce, blood products and samples at specific refrigeration and freezer temperatures. Due to the high value of many of these goods, Quality Assurance programs increasingly require that storage temperatures are to be verified several times per day and that records be maintained. The AS300 System will meet the alarm, monitoring and logging requirements.

Make an informed decision on what action to take based on instant readings from the equipment monitored. This can be sent at any time to your mobile telephone. This is particularly useful if an alarm occurs when the person is on call or out of the workplace. With no need for a dedicated telephone line or line rental the AS300RX is exceptionally easy to install.

## **AS300RX Receiver**

The AS300RX receiver will accept up to ten AS200TX transmitters, more receivers can be added to the system up to a maximum of nine, allowing alarm and monitoring of ninety pieces of equipment on one system. Additional systems can be added if the requirement is greater than ninety pieces of equipment. The AS300RX receivers uses the GSM network for communications, this allows the receivers to be located anywhere a network signal is present. This is particularly useful if customers equipment is distributed on several sites or different parts of the country. The distance between the AS300RX receiver and AS200TX transmitter is 200m, this distance will be greater or less, dependant on the building the equipment is installed in.

Three telephone numbers can be allocated in the AS300RX for each AS200TX transmitter on the system, this allows the person who is responsible for the equipment to be contacted immedi-

- Compact size. 170H x 85W x 35W (mm)
- Simple to use.
- 10 Transmitter channels, one transmitter per channel.
- Additional receivers can be added to expand system up to 90, AS200TX.
- Four menus available giving information on the system.
- Last 10 alarms stored.
- All alarms are time and date stamped.
- Acknowledged by mobile phone.
- Up to three telephone numbers can be allocated per AS200TX transmitter.
- Up to three telephone numbers can be allocated to the receiver.



## **AS200TX Transmitter**

The AS200TX transmitters will in the event of equipment failure, alarm and send data to the AS300RX receiver. The AS200TX transmitter will also display the alarm input and an audible alarm will be activated. The AS300RX receiver will displays the unit number and input of the activated transmitter. The alarm menu in the receiver will log the alarm time and date stamp it. Up to three telephone numbers will be texted giving details of the transmitter alarming and current temperature. The AS300RX receiver can be texted at any time to get details of current values in the AS200TX transmitter.

Up to three inputs are available on the transmitter. A common configuration would be to use input 1 to monitor the air temperature of the equipment, input 2 to monitor the sample temperature and input 3 to monitor the door open / closed condition. This configuration of the inputs will give maximum protection. All data will be recorded if web logging is employed on the system.



- Compact size. 110H x 65W x 27D (mm)
- Simple to use.
- Inputs available : Temperature, and Volt Free Contact.
- Temperature range -200°C to +100°C.
- Temperature transmitter resolution 0.1°C. Temperature transmitter accuracy better than 0.1°C across the full range.
- Front panel display of I/P1, temperature , time delay, unit number ID and alarm parameters.
- Temperature high and low alarms.
- Audible and visual indication of an alarm.
- Adjustable time delay for temperature alarm. (0 to 90 minutes)

Mains power failure alarm.

ately in case of an alarm.

- Alarm information sent to the receiver AS300RX, giving details of transmitter unit and fault.
- Status and alarm indication on transmitter.
  - Temperature, high & low alarms and delay time data transmitted every 15 minutes. This increases to every 10 minutes when in alarm..
- System self test.
- Rechargeable battery backup.
- Easy calibration procedure / Auto calibration reminder customer set. Settings ; None, 6 months or 1 year.
  - Designed bracket for easy installation and positioning.





Edit Reporting Window

BBBB Serial Number: D

Help

## Data Logging & Alarm Log for the AS300RX Receiver

The AS300RX can log data which can be sent to a server, this data can be down loaded to a computer and stored. The data collected can be analysed using our data analysis software. The data interval stored will be every 15 minutes this is the interval that the AS200TX transmits its data to the AS300RX receiver. An alarm log is also kept giving details of all telephone calls to contact and acknowledgement personal, which are time and date stamped.

Datadump Managing your remote data	2	ASPER S	INSTEMS	C		dump Iging your remote data	Logo
	24	d AL hini	···· 1 2 2 m				
Login		Received 2010-12-30	Device IMEI	Ident-Tele No.	Senal	Ham Nessage	Ferriware
Login	23	13:35:24	357820022514434	SNI 0784 2016257	1009	Geo: square ULT	2683U1V014
Username:	13	2010-12-12 12:49:31	357820022514434	SHI 0784 2016257	1009	GSN Atarm and Workloring System	2583U1V014
	田	2010-12-11	357820022514434	SNI 0784 2016257	1000	GSM Alarm and Monitoring System	2683U1V014
Password:	23	2010-12-10 12:40:38	357820022514434	586 0784 2016257	1000	GSM Alarm and Monitoring System	2683U1V014
	15	2010-12-09 12:52:55	357820022514434	BW 0794 2016257	1000	GBM Asarm and Monthoring Bystem	2563019014
Login	10	2010-12-08 12:49:54	357820322514434	SIN 0784 2016257	1002	GBIt Alarm and Wonitoring Byatem	2683011/014
	12	2010-12-07 10:43:01	357820022514434	SM 0784 2016257	1009	GSN Alarm and Workforing Bystem	2683U1V014
	15	2010-12-05 10:43:08	357820022514434	SBN 0784 2016257	1009	GBM Alarm and Monitoring Bystem	2683019014
	21	2010-12-05 10:43:18	357820022514434	SW 0784 2016257	1009	GSM Atarm and Monitoring Bystem	2683U1V014
	10	2010-12-04 10:43:43	357820022514434	SMI 0784 2016257	1009	GBN Alarm and Noreloring System	2683019014
sis - [10/11/2010 to 16/11/2010]	251	elected they and d	bairwas			٥	sentioad Help File Contact II

## Data Dump

Data from the AS300RX will be sent to the server once every 24 hours. The data can also be sent on demand by the user. To access the data a username and password is required. The data is organised into zip files and can be downloaded to the customers computer.

John H	amber. 0		Probe I		Pro	be 2			Trineda: 2	C - See Note: 110 March 1 Marc
First Re	ading: 10/11/2010 13:4	6:42	Max Temp (	°C): 21.6	Max	: Temp (°C): 2	1.7		Deatement 15	
Last Re	ading: 16/11/2010 15/2	E 10/11/2010 13:46:42 E 16/11/2010 15:29:43 ngs: 167 wember 2010 - 13:46:42 E ph. Alama me Probe 1 Actual of the second se	Min Temp ()	C): 11.1	Min	Min Temp ("C): 11.1			Description LS	Land Transfer (2445) 2022 2015 Weilling (2 12 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2
	Readings: 167		Average Ter	mp (°C): 16.1	Ave	arage Temp (°C):	10.0		P film	Description (1)     Descriptio
				10.1		nago romp ( o).	10.0		The: 24 B	ang 2007 🗰 2012 10 201 20 Antonio 2007 👘 2012 12 👘 😨 Bent Antonio 2017 Antonio
Filte									Jate On	C Real
Start:	10 November 2010 -	13:46:42 👙	End:	16 November 2	010 - 15:	29:43 🕂	Filter		-	for (2 hours 2017) [Both ] [10] [Instein] [Instein] [Instein] [Instein]
										1937
Data	1									
Data	Graph Alarms									N) N N N N N N N N N N N N N
	Time	Actual	Probe 1 Upper	Probe 1 Lower	Probe 1 Time	Probe 2 Actual	Probe 2 Upper	Probe 2 Lower		10
	Time		Limit ("C)	Limit ("C)	Delay	Temp ("C)	Limit ("C)	Limit ("C)		
•	10/11/2010 13:46:42	15.2	28.0	-199.0	1	14.9	99.0	-199.0		
	10/11/2010 14:47:05	15.9	28.0	-199.0	1	15.7	99.0	-199.0		
	10/11/2010 15:47:28	15.3	28.0	-199.0	1	15.0	99.0	-199.0		
	10/11/2010 16:47:52	14.8	28.0	-199.0	1	14.6	99.0	-199.0		
	10/11/2010 17:48:15	14.4	28.0	-199.0	1	14.4	99.0	-199.0		
	10/11/2010 18:48:39	13.8	28.0	-199.0	1	13.9	99.0	-199.0		41
	10/11/2010 19:49:02	13.3	28.0	-199.0	1	13.4	99.0	-199.0		
	10/11/2010 20:49:26	12.7	28.0	-199.0	1	12.7	99.0	-199.0		
	10/11/2010 21:49:49	12.5	28.0	-199.0	1	12.4	99.0	-199.0		
	10/11/2010 22:50:13	12.2	28.0	-199.0	1	12.2	99.0	-199.0		
	10/11/2010 23:50:36	12.1	28.0	-199.0	1	12.0	99.0	-199.0		
	11/11/2010 00:50:59	11.9	28.0	-199.0	1	11.9	99.0	-199.0		
	11/11/2010 01:51:23	11.9	28.0	-199.0	1	11.8	99.0	-199.0	10	
	11/11/2010 02:51:46	11.8	28.0	-199.0	1	11.8	99.0	-199.0	10	Data Analysis
	11/11/2010 03:52:16	11.6	28.0	-199.0	1	11.6	99.0	-199.0	10	Data analysis software allows the user to look at the data in numerical
	11/11/2010 04:52:33	11.3	28.0	-199.0	1	11.3	99.0	-199.0	10	
	11/11/2010 05:52:57	11.1	28.0	-199.0	1	11.1	99.0	-199.0	10	and graphical forms. Information at the top of the screen is displayed
	11/11/2010 06:53:20	11.4	28.0	-199.0	1	11.4	99.0	-199.0	10	
	11/11/2010 07:53:44	11.7	28.0	-199.0	1	11.5	99.0	-199.0	10	for the selected unit. Information includes, Serial Number, First Read-
	11/11/2010 08:54:09	12.3	28.0	-199.0	1	12.2	99.0	-199.0	10	ing Last Booding, Number of Boodings, May Tomp, Min Tomp and Av
	11/11/2010 09:54:31	12.5	28.0	-199.0	1	12.3	99.0	-199.0	10	ing, Last Reading, Number of Readings, Max Temp, Min Temp and Av-



12.5

12.

12/

28.0

28.0

199.0

12.6

99.0

-199.0

11/11/2010 09:54:31


11/11/2010 10:54:54

11/11/2010 11:55:17

The data in the columns displays information on Date, Time, Temperature Input 1, Temperature Input 2, High Alarm Set-point, Low Alarm Set-point and Delay Time. Graphs with scale breaks can be displayed by simply clicking on the graph tab.

erage Temp for I/P1 and I/P2.

The alarm log gives historical information of time, date, call direction, telephone numbers, alarm input and whether the call was successful or failed.



## Reporting

Print out of the data can be done in three ways.

- Transmitter Report prints the data ٠ column view.
- Graph Report prints the graph view. ٠
- Summary Report prints the "header" information for each transmitter and alarm log.



## CALIBRATION

A calibration service can be offered for all our products.

## SERVICE AND SUPPORT

We service and repair all of the electronic products we manufacture. All our goods carry a 1 year return to base warranty. If one of our products needs repair while under warranty, we will repair or replace it

### ALARM & MONITORING APPLICATIONS

- ♦ ULT FREEZERS
- ♦ LN2 STORAGE
- ♦ INCUBATORS
- ♦ ISOLATORS
- CLEAN ROOMS
- FREEZERS
- FRIDGES
- COLD ROOMS
- CULTURE ROOMS
- CLIMATE MONITORING
- CRYOGENICS
- COMPOST
- TEMPERATURE, -200°C to +100°C

## System Specification and Functions

## AS300RX RECEIVER

## Number of AS200TX Transmitters per AS300RX Receiver

10 AS200TX,
System expandable to 90 AS200TX

## Built in Communications

- 3 Telephone numbers allocated for each AS200TX transmitter.
- ♦ 3 Telephone numbers allocated per receiver.
- Web access for data and alarm history (Option)

#### **GSM Signal**

 Visual indication of signal strength. Audible and visual indication if network strength is lost or low.

#### Memory

Over 4000 readings

## Battery Back Up

Rechargeable Batteries

## **External Dimensions**

♦ 160H x 90L x 50D mm

## FUNCTIONS INITIATED FROM MOBILE PHONE

- Enable or Disable AS300RX.
- Setup Telephone Numbers, mobiles and land lines.
- Enquire Details from AS200TX transmitters Including Alarm Status, Temperature Readings, Alarm Log, Transmitter Signal Strength, GSM Signal Strength, Setup Values, Telephone Numbers and OGM Programmed, Amount of Credit in Sim Card.
- Acknowledge Alarms.
- Test Alarm Function.

## AS200TX TRANSMITTER

#### Inputs

- Temperature and Volt Free.
- Temperature Range +100°C to -200°C.
- Temperature resolution 0.1°C.
- Temperature Accuracy 0.1°C.

## Display

- Transmitter Unit Number
- Time Delay
- Temperature Input 1

## Battery Back Up

Rechargeable Batteries

- Alarms
- High and Low Alarms.
- Adjustable Time Delay For Temperature Alarms (0 to 90 mins).
- Power Fail Alarm.
- Status & Alarm information on transmitter.

## Data

 All data transmitted including temperature, high & low alarms and delay time every 15 minutes.

## External Dimension

110H x 65W x 27D mm





www.:www.ritemp.co.uk

contact@ritemp.co.uk

- eiver Passwords
  - Master and Acknowledger 4 digit password.
  - Lock switch on receiver, disabling change of settings by mobile phone.

## **Telephone Alarm Messages**

Up to three telephone numbers can be programmed into the AS300RX for each AS200TX on the system. A text message will be sent to a mobile or a voice message sent to a land line. An additional three telephone numbers can be programmed into the AS300RX for power fail or test alarms. The outgoing message (OGM) can be up to 80 characters long and the value and input alarming is tagged on to the end of the OGM.

#### Sim Card

 Low Credit Warning text message and visual, audible warning on GSM Unit. (Sim Credit Warning only available on O2 and Vodafone)

## Display

- Scrolling display of Sim Number, Sim Credit, Signal Strength. (Sim Credit only available on O2 and Vodafone)
- Four Menu's :-, System, Alarm Log, AS200TX Temperature Readings and Transmitter Information.