



Instructions for Use

AFO LIGHT

		3
EN	INSTRUCTIONS FOR USE	5
DE	GEBRAUCHSANWEISUNG	8
FR	NOTICE D'UTILISATION	11
ES	INSTRUCCIONES PARA EL USO	14
IT	ISTRUZIONI PER L'USO	17
DA	BRUGSANVISNING	20
SV	BRUKSANVISNING	23
NL	GEBRUIKSAANWIJZING	26
PT	INSTRUÇÕES DE UTILIZAÇÃO	29
JA	取扱説明書	32

ENGLISH

- a. Calf piece and velcro strap
- b. Elastic connection strap
- c. Posterior strut
- d. Thermoformable – tube
- e. Tapered footplate

FITTING THE AFO TO THE PATIENT

- Apply the J-shape spring to the patients' lower leg (**Figure 1**).
- Apply the calf piece (**Figure 2**).
- Verify height of the calf piece (**Figure 3**). To reduce the length of the posterior strut, remove calf piece and grind strut to an appropriate level (approximately at the thickest part of the calf).
Note: Do not grind the strut at the heel area as this may jeopardize product performance and strength.

WARNING

Make sure the strut doesn't reach the back of the knee in seated or slightly kneeled position and does not compress the perineal artery or nerve.

ASSEMBLY OF THE AFO

- Slide the calf-piece over the posterior strut.
- Apply slight tension on the elastic connection strap (**Figure 4**) and attach the PSA Hook.
- Connect the elastic connection strap using the Loop.
- Apply the tapered part (rounded edge) of the posterior strut pad to the distal part of the posterior strut and the entire length of the posterior strut covering the elastic connection strap (**figure 4**).
- Attach the calf pad to the calf-piece (**Figure 5**).
Note: The calf padding is slightly larger than the calf-piece to avoid contact between the piece and the leg.
- Insert the AFO into the shoe and allow the patient to walk around for a while to verify the fit of the AFO. Decrease or increase padding as needed.

Apply the shrink tube to the lower section of the strut to allow the user to don their shoe more easily while protecting the padding.

- Remove the upper part of the padding of posterior strut.
- Remove the elastic connection strap from the PSA Hook on the posterior strut.
- Remove the calf piece from the posterior strut.
- Slide the shrink tube over the posterior strut until it covers the heel area.
- Use a heat-gun to shrink the tube to proper size. Evenly distribute the heat to avoid burning and wrinkles (**Figure 7**). First apply the heat to the inside of the AFO, then to the outside of the AFO.
- Re-apply the calf piece, elastic connection strap and padding.

INSERTION INTO THE SHOE (**Figure 6**)

- If present, remove the insole from the shoe.
- Place the AFO inside the shoe.
- If necessary adjust the shape of the footplate (by grinding) to the contours of the shoe-sole (Use the insole as a reference for shaping).

Note: Do not grind the strut at the heel area as this may jeopardize product performance and strength.

- As applicable, re-insert the insole on top of the AFO-footplate.

DONNING THE AFO

- When donning the product, use the posterior strut of the AFO as a shoe-horn (**Figure 8**).
- Secure strap around the calf (**Figure 9**).
- If necessary, remove excessive strap material with a pair of scissors.

CUSTOMIZING THE AFO

- In case of point pressure at the heel-bone, it is advised to apply padding under the existing padding just proximal of the pressure point (Distal calf). Remove the padding and glue additional padding to the strut before re-attaching the padding.
- For patients with lower activity levels, where the calf piece doesn't move, it is advised to permanently fixate the calf-piece by application of strong glue (Loctite 495.) into the calf-piece before attaching it to the strut (**Figure 10**). The calf piece can also be riveted to the posterior strut. The elastic band and the micro-Velcro may be removed in this case.

GENERAL REMARKS AND ADDITIONAL ADVICE

1. To shorten the footplate, saw or grind off the appropriate length. Be certain to buff the edges to eliminate sharp edges. To lengthen the footplate, use foot insert of choice, and secure it to the foot plate.
2. To adjust the biomechanics of the foot, the use of a foot orthotic is advised.
3. The AFO cannot be reshaped by heating. Heating will cause delamination and negatively alter the dynamics of the orthosis.
4. The right footwear is critically important for proper function of the AFO. The following criteria should be met when selecting the right footwear to use with this AFO:
 - Firm heel counter – for proper control of the rear foot
 - Laced or Hook and Loop fastening
 - Removable insole – to allow space for the AFO footplate with modifications
 - Closed heel and toe box

AFO LIGHT INDICATIONS:

- Mild to moderate drop-foot secondary to neurological conditions

AFO LIGHT CONTRA-INDICATIONS:

- Excessive ankle instability or excessive deformities
- Severe spasticity
- Extreme activities
- Fluctuation Edema
- Foot/Lower Limb Neuropathy

WARNING.

1. The AFO Light is a supportive brace only and is not intended to prevent injury.
2. If you experience any unusual discomfort while using this product, consult your health care provider immediately.
3. Patients with foot neuropathy should not use this product.

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