


High Pressure Cryo-Cooler for X-Ray Crystallography (HPC-201)




ADC is offering a High Pressure Cryo-Cooler for preparing protein crystals. This device is based on a process developed by Cornell University scientists Prof. Sol M. Gruner and Dr. Chae Un Kim. This exciting new technology enables the simultaneous capture of both amplitude and phase information from single anomalous diffraction (SAD) of a cryocooled protein crystal, thereby providing sufficient data to solve the crystal structure of a protein with an unknown structure. Flash-freezing at atmospheric pressure requires the use of cryoprotectants. Finding the right cryoprotectant for a sample type can be a long trial-and-error process. The High Pressure Cryo-Cooler eliminates the need to

use cryoprotectants and produces superior results.

The scatter images below of a glucose isomerase crystal prepared at atmospheric pressure (left) and under high pressure (right) demonstrate the benefits of high-pressure cryocooling.



Resolution = 5.0 Å
Mosaicity = N/A



Resol. = 1.1 Å (1.3 Å for 3 crystals)
Mos. = 0.39° (0.48° for 3 crystals)

The high-pressure cryo-cooler is designed to hold 3 samples at a time. Crystal samples are picked up using a standard cryoloop. Cryoloops are mounted to heavy duty stainless steel tubing in the unit and are then ready to be pressurized and cooled. A high pressure oil pump provides helium gas to the samples. External controls allow the sample to be first pressurized and then cooled by a LN₂ bath. Once pressure is released the samples can be removed and handled like any other samples prepared by the conventional flash freezing.



126 Ridge Road, Lansing, NY, 14882
Tel: (607) 533-3531 • Fax: (607) 533-3618
adc@adc9001.com • www.adc9001.com

Features

Pressurizing Gas

Working Pressure

Cooling Fluid

Cryo Cooling Temp

Sample Capacity (per pressure & cooling cycle)

Process Time

ZEISS Microscope

LN₂ dewar

Helium

200 MPa

LN₂

77 K (-196 °C)

3

< 10 min (2 min for pump operation; ~ 5 min under pressure; 1min freezing)

StEREO Discovery.V8

Taylor-Wharton HC34

Connection Data

Voltages, frequencies

Power input, approx.

Oil pump pressure

LN₂ Dewar Holding

50/60 Hz

115/230 VAC

200 MPa

200 Days

