



Patient Fitting Instructions

- *Take your EZ Stride-ATB orthotic and slide it into your shoe*
- *Undo both straps*
- *Pull the anterior (front) section forward towards your shoe and “open” the EZ Stride up*
- *Slide your foot into your shoe*
- *As you slide your foot into the shoe, simultaneously slide your leg and calf into the posterior (back) section of the EZ Stride-ATB.*
- *Close the anterior (front) section by pulling it up towards your leg. Do this until the anterior section contacts your leg.*
- *The top of the anterior (front) section should slide into the posterior (back) section. There should be some overlap of the carbon fiber “wings” or extensions.*
- *Take each strap and thread it through it’s respective loop located on the anterior (front) section.*
- *Pull each strap through the loop, then out towards the side, and finally back towards the back of your leg. The EZ Stride-ATB should fit snugly.*
- *Lay each strap on itself after you have firmly pulled the strap tight. The hook in the end of the strap will secure to the loop of the strap.*
- *You are now ready to walk in your new EZ Stride-ATB!*