



GO!

ONE COMPONENT SELF-ETCHING
LIGHT CURED ADHESIVE



ONE STEP, ONE BOTTLE, SELF ETCHING ADHESIVE

A SINGLE COMPONENT, LIGHT CURED SELF-ETCHING, FLUORIDE RELEASING DENTINE AND ENAMEL BONDING SYSTEM

Go! is a super quick adhesive saving you valuable chair time: you don't have to shake the bottle, etch the dentine nor do you have to wash it off. You don't even need to rub or agitate it! Just one application is required. Once cured, you can immediately place your restorative over Go!. Your patient is then good to Go!. Try Go! on your pediatric patients because it is so fast to use.

BPA &
HEMA
FREE

7%


FILLED WITH
HYDROPHILIC AND
HYDROPHOBIC FILLERS



LOW FILM THICKNESS

Go!'s low film thickness provides an invisible margin between the restoration and the tooth.


FILM THICKNESS*



HIGH BOND STRENGTH

Go!'s high bond strength ensures the restoration will stay in place for many years to come.

MICROSHEAR BOND STRENGTH [MPa @ 24HRS]**



ULTRAFINE NANO FILLERS


ULTRAFINE NANO FILLERS

Go! is 7% filled with specially treated hydrophilic and hydrophobic fillers that change from a water loving to a water repelling state after polymerization. This provides a surface with no water residue that is perfect for bonding to the restorative.

SUPER HIGH FLUORIDE RELEASE

Go!'s high fluoride release assists in strengthening the tooth and preventing recurrent decay.

CUMULATIVE FLUORIDE RELEASE^{*}




QUICK COLOUR CHANGE

Initially, Go! is a dark yellow colour, but upon curing, it changes to clear. You know when you can go to the next step.


NO WATER RESIDUE

Unlike some 7th generation adhesives, there are no water droplets found at the Go!/ restorative interface to interfere with the adhesion.



Go! surface[^]

Microscopic image of Go! after polymerization.



G-Bond[^] surface[^]

Microscopic image of G-Bond after polymerization.

[^] Images courtesy of Dr G Knight, University of Adelaide

ADVANTAGES

No Bis-Phenol A or its derivatives

Etch, prime and bond in only 35 seconds

Quick - only one application, no scrubbing

Time efficient - no shaking of single dose or bottle required

High bond strength to cut enamel and dentine

Changes colour when completely cured

No more post-operative sensitivity

No water residue

Fluoride releasing

HEMA free

INDICATIONS

Direct filling restorations using light-cured composite and compomers

Composite repair of porcelain and hybrid ceramic restorations

Sealing exposed root surfaces causing hypersensitivity


Cavity sealing prior to placement of indirect restorations

Core build-ups using light cure or dual cure composites

Composite repair

INSTRUCTIONS


1 Clean and isolate tooth.



2 Remove excess water. Preferably keep moist, but a dry surface is ok.



3 Apply Go! to saturate all internal surfaces and leave in place for 20 seconds.




4 Blow firmly with dry, oil free air for 5 seconds. Leave surface glossy.




5 Light cure for 10 seconds.



6 Apply composite or compomer such as Aura, Aura Easy or Aura Bulk Fill according to manufacturer's instructions.



ORDER DETAILS



KITS

Go! 5mL Bottle Kit	8100400
1 x Go! 5mL bottle	
50 x Points disposable brush applicators	
1 x Mixing well	
Go! Single Dose Kit	8100401
50 x Go! 0.1mL Single Dose	
50 x Points disposable brush applicators	
1 x Single dose holder	

*SDI Test Data

**University of Melbourne

~G-Bond is not the registered trademark of SDI Limited.



YOUR SMILE. OUR VISION.

CE
0297

12-2021

M810400 G



+D036M8104001.

MADE IN AUSTRALIA
by SDI Limited
Bayswater, Victoria 3153
Australia
www.sdi.com.au

AUSTRALIA 1800 337 003
AUSTRIA 00800 0225 5734
BRAZIL 0800 770 1735
FRANCE 00800 0225 5734
GERMANY 0800 100 5759

ITALY 00800 0225 5734
NEW ZEALAND 0800 808 855
SPAIN 00800 0225 5734
UNITED KINGDOM 00800 0225 5734
USA & CANADA 1 800 228 5166