

Vision Tester VT1

Automated screening and visual function exploration device



Depisteo

Ideal for: * Occupational Health practitioners.

* FAA examiners

* School Nurses and pediatricians



Far and near vision

Without constraint or effort, the patient moves from far vision to near vision simply by lowering the regard.

The reading angle, through lower lenses, is adapted to wearers of progressive lenses, who are used to reading in this manner.







Revolution in Visual Exploration

Remotly controlled system



Remote control

No physical interaction with the device needed. The programmable remote control of the Depisteo Vision Tester VT1 is capable of initializing a quick series of routine tests or a full, variable series of in-depth tests.

Tests results are logged on a Patient Response Form.



Computer Interface

The Computer Interface System allows the use of a regular PC to run the Depisteo Vision Tester VT1. Driven by sophisticated yet user-friendly software, test results are viewed on your computer screen and printed with a standard printer.

The Depisteo Vision Tester VT1 is connected to the computer via the supplied USB lead.



Wide Range of Visual Screening Tests

A multitude of tests in a light-weight device!

The test band is composed of high-density photographic film, containing 16 different slides.

Looking into the Vision Tester, a sensation of quality is perceived, with the fineness of the optotype contours and the high density of the black in relation to the white.

Main tests*:

Test 1	Letters (Binocular)	
Test 2	Letters (Monocular)	
Test 3	Tumbling E (Monocular)	
Test 4	Tumbling E (Binocular)	
Test 5	Astigmatism Test (Binocular and Monocular)	
Test 6	Landolt ring (Binocular)	
Test 7	Landolt ring (Monocular)	
Test 8	Low Vision Test (Binocular)	
Test 9	Low Vision Test (Monocular)	
Test 10	Stereopsis Test (Binocular)	
Test 11	Phorias (Vertical and horizontal)	
Test 12	Fusion (Binocular)	
Test 13	Numbers colors test	
	(Binocular or Monocular)	
Test 14	Daisy colors test	
	(Binocular or Monocular)	
Test 15	ARMD Test	
	(Binocular or Monocular)	
Test 16	Contrast Test	
	(Binocular or Monocular)	


Silver Model:

With this model, all the above tests can be performed. The level of lighting can be lowered to perform visual acuity tests in "dusk" lighting. Tests are performed from as far as 20', or as close as 14".

Master Model:

The Master version includes the tests from the Silver model as well as the tests below:

Test 17** Intermediate vision test (24 inches)
Test 18 Vertical and horizontal visual field test



2

3

8

1			
1 0 * \(\text{\(\text{\) \exiting \eta}\}}}}}\end{\(\text{\(\text{\(\text{\) \eta}\}}}} \end{\(\text{\(\text{\(\text{\(\text{\(\text{\) \eta}\}}}}}} \end{\) \\ \) \\ \} \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\	2	3 □ O ★ Δ 6 □ O ★ Δ	1
1 7	2 5 8	3 6 9	1
	•	•	1
4	5	3 100 6 100 100 6 200 100 200 100 100 100 100 100 100 100 100	1
1 0 0 0 0 4 0 0 0 0	2 0 0 0 0 0 0 0	3 0 0 0 0 0 0 0	1
	·		1
1 2 3			1

^{*} The Band and the test list are an illustration of the tests available, customized tests are available on demand.

^{**} Other lenses are available on demand.

Visual Field Test



Horizontal and vertical visual field screening.

20 LED stimuli sweep the horizontal (temporal & nasal) and vertical visual fields.

Removable Headrest detection

To ensure the patient's face is in the correct position, the slides will only light up when the forehead is in the correct position, activating the pressure sensitive elastometer.

This part is easily removed for disinfection or interchanging.



Inclination adaptable to the patient

The Depisteo Vision Tester VT1 is inclined for ease of transport and for easy access to the connection leads.

It also adapts to each patient, allowing them to remain in a relaxed position throughout the test.









Features:

Depisteo LCC reserves the right Depisteo Vision Tester VT1.

Standard Accessories

- 1 user manual
- 1 power plug
- 1 headrest

Remote-control Version

- 1 remote control
- 1 pad of 100 patients response forms
- 1 table of correct responses

Computerized version

- 1 Depisteo software
- 1 USB lead

Optional Accessories

Transport case Disinfecting liquid or cloths Pads of 100 patients

The VT1 is automated and remotely operated. It can thus be placed in the ideal position in the test room, without the operator having to move around.

Amazingly quick and silent, the VT1 is operated with the fingertips, or with the click of a mouse. Tests can be adapted to the professional risks of each patient, by programming test protocols.

The lighting arrangement is composed of 8 white LEDs on the right side and 8 white LEDs on the left side, driven by a micro-controller. This lighting is homogenized by a diffuser, ensuring clear, powerful back-lighting.

The LEDS are low consumption, avoiding high voltage within the VT1. The absence of heat increases the life-span of the product and the display slides.

General features:

Silver Model: Far and near vision lenses

Master Model: Far, near and intermediate vision lenses,

vertical and horizontal visual fields

Remote control

Back-lit LCD remote control Driver Option:

Depisteo Software

VT1 software and USB lead. Driver Option:

In working position, maximum height of 19.7 " Height:

Width: Maximum width at base: 10.2 '

Without driver system: 10.3 lbs (without remote control) Weight: Luminous intensity: Maximum: 200 candelas/m2. Progressive increase in intensity.

Lighting by 16 white LEDs and diffuser system. Lighting system: Reactivity: Time between 2 neighbouring tests: 2 sec.

Time to go from one lens to another: 1.5 sec. Test time: Routine test: 2 min. In-depth test: 4 min. Visual acuity tests: Optotypes: letters, numbers, Landolt rings. Other tests: Contrast sensitivity, astigmatism, phorias,

stereopsis, fusion, ARMD, colors

(vertical and horizontal field option, Master model)

Luminous conditions: Day vision, low light vision, dusk vision.

Conformity:

Standards Test specifications defined according to data

of the EN ISO 8596 ophthalmic optics standard: measurement of visual acuity,

standardized optotypes and their presentation.

CE Marking: Class I medical devices

according to directive 93/42/EEC.

EN60601-1 standard. Safety:

Electromagnetic compatibility:

IEC 60601-1-2 standard.



3495 Piedmont Rd NE Building 11 suite 710 ATLANTA, GA 30305 **USA**



Distributed by:

www.depisteo.com

Tel (+1) 404-495-59-11

Fax (+1) 404-233-86-25

contact@depisteo.com