

Fitting guide for Pixi



- Hold the mask gently against the child's face, checking that the nose is fully enclosed by the cushion and the mouth is clear.



- Pull the headgear over the head and move the middle strap down until it sits across the middle of the back of the head.



- Keeping your hand on the frame to maintain correct headgear position, bring the lower headgear strap behind the back of the neck.



- Attach the lower headgear strap to the mask by holding the tab and hooking the quick-release clip onto the lower headgear arm.



- If necessary, adjust the straps by removing the tabs and gently tightening the headgear straps until comfortable. Do not overtighten.
- For best results, adjust the straps starting at the bottom strap and work upwards.




- Connect one end of the air tubing to the swivel and the other end to the therapy device.
- Turn on the therapy device.

Pixi tips



- When fully fitted, ensure the top strap is vertically aligned with the top of the ears, and not overtightened. The middle strap should sit across the back of the head, horizontally aligned with the top of the ears.



- If fitting the mask for the first time, we recommend the child be sitting in an upright position to ensure correct fit.

Note: For complete user instructions, refer to the Pixi Paediatric User guide.



Designed with care

ResMed has partnered with paediatric health experts to bring you Pixi™: a paediatric mask designed with children and their families in mind.

Engineered for little faces

Most masks used for children are small adult-sized masks—but not the Pixi mask. Specially designed using the head and facial characteristics of children aged between two and seven, Pixi takes into account the facial contours of a child's face.

Comfort, reassurance and acceptance

The headgear is designed to sit away from the child's eyes and ears to help reduce irritation and obstruction. Not only can the child see clearly while wearing the mask, but carers are also reassured and comforted by being able to see their child's eyes and more of their little face.

Welcome to Pixi's world

To help children overcome the psychological barriers to therapy, we've created a parent guide with tips from clinical specialists to support parents with their child's journey, a picture book, and colouring sheets of the Pixi character.

ResMed's custom paediatric solution

When you combine Pixi with ResMed's latest Astral™ or Stellar™ ventilators, or our VPAP™ S, VPAP ST or VPAP ST-A pressure support devices, you have an ideal paediatric solution designed to deliver more comfort and reliability.*



Pixi features and benefits

Soft, comfortable cushion
has a single-wall, thin silicone design that minimises pressure on children's soft faces, while providing a reliable seal

Comfortable Breath-O-Prene™ headgear materials
are lightly coloured for a more reassuring, unobtrusive appearance on the child's face

Side-facing back strap
ergonomically designed from the fitting adult's point of view

Two tubing positions
allow air tubing to be moved to either side of the cushion depending on your child's sleeping position

Flexible stabilisers
help the mask conform to each child's unique facial structure and shape for the most customised fit

Quick-release latch
makes it easy to take the mask off, while retaining desired headgear settings

Spring-flex lightweight tube
stretches and bends easily, so it reduces the chance of getting caught or dragging as your child moves around during sleep

In an emergency, release the latch even faster by pulling firmly on the lower headgear strap

Pixi components

Plug (attached to frame)

Tube assembly (elbow and tube)

Frame **61035**

Cushion **61033**

Headgear **61034**

Combine the Pixi mask with the latest Astral, Stellar or VPAP devices for a paediatric therapy system that delivers comfort, reliability and effective care.



**Stellar 100 not sold in ANZ

Product code

	ANZ/APAC
Pixi paediatric mask	61030

*Please refer to the relevant device User guide for indications for use with the mask.

Pixi™
Paediatric mask