

RehabPro E-Series Pulleys M70E Single & Double Wall Mount

User Manual



REHAB-PRO
www.rehabpropulleys.com

1. Inspection of Delivery

Before unpacking the Pulley please inspect the packaging for shipping damage. Observable damage must be noted **prior** to signing the bill of lading. A claim **cannot** be file unless damage is noted. Please contact RehabPro, Inc. or your local dealer with any concerns or questions about shipping damage.

2. Shipping, Content & Assembly

The RehabPro E-Series Pulleys (wall mounted) are shipped fully assembled in cardboard boxes bolted to wooden skids (below). M70E (single wall) is single stacked. M70E (double wall) is double stacked. Pulley accessories, if ordered with the pulley, are normally included inside the box. Please refer to the packing slip for content. Contact RehabPro or your local dealer immediately if ordered items are missing / content differs from attached packing slip.



The lower Pulley is bolted to the wood skid



**Tools needed to open
box & undo pulley**

There is no pulley assembly required.

3. Installation of Wall Mounts

Wall mounting pulleys require the following knowledge:


1. The base of the pulley must rest on the ground. Due to their design, they CANNOT be suspended on the wall.
2. Make sure your mounting wall reaches the ceiling. A wall ending at the false ceiling does not offer enough support to secure a pulley.
3. Identify the type of studs you have in your wall. Do you have wood studs or metal studs, or is your wall made of concrete? Different types of studs require different hardware to attach your pulley equipment safely.

RehabPro wall pulley attachment brackets:


RehabPro standard pulleys have two (2) attachment brackets; the center hole of the lower attachment bracket is 3.5" above the floor; the center hole of the upper attachment bracket is 80.5" above the floor. Each attachment bracket is 18" wide with two (2) holes for mounting hardware, 16.25" apart center to center.

Below, you will find three attachment suggestions based on the type of wall you plan on attaching your pulley to:

1. Drywall attached to wood studs. Attach two 1x4 boards horizontally in space, long enough to span two studs (4 studs for parallel pulleys), 2.5" lag bolts that go through the board, drywall and into the studs (see picture below). The center of the lower board must be 3.5" off the floor. The center of the upper board must be 80.5" off the floor. The pulley(s) can be anchored wherever along these boards using 1.5" lag bolts with a 5/16" flat washer.



2. Drywall attached to aluminum studs. Attach two 1x4 boards horizontally in space, long enough to span two studs (4 studs for parallel pulleys). Use 2.5” anchor bolts (see picture below) to sandwich the boards against the wall. Like attachment option # 1 the center of the lower board has to be 3.5” above the floor. The center of the upper board has to be 80.5” above the floor. The pulley(s) can be attached to these boards by 1.5” lag bolts with a 5/16” flat washer.



Drill your holes large enough for the anchors to pass through.

C. Pulleys attached directly to a concrete wall. Attach the pulley directly onto the concrete wall by placing correct hardware (concrete anchor bolts) at the two points of attachment for each attachment bracket 3.5” & 80.5” above the floor.

4. Utilization

1. Adjusting the weight

The RehabPro E-Series Pulley weight stacks are “magnetic, pin select”. The resistance is measured in pounds. The weight stack measures as follows:

M70E

Increment	Weight	Resistance	Increment	Weight	Total Weight
# 1	2.50 lbs.	2.50 lbs.	# 11	2.50 lbs.	27.50 lbs.
# 2	2.50 lbs.	5.00 lbs.	# 12	2.50 lbs.	30.00 lbs.
# 3	2.50 lbs.	7.50 lbs.	# 13	5.00 lbs.	35.00 lbs.
# 4	2.50 lbs.	10.00 lbs.	# 14	5.00 lbs.	40.00 lbs.
# 5	2.50 lbs.	12.50 lbs.	# 15	5.00 lbs.	45.00 lbs.
# 6	2.50 lbs.	15.00 lbs.	# 16	5.00 lbs.	50.00 lbs.
# 7	2.50 lbs.	17.50 lbs.	# 17	5.00 lbs.	55.00 lbs.
# 8	2.50 lbs.	20.00 lbs.	# 18	5.00 lbs.	60.00 lbs.
# 9	2.50 lbs.	22.50 lbs.	# 19	5.00 lbs.	65.00 lbs.
# 10	2.50 lbs.	25.00 lbs.	# 20	5.00 lbs.	70.00 lbs.

The RehabPro E-Series Pulleys (wall mount) come with 1x Pulley Handle per weight stack as standard accessories.


2. Adjusting the Angle of the Rope (based on the length – tension concept)

The upper rope adjustment bracket changes the angle of pull, allowing each exercise to be performed correctly from an anatomical and physiological point of view. The rule of thumb is as follows:

A – Check available range of motion.

B – Chose which exercise to be performed.

C – Pulley rope should be at 90 degrees with the lever arm at mid-range (mid-range may change as motion around the joint(s) improve).



Picture shows the upper rope adjustment bracket correctly attached to the rope adjustment bar. “Pull out” on the black spring-loaded adjustment pin in order to move the rope adjustment bracket up or down the bar.

Make sure the spring-loaded pin pops back into place (all the way) before attempting to exercise using the pulley.

3. Adjusting the Length of the Rope

The lower rope adjustment bracket changes the range of motion of the exercise that is being performed. The rule of thumb is as follows:

A – You may want to reduce range if the “initial” end range becomes painful during an exercise. Move the rope adjustment bracket “up” towards the upper adjustment bracket to allow more “slack”.

B – In the initial stages of rehab you may want to give the muscle group(s) a break between each repetition. Adjust the rope adjustment bracket up allowing the weights to hit the stack at the completion of each repetition before starting a new rep.



Picture shows the lower rope adjustment bracket correctly attached to the rope adjustment bar. “Pull out” on the black spring-loaded adjustment pin to move the rope adjustment bracket up or down the bar. Pull and rotate the spring-loaded pin head to allow the lower rope adjustment bracket to move together with the upper rope adjustment bracket.

Make sure you counter-rotate the spring-loaded pin head so it pops back into place (all the way) before attempting to exercise using the pulley.


5. Maintenance

Check the condition of the pulley daily prior to use:

1. Pay attention to any “loose” nuts or bolts. Re-tighten nuts and bolts after the pulley has “settled”, normally within a week or two after the pulley has been put into use. Check the pulley on a monthly basis after initial check.
2. Pay attention to the rope. On the average, the pulley rope lasts for 12 months with normal pulley use. Make sure the rope glides smoothly on the pulley wheels. AND, look for “fraying”. Change the rope when the outer layer has become “fuzzy”. You can contact RehabPro, Inc. directly or your local dealer for a rope replacement kit.

Rope Replacement Kit: Item # PR16316, 3/16” non-stretchable marine rope.

How to thread / replace the pulley rope:



1. Thread the rope through the hole in the center rod (from below – up) Make a simple knot to ensure that it cannot “pull out / up” through the hole).
2. Bring the rope up to the pulley wheel under the top box. Thread it around / over this wheel back to front.
3. Bring the rope all the way back down to the wheel on the lower glider. Thread the rope around this wheel back to front.
4. Bring the rope up to the upper glider. Thread the rope between the two pulley wheels, back to front.

Tying off the Rope at the hook:

- A. Pull the rope through the eyelet making a “U” & loop it around itself as shown in picture A.
- B. Loop the rope around itself, front to back 2 - 3 times & pull the rope through “the oval” created by the rope (back to front) as shown in picture B)
- C. Tighten the rope on itself by pulling the knot down to the eyelet as shown in picture C.



Picture A



Picture B



Picture C

D. Pull the shrink rubber over the knot and shrink it by using a heat gun or a hair blower as shown in pictures D & E.



Picture D



Picture E

3. Wipe off the guide rods with a clean rag weekly. After the wipe down spray **silicone lubricant** (NOT Silicone Grease) on a rag and wipe down the rods with the lubricant.

You may want to wipe down the rope adjustment bar also for a smoother glide of the rope adjustment assemblies.

RehabPro recommends Blaster silicone lubricant. It can be purchased at Home Depot.



6. Technical Information

Total weight:	150.00 lbs.
Width:	18"
Depth:	8"
Height:	80"
Weight Stack:	70.00 lbs.

7. Warranty

- RehabPro, Inc. offers a lifetime warranty on pulley frame.
- RehabPro, Inc. offers a 3- year warranty on all moving parts.
- RehabPro, Inc. offers a 30-day warranty on pulley rope, snap hooks, and pop-pins
(Items not made by RehabPro, Inc.)

(Cost of shipping replacement items are not covered under warranty)

8. Accessories

Item #	Accessory Description
R282H834	Standard Pulley Handle
2020	Padded Wrist / Ankle Strap
2021	Padded Shoulder / Knee Strap
2022	40" Padded Waist Strap
2024	Padded Multi-Purpose Strap
2025	Banana Sling
13107	Standard Triceps Bar

RehabPro, Inc.
28215 North 25th Dale
Phoenix, AZ 85085
Phone: 602-317-7207
E-mail: sales@rehabpropulleys.com
E-mail: kris@rehabpropulleys.com

