


Handicare 1000 Perch


Straight Stairlifts


Handicare 1000 Perch

Straight Stairlifts









Perch seat		Dimensions
Α	Minimum swivel radius*	-
В	Wall to inner seat back	2.95in
С	Seat depth	12.59in
D	Minimum folded width to footplate	15.15in
Е	Minimum open width to edge footplate	25.59in
F	Armrest width – external	21.65in
G	Armrest width – internal	18.30in
Н	Seat back height	17.91in
I	Footplate to seat height	26.57-34.44in
J	Minimum footplate height	2.75in
Minimum track intrusion into staircase		4.52in
Maximum weight limit		250lbs
Speed		0.07m/sec
Staircase angle		27 - 45°

^{*}Fixed seat

Some of the above dimensions depend on the characteristics of the staircase in the home. Please consult the sales manual for further information about these dimensions and features of this seat. All information is provided in good faith but does not constitute a guarantee. Handicare reserves the right to alter designs and specifications without prior notice.

HANDICARE 10888 Metro Court St. Louis, MO 63043 866.891.6502 handicareusa.com