



LIFE CARE HEARING CENTRE

2190 Warden Avenue (at Sheppard)
Scarborough, Ontario M1T 1V6

Tel.: (416) 496-3200
Fax: (416) 496-3833
lifecarehearingcentre@gmail.com

TIPS

Please clean hearing aid weekly.

RIGHT HEARING AID - **RED**

LEFT HEARING AID - **BLUE**

For behind the ear hearing aids



Click the button on the left ear to decrease the volume of both aids.

Click the button on the right ear to increase the volume of both aids

Hearing Aid Not Working?

- Change the battery when you do not hear anything
- Battery Flat side up (+/positive) on the battery door
- Look for wax buildup on the hearing aid and clean with the pick/brush
- Book an appointment at Life Care Hearing Centre if more help is needed

Call **(416)496-3200** for help, or to book an Appointment