

VORTEX SHAKER

Features:

General laboratory Vortex mixer.

Can be used for tube mixing and flask mixing by replacing head.

Push touches type with manual/continuous select switch.



Specifications:

Model		RIC-VM-B
Overall size(WXDXH)	mm	125X185X160
Speed range	rpm	3,000 Max.
Motor W		23
Material		ABS (acrylonitrile-butadiene-styrene)
Head		Rubber mixing cab (test tube Ø35mm & flask 250ml (max.))
Power source & consumption		AC 230V 50/60Hz 1phase
Net weight	kg	2.6



RELIABLE INSTRUMENTS CO. (ISO 9001 CO.)

Sales/Reg. Office: 301, 2nd Floor, Mittal Commercial Complex, N.I.T-5, Faridabad
Tel.: 91-129-4052611, 2410244, 6450611; Telefax: 91-129-2410219;
M: 91-9818047375
Email: ricindfbd@gmail.com; ricind@rediffmail.com; ricind@yahoo.co.in;
Web: www.ricind.com

