

PAIN RELIEF COLD LASER DEVICE

SEMICONDUCTOR
LASER
TREATMENT INSTRUMENT

GD-P

Application Range:

Cervical spondylosis, Slipped disc,
Arthritis, Shoulder of peri-arthritis,
Prostatitis, Bronchitis, Pelvic infection,
Facial paralysis, Back pain .



SPECIFICATION



ITEM	DESCRIPTION
Laser medium	GaALAs
Laser wavelength	808nm+650nm
650nm Output Laser	12 laser point for 650nm with 5mw
808nm Output Laser	1 laser point for 808nm with 150mw
Timing range	10-30 mins for adjustment
Environment temperature	-20°C-40°C
Relative humidity	≤85%
Heating	40°C-45°C
Power source use	Lithium ion battery with the capacity 1000m Ah
Packing	23.5x6.5x5cm

APPLICATIONS:BACK PAIN, NECK PAIN,INFLAMMATIONS,KNEE ARTHRITIS, RHEUMATOID ARTHRITIS, OSTEOARTHRITIS,ATHLETIC INJURIES, WOUNDS & ULCERS HEALING ,SHOULDER PAIN,TENNIS/ GOLFER'S ELBOW, REHABILITATION PAIN MANAGEMENT (SPRAINS,STRAINS, STIFFNESS) FRACTURES/WOUNDS.ETC.

ΟΡΑΤΑ ΑΠΟΤΕΛΕΣΜΑΤΑ ΣΕ ΔΕΚΑΠΕΝΤΕ ΛΕΠΤΑ ΜΟΝΟ!

ΕΙΔΙΚΕΣ ΤΙΜΕΣ ΓΙΑ ΙΑΤΡΟΥΣ ΚΑΙ ΦΥΣΙΟΘΕΡΑΠΕΥΤΕΣ!

ΤΗΛΕΦΩΝΗΣΤΕ ΓΙΑ ΜΙΑ ΕΠΙΔΕΙΞΗ ΧΩΡΙΣ ΚΑΜΙΑ ΥΠΟΧΡΕΩΣΗ!

AGENTS: MEGATRON LTD. TEL. 00357-99626137