

## Pertussis Diagnostics in Compliance with RKI Guidelines

# Bordetella pertussis ViraStripe® IgG: PT-100


Calibrated with International WHO Standards

According to new guidelines the Robert Koch Institute (RKI) recommends that a "once clearly increased value" for IgG antibodies to *Bordetella pertussis* proves recent infection (1).

The PT-100 band on the *Bordetella pertussis* ViraStripe® IgG is calibrated with International WHO Standards (2) and - when reactive - presents sup-

portive evidence for a recent contact with *Bordetella pertussis* (3).

Thus it is possible to confirm an infection with *Bordetella pertussis* by means of serology if vaccination has not been received 36 months prior to testing.



**Bordetella pertussis ViraStripe® IgG**