



Instruments, ActivaScrew™

| Screw Diameter (mm) | | | | |
|--|--------------|--------------|--------------|--------------|
| Thread Diameter | 2.0 | 2.7 | 3.5 | 4.5 |
| Drill Bit and Tap Diameter (mm) | | | | |
| Drill Bit | 1.5 | 2.0 | 2.5 | 3.2 |
| Tap Thread (AO) | 2.0 (Cortex) | 2.7 (Cortex) | 3.5 (Cortex) | 4.5 (Cortex) |
| Countersink | for 2.0 mm | for 2.7 mm | for 3.5 mm | for 4.5 mm |
| Drive Type | 2.5 Hex | 2.5 Hex | 2.5 Hex | 3.5 Hex |
| + Sleeve | for 2.7 mm | for 2.7 mm | for 3.5 mm | for 4.5 mm |

Instruments, ActivaScrew™ Cannulated

| Screw Diameter (mm) | | | |
|----------------------------|------------|------------|------------|
| Thread Diameter | 3.5 | 4.0 | 4.5 |
| K-wire | 1.25 | 1.6 | 1.6 |
| Drill Bit | 2.7 | 3.2 | 3.5 |
| Tap | 3.5 | 4.0 | 4.5 |
| Countersink | for 4.5 mm | for 4.5 mm | for 4.5 mm |
| Drive Type | 2.5 Hex | 3.5 Hex | 3.5 Hex |
| + Sleeve | for 3.5 mm | for 4.5 mm | for 4.5 mm |