



# TIMPAC HEALTHCARE

Lasers

Dermatoscopy

Radiofrequency

Cosmeceuticals

Diagnostics

Consumables

Liposuction



www.shoptimpac.in www.timpac.in



Timpac Healthcare, an ISO 9001:2008 certified organization offers the Expertise with an Experience of over 38 years in providing solutions to the Medical Fraternity and takes pride in achieving rapid growth whilst exercising ethical practices. Timpac has a strong sales and distribution network with a team of rich experienced exuberant youth, catering to over 6,500 happy customers across India and many abroad. Timpac takes pride in employing the largest pool of trained engineers as Service Experts in the Industry, to ensure that the sales process doesn't only result in a momentary gratifying experience, but a long lasting relationship, between the organization and the client. With the Corporate Head Office in New Delhi, branch offices in Mumbai and representatives in all major cities, Timpac's nation-

wide presence ensures its customers that we are never too far from them.

To gain the long lasting confidence of our clientele, our team follows the motto to offer quality in both sales and after sale service, with the focus owards providing World-class Solutions with the Latest Technology available worldwide, in accordance with the real needs of the customer.

Established as manufacturers of Physiotherapy Equipments, Timpac forayed nto Import and Marketing of Obesity Management technology in the Mid-90's, and then swiftly expanded it's trading arm to various branches of Aesthetics and Medicine, to become the leading company in the narrow niche that Timpac Specializes in. With Vast Experience, Strong Industry Relationships, Global Exposure and Strategic Alliances with Global Leaders, Timpac offers Complete Solutions and Implementation for:

- Dermatology Clinics
- Cosmetology Clinics/Departments
- Clinical Obesity Management Facilities
- Wellness Clinics & Medi-Spas
- Anti-Aging Clinics

• Timpac maintains it's over decade long leadership in the industry by consistently offering more to it's clients, whether it is through The International nstitute for Training and Research in Aesthetic Wellness (IITRAW), a training

ademy, with faculties from around the world, offering training programs for young and seasoned professionals. Each endeavor makes Timpac the most trusted resource for most of the leading Doctors, Institutions, and Corporate in the core areas of Timpac's operations.

Consulting and knowledge sharing are the core of Timpac's operations, utilizing every bit of the expertise gained by our experienced consultants in Conceptualization, Implementation, Training and Marketing of complete Projects, understanding the diversity and uniqueness of each project, with customized solutions, keeping in mind Vision of our Clients and the Satisfaction of their Clients.

TIMPAC HEALTHCARE PVT. LTD.

ACHIEVEMENT 2022

**BEST MEDICAL LASERS &** 

**EQUIPMENTS COMPANY** 

Ас

### **NEW IN 2023**

- Dermatoscope Dermlite DI5
- Dermatoscope Dermlite Lumio 2
- Timpac Freezpen USA
- Black Tree Hair Regrow Solution
- Layun Mono/Cog Threads
- Vip Italy Non Surgical Slimming
- Xtrac Excimer Laser USA TANK

- 577nm Pro Yellow Laser For Vascular Lesion & Melasma
- Monolith Asclepion Diode Laser 5000 Watt
- Amity L Hydrafacial
- Acne Scar / Keloid With Kentamed
- Kentamed RF 100/120 (4 Mhz)/ RF Surgery with Tightening Probes
- Timpac Excimer USA TANK
- Timpac Co2 Laser
- MCL 31 Dermatology, Gynaecology, Snoring
- BION PEELS & SERUMS

# **Timpac Healthcare - Contacts**

Corporate Office (Delhi):

C-168, 1st Floor, Naraina Industrial Area Phase-I New Delhi-110028 (India). Tel.: +91-11-45673300/20/40, +91-9811083717

Mumbai Office: Flat 4, Tulsi II, Four Bunglows Versova, Andheri (W) Mumbai - 400058

### info@timpac.in | www.timpac.in

Place Your Orders orders@timpac.in +91 9312234734

篇

Place Orders Online www.shoptimpac.in

Service Support service@timpac.in +91 9811097837 Customer Care wecare@timpac.in +91 9311758904

Account: Timpac Healthcare Private Limited

Bankers: ICICI Bank, E-47, Naraina Vihar, New Delhi - 110028

A/c#033505003785 | IFSC Code: ICIC0000335 | Account Type-Current

### Timpac WeCare - For you, WeCare!



# Enter Billion Dollar Slimming industry with Timpac

# Non-Surgical to Surgical Solution



Your way to
Perfect Body
with
Combination
Therapy



**CRYOSLIM** 



Vip LIPOLine DUO
Page No. 13



Complete project with consultancy and training





# World's First Variable Polarization

System In A Hybrid Dermatoscope

# **DermLite**

DL5

- •super practical with 100 mm ruler
- •16x white,
- •4x UV (365nm),
- •4x Pigmentboost,
- •1 torch LED.
- ø32mm 10x
- MCCTMadapter

Page No. 38





**RF100** 









**NEW GENERATION** of **RF** for DERMATOSURGERY & COSMETIC DERMATOLOGY

Page No. 35-36













Page No. 30-32









Page No. 7-8

# Contents...

S.No.	Particulars	Page No.
1	Asclepion MedioStar Monolith (Diode Hair Removal)	1-2
2	Asclepion NanoStar (Nano Second Q-Switched Nd:Yag Laser)	. 3 – 4
3	Asclepion PicoStar (Pico Second Q-Switched Nd:Yag Laser)	5 – 6
4	Asclepion Pro Yellow Laser (Vascular & Melasma Laser)	7 – 8
5	Asclepion MCL 31 Dermablate (Er:Yag Skin Lesions, Vaginal Tightening, Snoring)	9 – 10
6	Asclepion Bodylab (Weight Management & Body Shaping)	11 – 12
7	VIP Lipoline DUO (Non Surgical Slimming & Weight Loss)	13
8	Moller (Power Assisted LipoSuction)	14
9	UV Photo therapy	15
10	Excimer (Vitiligo Laser)	15-16
11	Hydrafacial	17
12	Microdermabrasion	18
13	Motorized Bed & Medical Devices Accessories	18
14	Economy Lasers Range	19-20
15	Explore (Microneedling)	21
16	Dana Hiqueen (HIFU Laser)	21
17	Cortex (Skin Ultrasound, DSM III Colormeter, Cryocan)	. 22-24
18	Layun (Thread Lift)	25-29
19	Bion (Peels & Serum)	30-32
20	Timpac Acugel (Ultrasound Assisted Dermaelectroporation)	33
21	Ecleris (Microdermabrasion with Dermaelectroporation)	33
22	PRP	. 33
23	Mesogun	34
24	Laser Accessories	. 34
25	Kentamed (Radiofrequency Surgery & Face Tightening)	35-36
26	Timpac Freezpen	. 37
27	Dermlite (Dermatoscopy DI5)	38 – 40

# New MEDIOSTAR State-of-the-Art Laser Workstation for

Medicine and Aesthetics



MeDioStar





THE DIODE LASER COMBINATION WITH 760 AND 1060 NM



# **APPLICATIONS**

	210, 1110110
0	Hair removal
•	Vascular treatments
•	Acne treatment
•	Pigment removal
•	Skin tightening
0	Skin whitening



### The MeDioStar features:

The new MeDioStar® offers an innovative, easy-to-maintain and expandable modular design with an intuitive user interface, efficient skin cooling and a two-handpiece system.

The new, lightweight and ergonomic Monolith handpieces are setting new standards in aesthetic laser medicine, ensuring maximum stability. The specially developed technology, design and treatment concept featuring a 360° contact skin cooling system as well as exclusively selected materials enables a long service life and absolute safety. There is a suitable solution for every need and numerous possible combinations thanks to the two-handpiece port system. The perfect tool for every user's day-to-day work.

The New MeDioStar® can be equipped with the following Monolith handpieces:

	Monolith S	Monolith M	Monolith L	Monolith XL
Spot size	1cm² (10x10 mm)	1.5cm² (15x10 mm)	3 cm <sup>2</sup> (30x10 mm)	10 cm <sup>2</sup> (31.5x31.5mm)
Power	max. 1,275W	max. 1,700W	max. 3,400W	max. 5,000W
Fluence	1-50 J/cm <sup>2</sup>	1-60 J/cm <sup>2</sup>	1-60 J/cm <sup>2</sup>	1-25 J/cm <sup>2</sup>
Pulse duration	3-400ms	3-400ms	3-400ms	3-400ms
Wavelength	810/940nm	810/940nm	810/940nm	810/940nm





The best choice for the quick treatment of dark hair types is the PowerLine ALX with a 3 cm<sup>2</sup> spot size and max. 5,000 watts of power.



The best choice for the quick treatment of dark skin types is the PowerLine YAG with a spot size of 3 cm<sup>2</sup> and a maximum output of 4,000 watts.



With a spot size of 1 cm<sup>2</sup> and a maximum output of 1,000 watts, the Monolith ALX is the simple and inexpensive way to get started with professional hair removal

### **Specifications**

- Maximum power combined with the shortest pulses.
- Unique wavelength mix 810/940 nm.
- Maximum speed for all skin types.
- Monolith handpieces with 360° contact cooling.
- Plug-and-play system for a fast set-up.
- Unique mechanics for easy maintenance.
- Ultra-sharp 10.1" LCD touchscreen display
- Consumable hygienic single use cap

Optional handpieces 755nm **ALEXANDRITE** 1.4cm<sup>2</sup>



Wireless Footswitch

# Man Star Family









# **QUALITY**



# STATE OF THE ART TECHNOLOGY

The NanoStar® Family comes with a variety of choices of flattop and fractional handpieces. The OptiBeam II technology guarantees a highly precise beam profile and the fractional handpieces expand the treatment options. *DF* stands for «Deep Fractional» and it's particularly indicated for skin resurfacing and scars procedures. *HC* stands for «High Coverage» and it's particularly indicated for the treatment of photoageing and pigmentations.

OptiBeam II

# Advantages

- Big and clear 10.4 inch LCD touch screen
- Big buttons and symbols that are clearly visible
- 습 Dark background for higher contrast and colour brilliance
- 🕒 Clear lines and rounded edges, similar to designed device





9 mmØ HC 8 mmØ DF 5x5 mm 4x4 mm 3x3 mm 2x2 mm

# TECHNICAL SPECIFICATIONS

Laser Nd:YAG, class 4

Wavelength 532 nm; 1064 nm

Frequency Max. 10Hz

Weight Approx. 150 kg

10.4" LCD Touchscreen

532nm:6 ns (QS)

6 ns + 6 ns (OP) 1064nm: 6 ns (QS)

6 ns + 6 ns (OP) 300 μs (PT)

532nm: 12.5 J/cm<sup>2</sup> (QS) 15 J/cm<sup>2</sup> (OP)

1064nm: 25 J/cm² (QS)

37.5 J/cm<sup>2</sup> (OP) 50 J/cm<sup>2</sup> (PT)

2x2;3x3;4x4;5x5mm<sup>2</sup>square

Fractional Ø 8 mm (DF) with typical cover rate 3% - 10% (depending on wavelength) Fractional Ø 9 mm (HC) with typical cover rate 20% - 40% (depending on wavelength)

Pulse

duration'

Max. Fluence

53 cm x 108.1 cm x 110.2 cm (W x D x H)

\* QS = Q-Switched; OP = Opti-Pulse; PT = Photo-Thermal(All specifications are subject to change without notice)



Display

Handpieces

Dimensions











# 




Skin rejuvenation


300ps

2.6 **GW** 

PULSE DURATION

POWER





# QUALITY

# STATE OF THE ART ENGINEERING

With the PicoStar® Asclepion has created a completely new handpiece. The ergonomically and lightweight design offers the physician an easy tool for his daily work. With special developed optics the laser creates a very precise beam profile.

MADE OF STAINLESS STEE

ZOOM ST HANDPIECE

Advantages

- Up to 4 wavelengths
- Highest energy and shortest pulses
- Homogeneous spot
- Automatic spot detection
- No pre-pulses
- No consumables




MODERN & ERGONOMIC DESIGN

DYE 585 HANDPIECE



DYE 650 HANDPIECE





# TECHNICAL SPECIFICATIONS

Laser	Nd:YAG, class 4		Spot Size	Zoom ST	2 - 5 mm Ø
Wayalan ath	Nd:YAG	532 nm, 1064 nm	Spot Size	DYE*	2.5 x 2.5 mm <sup>2</sup>
Wavelength	DYE*	585 nm, 650 nm	Frequency	0.5 - 10 Hz	
Pulse Duration	300 ps		Display	10.1" LCD with Touchscreen	
Energy	Max. 800 mJ		Dimensions	42 x 116 x 106 cm <sup>3</sup> (V	V x D x H)
Power	2.6 GW		Weight	approx. 160 kg	

(\* optional, availability expected Q1 2020)











(All specifications are subject to change without notice)



# SAY NO TO PULSED DYE LASER!

# **577nm Pro YELLOW LASER**



# Vascular Diode Laser multiple applications in dermatology

## The newest generation of "vascular" lasers

First argon laser for vascular treatments was introduced by Asclepion in 1984. The new **QuadroStarPRO** will be the first Table-Top-Laser for aesthetic laser treatment, available in four different wavelengths (green, yellow 940 & 980nm).

made in

Germany

- High performance and reliability
- Maximum for haemoglobin absorption (577 nm)
- Low melanin absorption for fewer side effects (577 nm)
- Scanner with integrated skin cooling available (532 nm & 577 nm)
- Homogeneous beam profile for maximum precision
- Small, compact, lightweight systems
- No consumables
- Transport case available

# APPLICATION

- Melasma
- Rosacea
- Hyper Pigmentation
- Port-wine stain

With iPad App



## **Pulse Dye Laser**

Dye has 6 months life Purpura is Produced Required Maintenance Not Portable No Fractional Mode

## **Yellow Laser**

Vs

No Dye required No Purpura No Maintenance Portable Fractional Mode



# Rosacea:

Typical immediate reaction using the sacnner

Lighten up of the scanned area

**INDIAN PATIENT FRIENDLY** 

NO CONSUMABLE | NO MAINTANANCE

### **QUADROSTAR PRO**



### **Application - Advantages of Yellow**

- Vessels heated to 70°C will get damaged
- Hot blood damages the inner wall of the vessel (intima)
- Very high absorption in blood causes heat very efficiently (with lower fluence)
- If less fluence is needed, shorter pulses are possible, which are additionally more effective for fine vessels

Using the pulse duration of milliseconds with 577nm YELLOW LASER there is no purpura as with the short pulsed Dye lasers

### **Specifications**

• Spot size: 1 mm, Area maximum 15 x 15 mm

Scan density: high

(100%; Spot by Spot; all shapes)


medium

(80%; distance 0.2mm; all shapes except point)

minimal

(60%; mild; distance 0.4mm; only big shapes)

Shapes

















Lentigo







# **MCL 31 Dermablate**

# THE WORLD'S BEST & MOST SUCCESSFUL

# Gold Standard

### INNOVATIONS FOR YOUR MEDICAL PRACTICE

The MCL31 Dermablate has a modern design and is equipped with latest laser technology. The system can be used in a wide variety of following specialist areas.

### **DERMATOLOGY**

Skin Resurfacing Periorbital Wrinkles Perioral Wrinkles Acne Scars **Exophytic Scars** 

### **GYNAECOLOGY**

**Vulvovaginal Atrophy** Dyspareunia Itching and burning Loss of elasticity Mild stress incontinence

### **ENT MEDICINE**

**Snoring** 

**SURGERY** 

Cutting

Xanthelasma








Advanced Skin Resurfacing

**ACNE SCARS** 













# MCL 31 Dermablate



# HANDPIECES UPGRADE

Juliet the feminine lase

45° Reflector, MicroSpot 9x9 mm, 169 MicroSpots 425 μm

Steri-Spot

optional: 4:5 mm optic adapter 360° cone handpiece



DISPOSABLE / REUSABLE HANDPIECE AVAILABLE



# Acne Scar & Vaginal Tightening

### **APPLICATIONS**


- Fractional Skin Resurfacing
- Epidermal and Dermal Lesions
- Acne Scar Treatment
- Vaginal Tightening (VT)

**VT** procedure takes **less than 15 minutes** and can be performed without need for anesthesia.

**The attenuation of the symptoms** can be observed immediately after the first treatment.








VAGINAL ATROPHY AND DRYNESS
VAGINAL RELAXATION SYNDROME (VRS)
STRESS URINARY INCONTINENCE (SUI)

# ROWED

# **ANTI SNORING TREATMENT**



ROMEO® REPRESENTS A FAST, GENTLE AND EFFECTIVE TREATMENT FOR PATIENTS WHO DO NOT WISH AN INVASIVE TREATMENT OR IN CASES WHERE OTHER TREATMENT METHODS ALREADY FAILED




VarioTEAM (XL)	СПТ	MicroSpot

# TIMPAC'S



WEIGHT

WAY OF WEIGHT MANAGEMENT & BODY SHAPING









# BodyLab









# **Treatment**

# areas









Arms

# TECHNICAL SPECIFICATIONS

 Technology
 Electromagnetic field

 Intensity
 Up to 2.5 T


 Frequency
 1 Hz to 150 Hz


 Pulse duration
 200 μs to 300 μs

 Display
 10.4" LCD with touchscreen

Dimensions 34 cm x 90 cm x 67 cm (W x H x D)

Weight Approx. 60 kg











Vip LIPOLine DUO

A MUST BUY for Non Surgical SLIMMING & WEIGHT LOSS Complete Body Shaping & Weight Management System

**Transion:** Slimming & firming **Isogei:** Quick and complete firming **Linfogei:** Lymphatic and venous

**Cavitation** 

The VIP Srl, an Italian company, manufactures and distributes equipment elettro estetiche computerized technology for building the centers of Aesthetics and Aesthetic Medicine. Founded in Turin 40 years ago, holds the record of the technology more 'advanced for the scientific treatment of the most' popular aesthetic and medical aesthetic blemishes.

Since 1970, 400,000 customers worldwide

# NON SURGICAL BODY CONTOURING ADVANCED TECHNIQUES

ADVANCED ISOTONIC ELECTRONIC EXERCISE

### **FATTY ABDOMEN AND THIGHS**







1<sup>st</sup> Session



**CRYOSLIM** Cryolilpolysis

Recognized as the best Quality/ Price Ratio

Sculpt your Body with Cryo Lipo No Consumable, Non Surgical Profitability reach in 1 to 4 months

# **CRYOSLIM Technology**

Unique device with 4 thermal extractors (Peltiers)



Optimal security thanks to healthy and controlled temperature



# Modular system for vibration-assisted liposuction





Vibration liposuction is a tried-and-tested method for fat removal. With this method, the suction cannula is vibrated axially which, in combination with the suction vacuum, cuts and removes the fatty tissue microscopically and gently. The vibrations enable the surgeon to work with greater precision and the separated fat cells are then easier to suction off. This technique offers advantages for both the surgeon and the patient.

- •Gentle and homogenous fat removal
- Time savings
- •Less force is required to move the cannula, the
- •surgeon is thus able to work with greater precision.





# **Liposat Basic**

The infiltration pump Liposat® power has been specially developed for tumescent local anaesthesia. The pump rate can be adjusted to the specific application. Once the set volume has been reached, the pump stops automatically.

## **Vibrasat Power**

The vibration drive Vibrasat® power makes it easier to suction away fat thanks to very fast and short longitudinal movements. The process is therefore less traumatic for patients and less tiring for the doctor. The ergonomic form of the hand grip makes working extremely pleasant.

### **Vacusat Power**

The surgical suction device Vacusat® power is an extremely powerful and untiring piece of equipment which is both quick and simple to operate. Thanks to its innovative insulation technology, it is audibly quieter when operated.







# **VITILIGO & PSORIASIS**



# The Excimer Lamp

# Targeted Phototherapy

Promising Solutions for Psoriasis, Vitiligo Leukoderma & Atopic Dermatitis.



# **Powerful Targeted Phototherapy 308nm**







Courtesy: Dr. Niteen Dhepe

- Monochromatic 308nm Ultra- Narrowband Excimer Light - Proven therapeutic wavelength
- 150 mW/cm<sup>2</sup> Irradiance
   Higher power for faster treatments
- Spot Size 18.9cm² (6.1cm x 3.1cm)
   -Large spot size enables faster, more convenient Tx
   -Targeted light minimizes exposure to healthy skin

### Users:

Dr. Niteen Dhepe, Pune PGI Chandigarh DR. Marwah, Mumbai

# **UV Phototherapy**



Model	No. of Lamps	Power KVA
TOJO 24A	24 UVA	6
TOJO 16NB	16 NBUVB	4
TOJO 24NB	24 NBUVB	6
TOJO 36NB	36 NBUVB	10
TOJO 12+12ANB	12UVA +12NBUVB	4
TOJO 24+12ANB	24UVA +12NBUVB	6
TOJO 18+18ANB	18UVA +18NBUVB	10

# **UV** Comb



# **308nm Excimer System**





# Atopic Dermatitis I Vitiligo I Psoriasis I Eczema

CI atoms and Xe atoms in the sealed handpiece form the unstable XeCI excimer under the action of a high frequency electric field, and this excimer is extremely unstable to a naturally radiate 308nm single wavelength excimer light in the process of transferring from quasi-molecular state into atomic state, which is mainly used for psoriasis, vitiligo and other diseases treatment.





The patented optical power of constant output technology: completely breaks the decreasing rule of optical power with long time use, and really realizes the constant output of optical power. During the effective period of the light source, the output optical power does not decrease, which can ensure the clinical effects.






### Double handpieces design:


- Two handpieces can be simultaneously perform treatment to improve treatment efficiency and utilization and meet clinical needs;
- Reduce the downtime possibility caused by handpiece life or failure
- · Small, flexible and Low cost
- Consumable control and remote control can support various business models
- 308nm monochromatic excimer light, more safe and effective.

### **TECHNICAL SPECIFICATIONS**

	5.1
Laser Type	Eximer System
Model	308
Wavelength	308nm monochromatic excimer light
Optical Power Density	50mw/cm <sup>2</sup>
Spot size	40*40 mm, 16 cm <sup>2</sup>
Energy Density	50-6000 mj/cm <sup>2</sup>
Display Size	8"
Dimension	548X325X192mm (L x W x H)
Power Requirement	100-240 V, 0.8 A, 50/60 Hz
Compliance	FDA / US 510 k*
	EC Medical Device Directive
	(MDD) 93/42/EEC (CE mark)

# **Excimer Light Theory**











Ice head



jet vacuum

0 - 8 Bar

Big Display 8" Screen

Vacuum Power 1 Bar Max.





Hydera Dermabrasion Handle






Facial Ultrasound Handle



Specifications





Electric Microneedle



Diamond Dermabrasion





Skin Scrubber



Liquid Bottle



PDT Light

Spraygun

Liquid Collection Bottle



## TIMPAC MOTORIZED FACIAL BED



### Beauty bed with 3 motors

High comfort bed equipped with 3 motors that control the height, the inclination of the back-rest and the inclination of the seat. Headrest with breather hole incorporated. PU upholstery of high quality and easy cleaning.

Up and Down 59 / 98 cm

Length: 184 x 60 cm (with armrests 85 cm)

Weight: 88 Kg

Other colors available on demand

# **ECLERIS**



## **MEDICAL DEVICES**





Disposable plastic hub needle Box of 100



# **Elliptical Punch**

Sterile, Disposable Excisional Scalpel

**Round Punch** (Biopsy)



Elliptical Incision Closes With NO "DOG EARS"







**PUNCH<sup>®</sup> KITS** 

A Complete, Disposable Kit for Skin Biopsy Procedures

# **CURETTAGE**



**Available In A Full Range Of Sizes** 

1mm 2mm 3<sub>m</sub>m 5mm 4mm

7<sub>mm</sub>

# SURGICAL INSTRUMENTS



The most-used dermatologic surgical instruments



# Lasers & Aesthetic Equipments AN ECONOMY RANGE FROM TIMPAC

# Nd Yag 3 IN 1 Laser Platform **Q-Switched**

IPL LASER (OPT + SHR) RF



Illuminant **IPL** 

Spectral Range 480/530/640-1200nm

Illuminance Area 15\*50mm **IPL** Energy 1-50|/cm<sup>2</sup> Adjustable pulse width: 2-99ms Pulse Number 1-5 1-50 **RF Energy** 

5Mhz(3 treating probes) RF Frequency Bipolar

Laser Wavelength 1064nm/532nm Q Switch Nd:YAG Laser Type

50 Laser Energy Laser Frequency 1-5Hz

Display Screen 8 itch color touch screen

Cooling system Water cooling +

semiconductor

refrigeration, air cooling

**Total Power** 2500W

AC100V-220V Voltage



**Q-switched** nd yag Laser With **Articulated** Arm

The wavelength 1064nm/532nm 100-2000mJ(1064nm) Energy range The screen 10.4 TFT touch screen

Spot size 2-10mm Repetition frequency 10Hz The pulse width 5-7ns

Continuous working time Can work continuously for

24 hours(when the room

Air cooling and water cooling **Cooling System** 

Temperature requirements  $10^{\circ}\text{C}{\sim}30^{\circ}\text{C}$ Humidity ≤85%

Voltage 220v±10%,50/60HZ

**INTENSE PULSE** 

LIGHT

RF + IPL E-LIGHT

## **Specifications**

Wavelength (Spectrum) 640-1200nm(HR) Energy Density (Fluence) 10-60J/cm2 Up to 50W **RF Energy** 15\*50mm2 (HR) Spot Size **Pulse Repetition Rate** 0.2Hz; 2Hz **Pulse Duration** 0.5-15ms 1-50ms Pulse Delay

Pulses 1-15 pulses IPL Peak Power 1400W RF Power Max. 100W

Cooling Continuous Crystal contact cooling (-4  $^{\circ}$ C  $^{\sim}$ 4  $^{\circ}$ C ) +

Air cooling + Closed water circulation cooling



# Lasers & Aesthetic Equipments AN ECONOMY PANCE EDOMATION

AN ECONOMY RANGE FROM TIMPAC

# Diode Laser 3 WAVELENGTH (755nm, 808nm And 1064nm)

## PAINLESS DIODE LASER HAIR REMOVAL

**Specifications** 





### Laser type Diode laser Laser wavelength 755nm,808nm&1064nm 10.4" color touch LCD display Display 1500W Output power Spot size 12 x 12 mm<sup>2</sup> Pulse width 8-685 ms adjustable Energy 2-120 J / cm<sup>2</sup> adjustable Frequency 0.5-10 Hz Cooling mode close-cycle water cooling water + air + semiconductor Cooling Temperature of probe 0~3°C 59\*59\*146cm Dimension of machine Dinension of package 70\*70\*126cm Net weight 45 kg Gross weight 61 kg AC 220V/50Hz;AC 110V/60Hz Voltage Packing Aluminum alloy case

# **Q-switched Nd:yag Laser**

System For Tattoo & Pigmentation



### **Technical Specification:**

Light Source	Picosecond Q-switch Nd:YAG laser
Controlling system	8.4 inch color touch screen
Wavelength	532nm, 755nm, 1064nm
Output Energy	Single Pulse: 532nm 400mJ; 1064nm 800mJ Single Pulse: 532nm 800mJ; 1064nm 1600mJ
Energy density(Max)	1064nm: 15.0; 532nm: 6
Pulse duration	2000ps
Repetition Rate	1-10Hz adjustable
Spot size	0-10mm adjustable
Cooling system	Closed cycle water to air heat exchanger
Weight	24kg

# Fractional CO<sup>2</sup> Laser

RF-activated CO2 Fractional laser Combines Ultra-pulse and fractional laser technology

10,600nm Wavelength Power Output 1-40W

UltraPulse; CW; Fractional Laser Type 0.12mm & 1.25mm (adjustable) Spot Size

Spot Density Up to 102400 dot Scan Size Up to 20\*20mm Scan Shapes

Square; Rectangle; Circle; Triangle; Rhombus; Ellipse; Line

Beam Patterns Standard; Random; Scatter Aiming Beam 5mW red diode laser, 635nm,

adjustable intensity, blink on/off

Beam Delivery Articulated 7-joint arm

Cooling

with up to 360° rotation Self-contained, closed cycle





# **EXPLORE**

## **Minimally Fractional RF Microneedeling**

4mm Needle Depth, 60 Watt, Maintenance Free

## **Invasive & Non-invasive**

The Harmony of 2 Way Hand-Piece



### **Gone with Explore**

- Aging skin (Lifting, Elasticity, Firmness, Anti-Wrinkle)
- Scar Treatment (Acne Scar, Normal Scar, Keloid)
- Acne, Goose bum skin, Striae distensae (Stretch Mark)

### Biometic RF - Only invasive RF can combine with stem cell or PRP

Power control - each 0.1 sec control Explore enables RF pulse time by	<b>SPECIFICATIONS EXPLORE</b>		
0.1sec and can freely control energy delivered area.	Туре	Electro-surgical Device	
energy delivered area.	Power	1.8W - 60W, OFF	
	ON TIME	1ms-120ms	
	OFF TIME	0.5sec - 3.4sec	
Time Power 30	Needle type	Invasive / Non-invasive	
THE RESERVE OF	Electrical	100-240VAC, 50Hz/60Hz	
	Dimension	330(W) X490(D) X 1175(H)mm	
Time Power 60	Weight	22kg	

### **Users:**

Dr. Punit Rekhi, Punjab, Dr. Ayush Agarwal, Bareilly Dr. Gambhir, (Plastic Surgeon), Delhi Dr. Manju Kesri, Ghaziabad, UP Dr. Leena Bhatnagar, Jaipur Dr. Phanisree, Hyderabad

Dr. Sandeep Gola, Hyderabad

Dr. Mahesh Jadwani, Nagpur Dr. Bandana Bindal, Gurgaon

# HI-QUEEN

# **High Intensity Focused Ultrasound**

Immediate visible effect just after one treatment.



Penetrate deep into the SMAS layer without any skin damage. This causes immediate skin tightening and lifting. With time more collagen and elastin are produced throughout all depths resulting in skin rejuvenation and a V shape face.



### **SPECIFICATIONS HIFU**

Energy Output	0.1-1.5 J/cm2
Cartridge Depth	4.5mm - 4Mz 3.0mm - 4Mz 1.5mm - 7Mz 13mm - 2Mz (Obesity)
Length	24 mm(Full)
Interval	12 mm (Half) 1.0, 1.5, 2.0mm
Weight	35kg

# **CORTEX** Technology

Users: AIIMS, New Delhi & AIIMS Rishikesh

**DermaLab® Combo** is a compact skin analyzer featuring touch screen operation making daily use easy and intuitive.

**DermaLab® Combo** is based on our long-term experience in developing and manufacturing scientific skin instrumentation. For the first time it is now possible to combine the classical topical measurements with high resolution cross-sectional in-depth imaging in one single instrument.

**DermaLab® Combo** offers high precision and reproducibility on all its measuring parameters making it ideal for comparison over time. A built-in client database makes it easy to recall previous recordings and compare with new measurements.

**DermaLab®** Combo easily integrates into the daily work routine, where it may provide supplementary scientific evidence in support of visual assessment of skin conditions, whether it is in a beauty clinic or salon. As such, the device is ideal for evaluating photoaging and the collagen improvement as a result of rejuvenating procedures as well as the recommendation of specific skin care products.

# **DERMALAB®** COMBO

Skin Test Center with Ultrasound

## **Introducing Low Cost Model**





**Ultrasound:** Ultrasound skin-age imaging clearly demonstrates the need for treatment and

documents treatment efficacy. Ultrasound may also be used for visualizing cellulite and the efficacy of cellulite treatment. Principle: Rotating 20 MHz X-tal, 3.7 mm penetration.











**Elasticity:** The visco-elastic properties of skin are assessed by <sup>20 years old 3</sup> means of a light weight suction cup. **Principle**: Stress/strain by applied suction cup.

30 years old 40 years old

years old Co

Elasticity

**Moisture:** Conductivity is a well-known indicator for the water binding capacity of Stratum Corneum the skin hydration level. Pin probe design facilitates dry skin and scalp applications.

Principle: Conductance, single frequency.



Moisture

**Sebum:** Sebum is objectively assessed by measuring the translucency of an oil absorbing microporous film material. Offers on-screen step-by-step guidance and strip offset calibraton. **Principle:** Translucency of sebum collecting material.



Sebum

**Pigmentation:** Calculates the melanin and erythema index to provide information about the pigmentsation level and reaction to irritants and allergy. **Principle:** Colour detection, reflectance spectrophotometry.



Pigmentation

**Transepidermal Water Loss:** Transepidermal Water Loss (TEWL) is an established and widely accepted method to accurately assess the skin barrier function. **Principle:** Open chamber, vapour diffusion gradient.



TEWL

**Video-scope:** The video-scope is an electronic video microscope with 20-50X magnification. With built-in illumination it provides high quality close-up images of skin conditions, which may easily be captured and stored in digital form. **Principle:** CCD, built-in illumination.



Videoscope

## **DERMASCAN® C USB**


Skin Ultrasound

# **CORTEX** Technology

# **DSM III**COLORMETER

2D, Ultra Compact, Ultrasound skin scanner with advanced skin image analysis software











dermis, age 44.





**AVAILABLE WITH ONE PROBE ONLY** 

### **PROBE SELECTION:**

- 10 MHz, 120 x 1000 pm, 30 mm penetration.
- 20 MHz, high definition, 60 x 150 pm, 13 mm penetration.
- 20 MHz, high sensitivity, 60 x 200 pm, 23 mm penetration.
- 50 MHz, 25 x 60 pm, 3 mm penetration.

### **SYSTEM SPECIFICATIONS**

- Swept gain control, level and profile.
- True broadband system, 12 bit A/D conv.
- · Cineloop recording, 100 images.
- Full-featured image analysis package.
- Four simultaneous images.
- Autosave function.
- Real-time scanning with scroll and zoom.
- Meets international regulatory requirements.

### DSM III Skin Colormeter

combines advanced color sensing technology, high sensitivity and durability with extreme flexibility and ease of operation. Proven and highly sensitive technologies for skin colour detection are built into a handheld and light weight device to form a fully integrated instrument, which is operated by just one button. To facilitate correct positioning on the skin a guiding light illuminates the target area.

Measurements are presented on the integrated display, and may also be sent to a PC or tablet for further spreadsheet processing. For this purpose the DSM III features wireless Bluetooth communication.

It offers read-out of erythema and melanin (E&M) based on the absorbance characteristics of human skin as well as CIEL\*a\*b\* values.



### **Specifications**

- Light source: 2 x high intensity white LED
- Colorsensor: 4 x 16 active elements
- Display: LC-display
- Read-out: E&M (erythema & melanin), CIEL\*a\*b\* (optimised for white LED illumination)
- Operation: One-button control
- Built-in Bluetooth for wireless communication with tablet/ laptop
- Battery: 1 pcs. rechargeable, standby time: 1 month
- Charger: Wireless induction
- Calibration: White calibrator included
- Operation temp.: 10-40 Deg. C

# DermaLab USB single parameter units



# **CORTEX** Technology

# **CRYO PRO**

## Liquid Nitrogen Cryosurgical Equipment

At -196 deg.C liquid nitrogen is the most efficient coolant to be used for cryosurgical treatment. Our CryoPro treatment units as well as the comprehensive accessory program are designed to provide reliable and trouble free operation even after extensive freezing due to the use of high quality materials and careful quality control in all production stages.

In addition to the treatment units we also supply liquid nitrogen storage containers and withdrawal devices to conveniently fill the CryoPro units.

The CryoPro units are available in sizes of 0.5 I (Maxi) and 0.35 I (Mini) and come with a three year world wide warranty against manufacturing and material defects.





### **ACCESSORIES:**

Closed contact	sharp pointed	1 mm		2 mm	3 mm	4 mm
probes:	5 mm	6 mm	8 mm	20 mm	30 mm	cervical probes
Open spray tips:	std. tips A/B/C/D	bent spray	straight spray			
LN2 storage:	AL dewars	CD dewars	withdrawal tube	roller Base		
Misc. accessories:	lexan cones	neoprene cones	carrying bag			

Contact probes: Miniature sized, gold plated. For tiny to small lesions - e.g. small flat warts, lesions near the eye or on the eyelid.













Contact probes: Medium sized, gold plated. For small to medium size lesions - e.g. flat warts, small benign tumors, lesions near the eye.



4 mm

10 mm







Contact probes: Large, gold plated. For medium size to large lesions - e.g. flat warts, benign tumors.









Standard spray tips: A - heavy spray for large lesions and superficial desquamation, B & C - medium output for general application, D - fine spray for small lesions.










A aperture (Ø 1.0)



C aperture (Ø 0.55)

ALCOHOL: N

Special spray tips:









Bent spray extension

Straight spray extension

Soft peel desquamation tip

Luer lock adapter (KELOID & MORE)

Special attachments & accessories:









Open flow restrictor

Lexan cones









# ReNew by















- Basic lifting thread
- Neck, buttocks, forearms, and abdomen wrinkle care
- Elasticity restoration
- Collagen production

### Specification

specification								
Needle	Needle	Thread	LICD		Thread Optior	า	Pack	aging
Gauge	Length	Length	USP	PDO	PLLA	PLCL	Pack	Вох
25G	60	90	5-0	0	0	0	20ea	5packs
27G	40	50	5-0	0	0	0	20ea	5packs
29G	30	30	6-0	0	0	0	20ea	5packs
29G	40	50	6-0	0	0	0	20ea	5packs
29G	50	70	6-0	0	0	0	20ea	5packs
30G	30	30	6-0	0	0	0	20ea	5packs
30G	40	50	6-0	0	0	0	20ea	5packs



# ReNew by Mono

Application









- Wrapped in a tornado shape
- Increasing the contact surface on insertion
- Undereye wrinkles and crow feet

Needle	leedle Needle Thre		LICE		Thread Option	Packaging		
Gauge	Length	Length	USP	PDO	PLLA	PLCL	Pack	Вох
25G	60	90	5-0	0	0	0	20ea	5packs
27G	40	50	5-0	0	0	0	20ea	5packs
29G	30	30	6-0	0	0	0	20ea	5packs
29G	40	50	6-0	0	0	0	20ea	5packs
29G	50	70	6-0	0	0	0	20ea	5packs
30G	30	30	6-0	0	0	0	20ea	5packs
30G	40	50	6-0	0	0	0	20ea	5packs















# RENEW EYE

### Application








- Micro U type cannula
- Produces collagen that fills out wrinkles on insertion
- Improving the contour of the eye and the elasticity of the skin

### Specification

Needle	Needle	Thread	LICD		Thread Optior	1	Pack	aging	No adla Tana
Gauge	Length	Length	USP	PDO	PLLA	PLCL	Pack	Вох	Needle Type
G			-	М	М	М	ea	packs	
G			-	М	М	М	ea	packs	

# ReNew by Cog

# CORSET





### Application







- Creating a V-line on the face
- Patented asymmetric structure prevents sagging after the procedure
- Removes wrinkles and provides long-lasting skin elasticity
- Thread is 180° cog, and the retention part (rear half) is 360° dual cog
- W type cannula minimizes tissue damage during insertion

### Specification

	aging	Pack	Thread Option Pa		Thread	Needle	Needle		
Needle Type	Box	Pack	PLCL	PLLA	PDO	USP	Length	Length	Gauge
	packs	ea	В	В	В	-			G
	packs	ea	В	В	В	-			G
	packs	ea	В	В	В	-			G
		ļ.							











Fixation part

# ReNew by Cog

# CORSET

### Application









Maximized lifting eect by "Three-part" asymmetric structure


Retention part

- Unique structure never seen before, divided into three parts to grip, retain and fix
- Anchoring with 180 ° cog, the Retention part (middle) is the 90 ° dual cog
- Fixation part (rear) is the 90 ° cog to fix the thread once more
- Combination with CORSET Black, our other patented flagship product

### Specification

۰	Needle	Needle	Thread	LICD	_	Thread Optior	า	Pack	aging	N II. T
٠	Gauge	Length	Length	USP	PDO	PLLA	PLCL	Pack	Вох	Needle Type
	G		-	D		D	D	ea	packs	





### Application





- Patented asymmetric structure
- Embedding cog that runs in a 360° spiral
- Black raises sagging cheeks and improves facial lines, and is more effective when treated with CORSET
- divided into a grip part that fixes it and a retention part that maintains the shape
- Ideal angle, shape, and size
- More effective and lasts longer

Specification

Needle	Needle	Thread	1100		Thread Option	ı	Pack	aging	
Gauge	Length	Length	USP	PDO	PLLA	PLCL	Pack	Вох	Needle Type
G				G	G	G	ea	packs	
G		-	G		G	G	ea	packs	














# ReNew by Cog DOUBLE-LOC



### Application







Gear wheel like eect


- DOUBLE-LOC is effective for thin skin such as nasolabial folds and forehead wrinkles.
- Inside and outside of the needle are designed differently
- Inside of the needle is made of bi-directional cogs facing each other.
- Outside of the needle is made of densely coiled cogs.

### Specification

Needle	Needle	Thread	1100		Thread Option	1	Pack	aging	
Gauge	Length	Length	USP	PDO	PLLA	PLCL	Pack	Вох	Needle Type
G			-	F	F	F	ea	packs	
G			-	F	F	F	ea	packs	
			•				ļ.		



# ReNew by Cog Plus ® DOUBLE-LOC



### Application





# The specially designed dual needle thread Removes the forehead and eye tail wrinkles

- DOUBLE-LOC Plus has press made rigid 180° molding cog giving it tensile, traction, and durability
- Optimized for thin skin layer treatment such as forehead and chin
- Round Blunt needle minimizes tissue damage and reduces bruising and swelling

### Specification

	Needle Gauge	Needle	Thread	USP		Thread Option	Packaging		
		Length	Length		PDO	PLLA	PLCL	Pack	Вох
	G		-	а		×	а	ea	packs













# ReNew by Cog



Application









- The KO thread has the scaffold structure.
- for Nose treatment (dorsum, tip and wings of the nose)
- Cogs on the inside and outside of the needle are designed in the opposite direction,
- Effective for fixation and traction
- L cannula minimizes tissue damage, bleeding, bruising and swelling.

### Specification

	Needle	Needle	Thread	LICD		Thread Optior	ı	Pack		
ı	Gauge	Length	Length	USP	PDO	PLLA	PLCL	Pack	Box	
	G				а	a	a	ea	packs	
	G				а	а	а	ea	packs	
	G			-	a	а	а	ea	packs	
	G			-	a	a	a	ea	packs	1



# RENEW BY









- Bi-directional structure are essentially monothreads with barbs to hook to the underside of the skin.
- effective thread to lift jawline and restore the young and slimming look.
- L cannula provides the advantage of smooth insertion into skin
- Making the whole lifting procedure totally safe without any bleeding, bruising and swelling.

### Specification

Needle Gauge	Needle Length	Thread Length	USP		Thread Optior	Packaging		
				PDO	PLLA	PLCL	Pack	Box
G			-	а	а	а	ea	packs
G			-	a	a	а	ea	packs
G			-	a	a	a	ea	packs













SKIN REJUVENATION SYSTEM
with Glycolic Acid Peel

20%, 35%, 50% & 70%

**PRE PEEL SKIN CLEANSER** 





# **SKIN NEUTRALISER**



**HYDRAFACIAL REJUVENATION SYSTEM** 







LOTION

SKIN PURE LOTION (BHF 2)

LOTION

### TO BE USED WITH ACUGEL & RF TIGHTENING



SKIN DEEP **BRIGHTENING GEL** 




BRIGHTENING MASK

**Deep Cleaning & Depigmenting System** 





## Salicylic Acid Peel





FACE PURIFYING SERUM Salicylic Acid 20%

### Lactic Acid Peel



BRIGHTENING GEL Lactic Acid 20%, Citric Acid 10%



BRIGHTENING GEL Lactic Acid 10%, Citric Acid 10%

### TCA Peel







SALAZ Peel



### SM Peel





### YELLOW Peel



Launching New Combination Peel Lightening & Anti-Aging serum soon...





# Nourish Well, Grow More



# **TOPICAL**MESO HAIR AQUASOL



Peptides
Growth Factors
Amino Acid
Vitamin
Trace elements



### **HAIR SERUM**

Mulberry Root Extract Biotinoyl Tripeptide Amino Acid Cocktail

**CAN BE USED AS MESOTHERAPHY** 

## **CARBON PEEL**



TO BE USED WITH Q-SWITCHED ND:YAG

### **CREAMS**



**BRIGHTENING CREAM** 



EYE REJUVENATING CREAM

### **LOTION**



ULTRA SOOTHING LOTION



SOFT HEAL LOTION

### **WASH**



FOAMING EXFOLIATING WASH



CLARIFYING SKIN WASH

# **ACU-GEL**

### **Ultrasound Assisted Dermaelectroporation**

• Non-Invasive • Painless • Inexpensive

NEW Application Post Q-Switched Nd:YAG for

**TATTOO** and Pigmentation



TIMPAC ACUGEL is a portable & highly personalized machine. It is completely solid state & very easy to use. It uses a unique form of electroporation alongwith electro mechanical energy which penetrates the tissues very deeply. The energy released into the tissue vibrates mechanically at high frequency for a very small range of movement, at cellular level, the molecules are oscillated about million times a second but for a distance of few millimicrons. The result is electroporation alongwith cellular massage for the cells, improvement in blood circulation, reduction edema. Resultant synergy works much better than electroporation or ultrasound used independently.

#### **APPLICATIONS**

- Under Eye Dark Circles,
- Under Eye Bags,
- Facial Rejuvenation
- Lymphatic Drainage
- Penetration of Ingredients to deep dermalayers (electroporation)
- · Post Procedure erythema







Courtesy:
Dr. Aashish Khurana
New Delhi







Lightening products
Use with Acugel

Brightening Gel

**Brightening Mask** 



#### TYROSINASE INHIBITORS WITHOUT HYDROQUINONE

Brightening Gel Use with Acugel
Brightening Mask Post Acugel Treatment

Made in France I-II Level Lightening In 20 min.

**Composition:** Bearberry Extract, Mulberry Extract, Licorice Extract, Grapefruit Extract, Scutellaria Extract

# ECLERIS

# **DERMA-ELECTROPORATION**Needle-FreeMesotherapy



# Platelet Rich Plasma AUTOLOGOUS CELLULAR REJUVENATION



## **MESO DEVICES**

### **Painless Mesoguns**



#### THE BENCHMARK FOR MESOGUNS






### **DERMAROLLER**






## LASER ACCESSORIES®



#### LASER GOGGLES



19row







# SURGITRON® EMC

**Advance 4.0 MHZ Radiowave Technology** 









# SURGITRON DUAL RF

**Advance 4.0 MHZ Radiowave Technology** 



Surgitron® Dual Frequency 4.0 MHz Radiowave Technology

## RF-B









# **RF Surgery Unit**



#### RF-B - 4MHZ MONOPOLAR & BIPOLAR

for Dermatology, Plastic, ENT, Aesthetic Surgery, Ophthalmology, Neurosurgery, Gynaecology

- Timer can be used in each mode to secure reproducibility of certain procedures
- Inbuilt Cooling

#### Users Dr. J. B Ratti, NEW DELHI

Dr. Prabha Eye Hospital, BANGALORE

Amardeep Diya Eye Care Private Limited, KOLLAM

Pandurang Shasri Athavale Education Trust, MUMBAI

Medical Research Foundation (Sankara Nethralaya), CHENNAI

Aravind Eye Hospitals, Tirunelveli

Sir Ganga Ram Hospital, NEW DELHI

Technical Characteristics			
Mode	Maximal Output Power	Crest Factor	
Fine CUT	130 W / 150 Ω	1.5	
Heavy CUT	130 W / 450 Ω	1.6	
Bi CUT	130 W / 150 Ω	1.5	
Soft COAG	130 W / 50 Ω	1.5	
Forced COAG	100 W / 600 Ω	3.8	
Bi COAG	130 W / 50 Ω	1.5	
Operating Frequency [MHz]	4		
Modulation Frequency [kHz]	33		
Patient Safety	NEMSY (Neutral Electrode Monitoring System)		
Activation	Twin Footswitch and Hand-switch Pencil		
Duty Cycle	10/30		
Safety Class/MDD Class	I/IIb		
Applied Parts Type	BF		
Dimensions [mm]	280x430x150		
Weight [kg]	12.7		



### **NEW GENERATION** of **RF** for DERMATOSURGERY & COSMETIC DERMATOLOGY

4 Mhz

100 Watt

Warranty

RF100 - 4MHZ MONOPOLAR & BIPOLAR with Face & Body Tightening Electrodes\* & Insulating Needle\* for White Hair

More than 200 Machines Sold

Partial list of Users of Kentamed

DR. MOHAN SINGH, PUNJAB DR. P.S ANEJA, JALANDHAR DR. MAHESH JADWANI, NAGPUR

DR. SUBODH BANSAL (ABHAR) DR. DEEPAK SHARMA (REWARI)

DR. K.S. KUMAR, BANGALORE DR. VIJAY GOYAL MONGA, PUNJAB

EYE INSTITUTE, HYDERABAD MAX HOSPITAL, DELHI

MAX HOSPITAL, DELHI ESI HOSPITAL, DELHI GOVT. MEDICAL COLLEGE, LATUR DR. ANUJ PALL, GURGAON DR. DHANANJAY CHAVAN - PUNE

DR. BALACHANDRA S ANKAD - KARNATAKA DR. BINOD KUMAR THAKUR - ASSAM DR. HARSHIT RANPARA - RAJKOT

DR. DEBJIT KAR - ORISSA DR. ABHISHEK KUMAR SINHA - PATNA

DR. VIMALA SKIN HAIR & LASER CLINIC - HYDRABAD DR. PANJABRAO ALIAS BHAUSAHED - AMRAVATI

SIR GANGA RAM HOSPITAL, DELHI



Users: Dr. Bipin Despande, Pune

**( (** 0044

#### **6 Modes of Operation**

- No Carbonization
- Heavy Cutting Cut with Coagulation
- **Soft Coagulation** • Forced Coagulation - Fulguration
- Very Fine Cutting
- Insulated Needle for White Hair

#### Radiofrequency by Dr. Bipin Deshpande





Scarless Excision of Mole





Virtually scarless excision





Scarless Excision of Actinic Keratosis





**Excellent Removal** of Warts

Technical Characteristics			
Mode	<b>Maximal Output Power</b>	Crest Factor	
Fine CUT	70 W	1.5	
Heavy CUT	100 W	1.5	
Soft COAG	50 W	2.0	
Forced COAG	50 W	2.8	
Bi COAG	80 W	1.5	
Micro Bi COAG	35 W	1.5	
Operating Frequency [MHz]	4		
Patient Safety	NEMSY (Neutral Electrode Monitoring		
	System)		
Activation	Twin Footswitch and Hand-switch Pencil		
Duty Cycle	10/30		
Safety Class/MDD Class	I/IIb		


#### Standard Set of Accessories for Dermatology & Plastic Surgery

Applied Parts Type

Dimensions [mm]

Weight [kg]

- 1 reusable finger-switch handle with 3 m cable
- 5 monopolar electrodes
- · Reusable earth cable
- 1 bipolar cable 3 m
- 1 bipolar forceps • 1 twin footswitch
- Face tightening electrode (optional)
- Body tightening electrode (optional)
- Insulating Needle for White Hair (optional)





Other Accessories – Available Upon Request

BF

280x150x280

8.6

# Timpac Freezpen Made in Belgium



Angioma
Milii
Acne (scars)
Acrochordon
Actinic cheilitis
Actinic keratosis
Angiofibromas
Callosity

Dermatofibroma Flat warts

Horns Leukoplakia

Leuropiaria Tilomi

Syringoma

Stuccocheratosi

Vulgar warts

Warty carcinoma

Lentigo solar

Milium colloid

Mucoid cyst

Plantar warts

Rosacea

Piercing elastosis

Lichen planus

Lichen sclerosus

Lichen simplex

Capillary hemangioma

Cutaneous leishmaniasis

Cutaneous sarcoidosis

Granuloma fissuratum

Herpesvirus infections

Infectious mollusk

Keratoacanthoma

# What Indications can be treated?

**Approved** 

Warty epidermal nevus Xantomi Kaposi Ano-genital warts Annular granuloma Nevus sebaceous Nodular prurigo Porokeratosis linear Pruritus anirhinophyma Pyogenic granuloma Dermatosis papulosa nigra Facial eosinophilic granuloma Fibrous papule of the face Filiform and typed warts Angiocheratoma of Mibelli and the mammary areola Melanotic macules of the lips Perforating dermatoses

Seborrheic keratosis

### **Applicators**



Angiocheratoma of the scrotum and vulva
Dermatoses due to mechanical traumas
Fixed discoid erythematosus lupus
Idiopathic guttata hypomelanosis
Keloid and hypertrophic scars
Papillary hyperplasia of the palate
Nevoid hyperkeratosis of the nipple
Nodular chondrodermatitis of the helix
Viral lesions in HIV positive and negative patients

cryosurgery made simple

# **SKIN DIAGNOSTICS**



# **DermLite**

# Accurate Diagnostic Leads To Precise Treatment

# DL5

World's First Variable Polarization

System In A Hybrid Dermatoscope



- Rechargeable, includes charging base
- super practical with 100 mm ruler
- 16x white,
- 4x UV (365nm),
- 4x Pigmentboost,
- 1 torch LED.
- ø32mm 10x
- MCC<sup>™</sup> adapter

# DL4

**World's Best-Selling Dermatoscope** 



- Polarized & Non-Polarized
- 18 white LEDs (12 Polarized,
   6 Non-Polarized), 6 Orange
   Polarized LEDs
- IceCap™
- Pigment Boost™ Plus
- ø30mm 10x

# **SKIN DIAGNOSTICS**



# **DermLite**

# Accurate Diagnostic Leads To Precise Treatment

### **DL200 HYBRID**



#### **Pocket Dermatoscope**

- Polarized & Non-Pol.
- 21 white LEDs (15 pol., 6 non-pol.)
- IceCap™
- PigmentBoost™
- ø25mm 10x
- Longer Battery Life
- Infection Control Caps

## **DL3N**

- 25 mm four-element lens,
- 28 high-powered LEDs and the all-new
- PigmentBoost™ illumination.
- Polarized & Non-Polarized.
- 18 white LEDs (12 pol., 6 non-pol.), 10 orange pol. LEDs
- IceCap™
- PigmentBoost™
- ø25mm 10x



### **Foto X Plus**

- 16.1MP Four Thirds Live MOS Sensor
- Dual Quad-Core TruePic VIII Processor
- UHD 4K30 Recording; Clips Video Mode
- 2.36m-Dot 1.23x Electronic Viewfinder
- 3.0" 1.04m-Dot Tilting Touchscreen LCD
- In-Body 5-Axis Image Stabilization
- 121-Point Contrast-Detect AF System
- Up to 8.6 fps Shooting and ISO 25600
- Digital ED 14-42mm f/3.5-5.6 EZ Lens
- 16GB SDHC Memory Card



## **Foto II Pro**

- Works with any Nikon or Canon full-frame or 2/3 frame DSLR camera (with interchangeable lenses) out of the box
- Unsurpassed edge-to-edge clarity and resolution
- Toggle between polarized and non-polarized illumination via button push or automatically via flash cable
- On-board rechargeable battery
- Integrated zoom with secure zoom lock
- Removable magnetic spacer and removable glass contact plate



# **SKIN DIAGNOSTICS**



# DermLite

# Which Dermlite Suits Your Practice?





### Universal Adapte Smartphone



MCC CONNECT CLAMP



**Iphone Adapter** 









Page No. 25-29

## **MAXIMIZE**

YOUR RESULTS WITH COMBINATION OF

# **KENTAMED & FREEZPEN**



Page No. 35





# **Lasers & Aesthetic Equipments**

AN ECONOMY RANGE FROM TIMPAC



RF + IPL E-LIGHT



Nd Yag 3 IN 1 Platform Q-Switched IPL LASER (OPT + SHR) | RF



**Diode Laser** 



Q-Switched Nd: Yag with Articulated Arm



Fractional CO2 Laser Q-Switched Nd: Yag Laser



Page No. 18-19



Page No. 20-24



Glycolic Acid | Salicylic Acid | Yellow Peel













# Skin & Body Solutions

Non-Surgical to Surgical

Lasers | Dermlite Dermatoscopy | Radiofrequency
Diagnostics | Consumables | Liposuction
Bion Peel & Serum





knowledge thru practice a unit of joto to to the second results of the second results of

International Institute of Training & Research

Cosmetic Dermatology • Aesthetic Surgery • Wellness

COURSES ON COSMETIC DERMATOLOGY
AESTHETIC MEDICINE AND WEIGHT MANAGEMENT \*

### One Stop Shop For Medical Lasers & Cosmeceuticals Trusted by Over 6000 Doctors



























Service Helpline service@timpac.in +91 9811097837

Place Your **Orders** orders@timpac.in +91-9312234734



Customer Care
wecare@timpac.in
+91-9311758904

Account: Timpac Healthcare Private Limited Bankers: ICICI Bank, E-47, Naraina Vihar New Delhi - 110028 A/c # 033505003785 IFSC: ICIC0000335 Business Hours on weekdays: 9.30 am to 6.00 pm.



C-168, 1st Floor, Naraina Industrial Area, Phase -I, New Delhi - 110028, INDIA Tel: +91-11-45673300/20/40, wecare@timpac.in | +91 9811083717