


Implant Catalogue

MENTOR® BREAST IMPLANTS

www.mentormedical.co.uk
Customer Service 0800 973 082
ccontacts@medgb.jnj.com



Our commitment to you: the highest standard of quality, proven safe and effective clinical outcomes and differentiated service, everyday.

MENTOR is committed to supporting and partnering with you and we provide you the broadest range of products and services.

- We are dedicated to partnering and providing the expertise needed to help you differentiate your practice.
- Mentor Professional Education provides you the means to achieve excellence through better patient outcomes.
- From our customer service to our policies and procedure, we aim to exceed your expectations in working with Mentor.

MENTOR Core Study: The latest data.

LOWEST	LOW	NO
Reported incidence of the following complications at 10 years ^{1*}	Complication risk in over 9000 SILTEX Texture Implants ²	Silicone globules left behind ³
3.8% Capsular contracture in primary augmentation ¹	0% Reported double capsule ¹	0 out of 5 Capsules had silicone globules in a published study ³
1.5% Rotation in primary augmentation ¹		

*Not a head to head comparison. Based on 10 year core clinical study reported data of CPG™ Gel Breast Implants, NATRELLE™ 410 TruForm™ 3 Gel Breast Implants, NATRELLE™ Round TruForm™ 1 Gel Breast Implants, and MemoryGel™ Breast Implants. These core clinical studies DID NOT include NATRELLE™ Inspira Round TruForm™ 1 Gel Breast Implants, NATRELLE™ 410 TruForm™ 2 Gel Breast Implants, or NATRELLE™ 510 Dual Gel Anatomical Breast Implants. Only clinical studies with 10 year data were considered.

1. Mentor Worldwide, LLC. MemoryShape® Post-Approval Cohort Study (formerly Contour Profile Gel Core Study) Annual Report. June 2014.

2. Mentor's Core and Continued Access studies for round and shaped breast implants.

3. Danino, A.M. Basmacioglu, P, Saito, S, et al. PRSJ, 2001; 108:2047-2052. All other companies' brand names mentioned herein are the trademarks of their respective owners.

Limited warranty plan for MENTOR® Silicone Gel-filled Breast Implants.*

**You can't put a price on peace of mind.
So neither do we.**

PatientCare Promise means extra peace of mind for your MENTOR® Silicone Gel-filled Breast Implants.

- Free automatic enrollment.
- Lifetime product replacement in case of confirmed rupture.
- In the case of confirmed rupture, up to €1,000 of financial assistance, for a period of 10 years from the date of its implantation.
- Free product replacement in the case of confirmed capsular contracture Baker grade III and IV in augmentation surgery, for a period of 10 years from the date of its implantation.

MemoryGel™ Breast Implants.

MemoryGel™ Breast Implants Product Scale:

Product Profile	Cohesive I™	Cohesive II™	Cohesive III™
Round Moderate	●	●	
Round Moderate Plus	●	●	
Round High	●	●	
Round Ultra High	●		
CPG™ Gel Breast Implant			●
Round BECKER™ 25 Expander/Breast Implant	●		
Round BECKER™ 50 Expander/Breast Implant	●		
CONTOUR PROFILE™ BECKER™ 35 Expander/Breast Implant		●	

All Cohesive, all the time.™

Advanced MemoryGel™ Breast Implants.
The more evolved implant.



All MENTOR® Silicone Gel-Filled Breast Implants contain gel that is cohesive, safe and aesthetically pleasing.

Cohesive I™

The standard cohesive level gel used in MENTOR® Breast Implants. This is the softest gel and has been preferred for years in textured and smooth Round Low, Moderate, Moderate Plus, High, and Ultra High Profile Gel Breast Implants as well as the BECKER™ Expander/Breast Implants.

Cohesive II™

A slightly firmer gel, for surgeons wanting a firmer feeling implant. This gel is used in textured Round Moderate, Moderate Plus, High Profile Gel Breast Implants and CONTOUR PROFILE™ BECKER™ Expander/Breast Implants.

Cohesive III™

Mentor's most cohesive gel, providing shape retention with a pleasing level of firmness for optimal aesthetic results. This gel is used in the MENTOR® CPG™ Gel Breast Implants.

Round

Siltex™ Round Gel Breast Implants

Cohesive I™

Cohesive II™

CPG

CPG™ Gel Breast Implants, Cohesive III™

Expanders

CONTOUR PROFILE™ Tissue Expanders

BECKER™ Expander/Breast Implants

Accessories

Other

Smooth Implants

Saline Products

Other Expanders

Sizers

Surgical Resources

CPX-CPG™ Stage 1 to Stage 2 process

BodyLogic™ – Round Implants by base width

Important Safety Information

Sizing Guide BodyLogic™

**SILTEX™ Round Breast Implants,
Moderate Profile, Cohesive I™**


Vol cc	Diameter cm	Proj cm	Order Catalog #	Restorable Gel Sizer*
100	8.8	2.5	354-1007	
125	9.3	2.8	354-1257	RSZ-7125
150	10.2	2.7	354-1507	RSZ-7150
175	10.7	2.8	354-1757	RSZ-7175
200	11.2	2.8	354-2007	RSZ-7200
225	11.4	3.0	354-2257	RSZ-7225
250	11.5	3.3	354-2507	RSZ-7250
275	12.4	3.4	354-2757	RSZ-7275
300	12.6	3.5	354-3007	RSZ-7300
325	12.9	3.6	354-3257	RSZ-7325
350	13.4	3.7	354-3507	RSZ-7350
375	13.4	3.8	354-3757	RSZ-7375
400	13.5	3.9	354-4007	RSZ-7400
450	13.9	4.1	354-4507	RSZ-7450
500	14.2	4.3	354-5007	RSZ-7500
550	14.8	4.4	354-5507	RSZ-7550
600	15.4	4.5	354-6007	RSZ-7600
700	16.8	4.3	354-7007	RSZ-7700
800	17.2	4.6	354-8007	RSZ-7800

*Restorable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

Round Gel Breast Implants, Cohesive I™

**SILTEX™ Round Breast Implants,
Moderate Plus Profile, Cohesive I™**

Vol cc	Diameter cm	Proj cm	Order Catalog #	Resterilizable Gel Sizer*
100	8.1	2.7	354-1001	
125	8.8	2.9	354-1251	RSZ-1251
150	9.4	3.0	354-1501	RSZ-1501
175	10.0	3.2	354-1751	RSZ-1751
200	10.5	3.3	354-2001	RSZ-2001
225	10.9	3.5	354-2251	RSZ-2251
250	11.3	3.6	354-2501	RSZ-2501
275	11.7	3.7	354-2751	RSZ-2751
300	12.0	3.7	354-3001	RSZ-3001
325	12.3	3.8	354-3251	RSZ-3251
350	12.6	3.8	354-3501	RSZ-3501
375	12.9	3.9	354-3751	RSZ-3751
400	13.2	4.0	354-4001	RSZ-4001
450	13.7	4.1	354-4501	RSZ-4501
500	14.1	4.2	354-5001	RSZ-5001
550	14.4	4.4	354-5501	RSZ-5501
600	14.7	4.5	354-6001	RSZ-6001
700	15.7	4.8	354-7001	RSZ-7001
800	16.6	5.0	354-8001	RSZ-8001

Round Gel Breast Implants, Cohesive I™

**SILTEX™ Round Breast Implants,
High Profile, Cohesive I™**

Vol cc	Diameter cm	Proj cm	Order Catalog #	Resterilizable Gel Sizer*
125	8.4	3.6	354-4125	RSZ-1254
150	8.9	3.8	354-4150	RSZ-1504
175	9.4	4.0	354-4175	RSZ-1754
200	9.9	4.1	354-4200	RSZ-2004
225	10.2	4.3	354-4225	RSZ-2254
250	10.5	4.5	354-4250	RSZ-2504
275	10.9	4.6	354-4275	RSZ-2754
300	11.1	4.7	354-4300	RSZ-3004
325	11.5	4.8	354-4325	RSZ-3254
350	11.7	4.9	354-4350	RSZ-3504
375	12.0	5.0	354-4375	RSZ-3754
400	12.3	5.1	354-4400	RSZ-4004
425	12.5	5.2	354-4425	RSZ-4254
450	12.7	5.2	354-4450	RSZ-4504
500	13.2	5.4	354-4500	RSZ-5004
550	13.5	5.6	354-4550	RSZ-5504
600	14.0	5.7	354-4600	RSZ-6004
650	14.3	5.8	354-4650	RSZ-6504
700	14.7	6.0	354-4700	RSZ-7004
800	15.4	6.3	354-4800	RSZ-8004

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

Round Gel Breast Implants, Cohesive I™

**SILTEX™ Round Breast Implants,
Ultra High Profile, Cohesive I™**

Vol cc	Diameter cm	Proj cm	Order Catalog #	Resterilizable Gel Sizer*
135	8.0	4.3	354-5135	RSZ-5135
160	8.4	4.5	354-5160	RSZ-5160
185	8.6	4.6	354-5185	RSZ-5185
215	8.9	4.7	354-5215	RSZ-5215
240	9.2	4.9	354-5240	RSZ-5240
270	9.5	5.0	354-5270	RSZ-5270
295	9.9	5.2	354-5295	RSZ-5295
320	10.0	5.3	354-5320	RSZ-5320
350	10.3	5.4	354-5350	RSZ-5350
375	10.5	5.5	354-5375	RSZ-5375
400	10.8	5.6	354-5400	RSZ-5400
430	11.1	5.8	354-5430	RSZ-5430
455	11.3	5.9	354-5455	RSZ-5455
480	11.6	6.0	354-5480	RSZ-5480
535	12.1	6.3	354-5535	RSZ-5535
590	12.6	6.5	354-5590	RSZ-5590

Round Gel Breast Implants, Cohesive II™

**SILTEX™ Round Breast Implants,
Moderate Profile, Cohesive II™**

Vol cc	Diameter cm	Proj cm	Order Catalog #	Resterilizable Gel Sizer*
100	9.1	2.6	354-1009	
125	9.3	2.8	354-1259	RSZ-7125
150	10.1	2.4	354-1509	RSZ-7150
175	10.5	2.9	354-1759	RSZ-7175
200	11.2	2.8	354-2009	RSZ-7200
225	11.5	3.0	354-2259	RSZ-7225
250	11.6	3.5	354-2509	RSZ-7250
275	12.6	3.1	354-2759	RSZ-7275
300	12.6	3.6	354-3009	RSZ-7300
325	13.1	3.6	354-3259	RSZ-7325
350	13.4	3.5	354-3509	RSZ-7350
375	13.4	3.9	354-3759	RSZ-7375
400	13.7	3.8	354-4009	RSZ-7400
450	13.7	4.6	354-4509	RSZ-7450
500	14.2	4.8	354-5009	RSZ-7500
550	14.5	4.5	354-5509	RSZ-7550
600	15.3	4.6	354-6009	RSZ-7600
700	16.3	5.0	354-7009	RSZ-7700
800	16.7	5.4	354-8009	RSZ-7800

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

Round Gel Breast Implants, Cohesive II™

**SILTEX™ Round Breast Implants,
Moderate Plus Profile, Cohesive II™**

Vol cc	Diameter cm	Proj cm	Order Catalog #	Resterilizable Gel Sizer*
125	8.8	2.9	324-5125	RSZ-1251
150	9.4	3.0	324-5150	RSZ-1501
175	10.0	3.2	324-5175	RSZ-1751
200	10.5	3.3	324-5200	RSZ-2001
225	10.9	3.5	324-5225	RSZ-2251
250	11.3	3.6	324-5250	RSZ-2501
275	11.7	3.7	324-5275	RSZ-2751
300	12.0	3.7	324-5300	RSZ-3001
325	12.3	3.8	324-5325	RSZ-3251
350	12.6	3.8	324-5350	RSZ-3501
375	12.9	3.9	324-5375	RSZ-3751
400	13.2	4.0	324-5400	RSZ-4001
450	13.7	4.1	324-5450	RSZ-4501
500	14.1	4.2	324-5500	RSZ-5001
550	14.4	4.4	324-5550	RSZ-5501
600	14.7	4.5	324-5600	RSZ-6001
700	15.7	4.8	324-5700	RSZ-7001
800	16.6	5.0	324-5800	RSZ-8001

Round Gel Breast Implants, Cohesive II™

**SILTEX™ Round Breast Implants,
High Profile, Cohesive II™**

Vol cc	Diameter cm	Proj cm	Order Catalog #	Resterilizable Gel Sizer*
125	8.4	3.6	324-4125	RSZ-1254
150	8.9	3.8	324-4150	RSZ-1504
175	9.4	4.0	324-4175	RSZ-1754
200	9.9	4.1	324-4200	RSZ-2004
225	10.2	4.3	324-4225	RSZ-2254
250	10.5	4.5	324-4250	RSZ-2504
275	10.9	4.6	324-4275	RSZ-2754
300	11.1	4.7	324-4300	RSZ-3004
325	11.5	4.8	324-4325	RSZ-3254
350	11.7	4.9	324-4350	RSZ-3504
375	12.0	5.0	324-4375	RSZ-3754
400	12.3	5.1	324-4400	RSZ-4004
425	12.5	5.2	324-4425	RSZ-4254
450	12.7	5.2	324-4450	RSZ-4504
500	13.2	5.4	324-4500	RSZ-5004
550	13.5	5.6	324-4550	RSZ-5504
600	14.0	5.7	324-4600	RSZ-6004
650	14.3	5.8	324-4650	RSZ-6504
700	14.7	6.0	324-4700	RSZ-7004
800	15.4	6.3	324-4800	RSZ-8004

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.


Not all units will conform exactly to the dimensions noted above.

CPG™ Breast Implant Naming Legend.

Cohesive Level	Implant Height	Implant Projection
III 3	Low 1	Moderate 1
	Medium 2	Moderate Plus 2
	Tall 3	High 3


Lower Pole Arc (LPA) =

Length measured along the lower pole of the implant running from the bottom of the implant (Point B) up to the midway point (Point A) including a 0,5 cm of tissue thickness.



**CPG™ Implants 311, Cohesive III™,
Low Height, Moderate Projection**


Vol cc	Width cm	Height cm	Proj cm	Lower Pole Arc A-B cm	Order Catalog #	Resterilizable Gel Sizer*
120	9.5	8.5	3.2	6.9	334-0951	RSZ-0095
140	10.0	8.9	3.4	7.2	334-1001	RSZ-0100
160	10.5	9.4	3.6	7.5	334-1051	RSZ-0105
180	11.0	9.8	3.7	7.9	334-1101	RSZ-0110
210	11.5	10.3	3.9	8.2	334-1151	RSZ-0115
235	12.0	10.7	4.1	8.5	334-1201	RSZ-0120
270	12.5	11.2	4.2	8.8	334-1251	RSZ-0125
300	13.0	11.6	4.4	9.2	334-1301	RSZ-0130
335	13.5	12.0	4.6	9.5	334-1351	RSZ-0135
375	14.0	12.5	4.7	9.8	334-1401	RSZ-0140
415	14.5	12.9	4.9	10.1	334-1451	RSZ-0145
460	15.0	13.4	5.1	10.5	334-1501	RSZ-0150
510	15.5	13.8	5.2	10.8	334-1551	RSZ-0155
560	16.0	14.3	5.4	11.1	334-1601	RSZ-0160
615	16.5	14.7	5.6	11.4	334-1651	RSZ-0165

**CPG™ Implants 321, Cohesive III™,
Medium Height, Moderate Projection**


Vol cc	Width cm	Height cm	Proj cm	Lower Pole Arc A-B cm	Order Catalog #	Resterilizable Gel Sizer*
120	9.0	8.5	3.3	7.0	354-0908	RSZ-0908
135	9.5	8.9	3.5	7.3	354-0958	RSZ-0958
155	10.0	9.4	3.7	7.7	354-1008	RSZ-1008
180	10.5	9.9	3.8	8.0	354-1058	RSZ-1058
215	11.0	10.3	3.9	8.2	354-1108	RSZ-1108
245	11.5	10.8	4.0	8.6	354-1158	RSZ-1158
280	12.0	11.3	4.2	8.9	354-1208	RSZ-1208
315	12.5	11.8	4.4	9.2	354-1258	RSZ-1258
355	13.0	12.2	4.6	9.6	354-1308	RSZ-1308
395	13.5	12.7	4.7	9.9	354-1358	RSZ-1358
440	14.0	13.2	4.9	10.3	354-1408	RSZ-1408
480	14.5	13.6	5.0	10.5	354-1458	RSZ-1458
530	15.0	14.1	5.2	10.9	354-1508	RSZ-1508
640	16.0	15.0	5.6	11.5	354-1608	RSZ-1608
775	17.0	16.0	5.9	12.2	354-1708	RSZ-1708


*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

**CPG™ Implants 331, Cohesive III™
Tall Height, Moderate Projection**


Vol cc	Width cm	Height cm	Proj cm	Lower Pole Arc A-B cm	Order Catalog #	Resterilizable Gel Sizer*
125	9.0	9.2	3.2	7.3	334-0903	RSZ-0903
150	9.5	9.7	3.3	7.7	334-0953	RSZ-0953
175	10.0	10.2	3.4	8.1	334-1003	RSZ-1003
200	10.5	10.7	3.5	8.4	334-1053	RSZ-1053
230	11.0	11.3	3.6	8.8	334-1103	RSZ-1103
265	11.5	11.8	3.7	9.1	334-1153	RSZ-1153
300	12.0	12.3	3.9	9.5	334-1203	RSZ-1203
340	12.5	12.8	4.0	9.9	334-1253	RSZ-1253
380	13.0	13.3	4.1	10.2	334-1303	RSZ-1303
425	13.5	13.8	4.3	10.6	334-1353	RSZ-1353
475	14.0	14.3	4.5	11.0	334-1403	RSZ-1403
530	14.5	14.8	4.6	11.3	334-1453	RSZ-1453
585	15.0	15.3	4.8	11.7	334-1503	RSZ-1503
645	15.5	15.9	5.0	12.1	334-1553	RSZ-1553

**CPG™ Implants 312, Cohesive III™
Low Height, Moderate Plus Projection**


Vol cc	Width cm	Height cm	Proj cm	Lower Pole Arc A-B cm	Order Catalog #	Resterilizable Gel Sizer*
125	9.0	8.0	3.8	7.2	334-0907	RSZ-0907
145	9.5	8.4	4.0	7.5	334-0957	RSZ-0957
170	10.0	8.8	4.2	7.9	334-1007	RSZ-1007
195	10.5	9.3	4.4	8.3	334-1057	RSZ-1057
225	11.0	9.7	4.7	8.6	334-1107	RSZ-1107
255	11.5	10.2	4.9	9.0	334-1157	RSZ-1157
290	12.0	10.6	5.1	9.3	334-1207	RSZ-1207
330	12.5	11.1	5.3	9.7	334-1257	RSZ-1257
370	13.0	11.5	5.5	10.0	334-1307	RSZ-1307
415	13.5	12.0	5.7	10.4	334-1357	RSZ-1357
465	14.0	12.4	5.9	10.8	334-1407	RSZ-1407
515	14.5	12.8	6.1	11.1	334-1457	RSZ-1457
570	15.0	13.3	6.3	11.5	334-1507	RSZ-1507
690	16.0	14.2	6.8	11.2	334-1607	RSZ-1607

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

**CPG™ Implants 322, Cohesive III™,
Medium Height,
Moderate Plus Projection**



Vol cc	Width cm	Height cm	Proj cm	Lower Pole Arc A-B cm	Order Catalog #	Resterilizable Gel Sizer*
140	9.0	8.5	3.8	7.6	334-0905	RSZ-0905
165	9.5	8.9	4.0	8.0	334-0955	RSZ-0955
195	10.0	9.4	4.2	8.4	334-1005	RSZ-1005
225	10.5	9.9	4.4	8.8	334-1055	RSZ-1055
255	11.0	10.3	4.7	9.2	334-1105	RSZ-1105
295	11.5	10.8	4.9	9.5	334-1155	RSZ-1155
330	12.0	11.3	5.1	9.9	334-1205	RSZ-1205
375	12.5	11.8	5.3	10.3	334-1255	RSZ-1255
420	13.0	12.2	5.5	10.7	334-1305	RSZ-1305
475	13.5	12.7	5.7	11.1	334-1355	RSZ-1355
525	14.0	13.2	5.9	11.4	334-1405	RSZ-1405
585	14.5	13.6	6.1	11.8	334-1455	RSZ-1455
650	15.0	14.1	6.3	12.2	334-1505	RSZ-1505

**CPG™ Implants 332, Cohesive III™,
Tall Height,
Moderate Plus Projection**



Vol cc	Width cm	Height cm	Proj cm	Arc Length cm	Order Catalog #	Resterilizable Gel Sizer*
145	9.0	9.4	3.8	8.1	334-0909	RSZ-0909
175	9.5	9.9	4.0	8.5	334-0959	RSZ-0959
205	10.0	10.4	4.2	8.9	334-1009	RSZ-1009
235	10.5	10.9	4.4	9.3	334-1059	RSZ-1059
270	11.0	11.5	4.7	9.7	334-1109	RSZ-1109
305	11.5	12.0	4.9	10.1	334-1159	RSZ-1159
350	12.0	12.5	5.1	10.5	334-1209	RSZ-1209
395	12.5	13.0	5.3	10.9	334-1259	RSZ-1259
445	13.0	13.5	5.5	11.3	334-1309	RSZ-1309
495	13.5	14.1	5.7	11.7	334-1359	RSZ-1359
555	14.0	14.6	5.9	12.1	334-1409	RSZ-1409
615	14.5	15.1	6.1	12.6	334-1459	RSZ-1459
680	15.0	15.6	6.3	13.0	334-1509	RSZ-1509

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

**CPG™ Implants 313, Cohesive III™,
Low Height, High Projection**


Vol cc	Width cm	Height cm	Proj cm	Lower Pole Arc A-B cm	Order Catalog #	Resterilizable Gel Sizer*
130	9.0	8.1	4.4	7.6	334-0906	RSZ-0906
155	9.5	8.6	4.5	7.8	334-0956	RSZ-0956
180	10.0	9.0	4.6	8.3	334-1006	RSZ-1006
210	10.5	9.5	4.8	8.7	334-1056	RSZ-1056
240	11.0	9.9	4.9	9.1	334-1106	RSZ-1106
270	11.5	10.3	5.0	9.5	334-1156	RSZ-1156
310	12.0	10.8	5.2	9.9	334-1206	RSZ-1206
350	12.5	11.2	5.4	10.2	334-1256	RSZ-1256
395	13.0	11.7	5.6	10.6	334-1306	RSZ-1306
440	13.5	12.2	5.8	11.0	334-1356	RSZ-1356
490	14.0	12.6	6.1	11.4	334-1406	RSZ-1406
545	14.5	13.1	6.3	11.7	334-1456	RSZ-1456
605	15.0	13.5	6.6	12.1	334-1506	RSZ-1506

**CPG™ Implants 323, Cohesive III™,
Medium Height, High Projection**


Vol cc	Width cm	Height cm	Proj cm	Lower Pole Arc A-B cm	Order Catalog #	Resterilizable Gel Sizer*
165	9.0	8.5	4.6	8.1	334-0902	RSZ-0902
195	9.5	8.9	4.8	8.5	334-0952	RSZ-0952
225	10.0	9.4	5.1	8.9	334-1002	RSZ-1002
260	10.5	9.9	5.3	9.3	334-1052	RSZ-1052
300	11.0	10.3	5.6	9.7	334-1102	RSZ-1102
345	11.5	10.8	5.8	10.1	334-1152	RSZ-1152
390	12.0	11.3	6.0	10.5	334-1202	RSZ-1202
440	12.5	11.8	6.2	10.9	334-1252	RSZ-1252
495	13.0	12.2	6.5	11.3	334-1302	RSZ-1302
555	13.5	12.7	6.7	11.8	334-1352	RSZ-1352
620	14.0	13.2	6.9	12.2	334-1402	RSZ-1402
685	14.5	13.6	7.1	12.6	334-1452	RSZ-1452

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

CPG™ Gel Breast Implants, Cohesive III™

CPG™ Implants 333, Cohesive III™
Tall Height, High Projection

Vol cc	Width cm	Height cm	Proj cm	Lower Pole Arc A-B cm	Order Catalog #	Restorable Gel Sizer*
180	9.0	9.4	4.6	8.4	334-0904	RSZ-4090
215	9.5	9.9	4.8	8.8	334-0954	RSZ-4095
250	10.0	10.4	5.1	9.3	334-1004	RSZ-4100
290	10.5	10.9	5.3	9.7	334-1054	RSZ-4105
330	11.0	11.5	5.6	10.1	334-1104	RSZ-4110
380	11.5	12.0	5.8	10.5	334-1154	RSZ-4115
430	12.0	12.5	6.0	11.0	334-1204	RSZ-4120
485	12.5	13.0	6.2	11.4	334-1254	RSZ-4125
545	13.0	13.5	6.5	11.8	334-1304	RSZ-4130
610	13.5	14.1	6.7	12.2	334-1354	RSZ-4135
680	14.0	14.6	6.9	12.7	334-1404	RSZ-4140
755	14.5	15.1	7.1	12.1	334-1454	RSZ-4145

*Restorable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

BufferZone™ Self-sealing patch and posterior Dacron® Patch.

- Protecting the chest wall whilst driving lower pole expansion to mimic the shape of the natural breast.

Improved stronger magnet.*

- Easier to locate injection dome for larger patients.

Enhanced smooth, flat injection dome and shell pliability.*

- Helps ease insertion and removal through incision.

Better more pliant proprietary BufferZone™ Patch.*

- Enhanced patient comfort.

Patented SILTEX® Texture.

- Prevents the pooling or sequestration of body fluids.

Implant matching.

- Matching range of definitive implants to be used after tissue expansion.



*Compared with the MENTOR® CPX2™ and CPX3™ Breast Tissue Expanders.

CONTOUR PROFILE™ Tissue Expanders

**CPX™ 4 Tissue Expander,
Low Height, Style 8100***

Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
250	11.4	8.1	6.1	354-8121
350	12.7	9.4	6.5	354-8122
450	14.0	10.2	7.1	354-8123
550	15.0	10.9	7.4	354-8124
650	15.7	11.2	7.9	354-8125
750	16.5	11.9	8.1	354-8126

**CPX™ 4 Tissue Expander,
Medium Height, Style 8200***

Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
275	10.7	9.3	6.2	354-8221
350	11.7	10.0	6.6	354-8222
450	12.7	10.8	7.0	354-8223
550	13.5	11.7	7.4	354-8224
650	14.6	12.6	7.6	354-8225
800	15.6	13.3	8.0	354-8226

CONTOUR PROFILE™ Tissue Expanders

**CPX™ 4 Tissue Expander,
Tall Height, Style 8300***

Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
250	10.1	10.7	5.6	354-8321
350	11.3	11.8	6.0	354-8322
450	12.3	12.9	6.5	354-8323
550	13.2	13.8	6.9	354-8324
650	14.0	14.6	7.3	354-8325
750	14.6	15.3	7.6	354-8326
850	15.4	15.9	7.9	354-8327

*Includes CENTERSCOPE™ Magnetic Detection Device.

PLEASE NOTE: CONTOUR PROFILE™ Tissue Expanders are not MRI compatible.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

CONTOUR PROFILE™ Tissue Expanders

CPX™ 4 with Suture Tabs Tissue Expander, Low Height, Style 9100*



Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
250	11.4	8.1	6.1	354-9121
350	12.7	9.4	6.5	354-9122
450	14.0	10.2	7.1	354-9123
550	15.0	10.9	7.4	354-9124
650	15.7	11.2	7.9	354-9125
750	16.5	11.9	8.1	354-9126

CPX™ 4 with Suture Tabs Tissue Expander, Medium Height, Style 9200*



Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
275	10.7	9.3	6.2	354-9221
350	11.7	10.0	6.6	354-9222
450	12.7	10.8	7.0	354-9223
550	13.5	11.7	7.4	354-9224
650	14.6	12.6	7.6	354-9225
800	15.6	13.3	8.0	354-9226

CONTOUR PROFILE™ Tissue Expanders

CPX™ 4 with Suture Tabs Tissue Expander, Tall Height, Style 9300*



Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
250	10.1	10.7	5.6	354-9321
350	11.3	11.8	6.0	354-9322
450	12.3	12.9	6.5	354-9323
550	13.2	13.8	6.9	354-9324
650	14.0	14.6	7.3	354-9325
750	14.6	15.3	7.6	354-9326
850	15.4	15.9	7.9	354-9327

*Includes CENTERSCOPE™ Magnetic Detection Device.

PLEASE NOTE: CONTOUR PROFILE™ Tissue Expanders are not MRI compatible.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

Becker™ Expander/Breast Implants

**SILTEX™ Round BECKER™ 25
Expander/Implants, Cohesive I™**

(25% Silicone Gel, Cohesive I,™ in outer lumen, 75% Saline in inner lumen)



Vol cc	Gel Vol cc	Temporary Over Expansion Volumes		Total Saline cc	Total Gel Saline cc	Diameter cm	Proj cm	Order Catalog #
		Max Saline cc	Total Gel Saline cc					
150	40	185	225	85-150	125-190	10.7	2.4	354-1500
200	50	250	300	125-200	175-250	11.6	2.7	354-2000
250	60	315	375	165-255	225-315	12.1	3.0	354-2500
300	75	375	450	200-300	275-375	12.9	3.0	354-3000
350	90	435	525	235-350	325-440	13.5	3.0	354-3500
400	100	500	600	275-400	375-500	13.7	3.4	354-4000
500	125	625	750	350-500	475-625	14.7	3.7	354-5000
600	150	750	900	425-600	575-750	15.6	4.2	354-6000
700	175	875	1050	500-700	675-875	16.8	4.3	354-7000
800	200	1000	1200	575-800	775-1000	17.2	4.6	354-8000

Becker™ Expander/Breast Implants

**SILTEX™ Round BECKER™ 50
Expander/Implants, Cohesive I™**

(50% Silicone Gel, Cohesive I,™ in outer lumen, 50% Saline in inner lumen)



Vol cc	Gel Vol cc	Total Saline cc	Total Gel Saline cc	Diameter cm	Proj cm	Order Catalog #
300	150	150-200	300-350	12.9	3.0	354-1515
400	200	200-300	400-500	13.9	3.4	354-2020
500	250	250-350	500-600	14.7	3.7	354-2525
600	300	300-425	600-725	15.6	4.2	354-3030
700	350	350-500	700-850	16.0	5.0	354-3535

Note: Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

SILTEX™ CONTOUR PROFILE™ BECKER™ 35 Expander/Implants, Cohesive II™

(35% Silicone Gel, Cohesive II™, in outer lumen, 65% Saline in inner lumen)



		Temporary Over Expansion Volumes								
Vol cc	Gel Vol cc	Max Saline cc	Total Gel Saline cc	Total Saline cc	Total Gel Saline cc	Width cm	Height cm	Proj cm	Order Catalog #	
145	50	120	170	85-95	135-145	9.5	8.3	4.2	324-0955	
195	70	155	225	110-125	180-195	10.5	9.1	4.6	324-1055	
255	90	205	295	145-165	235-255	11.5	10.0	5.1	324-1155	
290	100	235	335	170-190	270-290	12.0	10.5	5.3	324-1205	
325	110	270	380	190-215	300-325	12.5	10.9	5.5	324-1255	
365	125	300	425	215-240	340-365	13.0	11.3	5.7	324-1305	
400	130	335	465	240-270	370-400	13.5	11.8	5.9	324-1355	
460	160	375	535	270-300	430-460	14.0	12.2	6.2	324-1405	
565	195	460	655	330-370	525-565	15.0	13.1	6.6	324-1505	
685	240	555	795	400-445	640-685	16.0	13.9	7.0	324-1605	

Injection Port Detector

Description

CENTERSCOPE™ Magnetic Detection Device

Order
Catalog #

350-4402

Injection Dome

Description

Dome Pack Set

Order
Catalog #

350-DOMPK

Can be used with SPECTRUM™ Post-Operatively Adjustable Saline Breast Implant, Round BECKER™ Expander/Breast Implant, and Smooth Tissue Expanders.

Replacement Butterfly Needle

Description

Winged Infusion Set

Order
Catalog #

7B3050

Note: Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

Round Gel Breast Implants, Cohesive I™

Smooth Round Breast Implants, Moderate Profile, Cohesive I™



Vol cc	Diameter cm	Proj cm	Order Catalog #	Resterilizable Gel Sizer*
125	10.0	2.2	350-7125 BC	RSZ-7125
150	10.3	2.3	350-7150 BC	RSZ-7150
175	11.2	2.4	350-7175 BC	RSZ-7175
200	11.7	2.5	350-7200 BC	RSZ-7200
225	12.2	2.6	350-7225 BC	RSZ-7225
250	12.3	2.8	350-7250 BC	RSZ-7250
275	13.2	2.9	350-7275 BC	RSZ-7275
300	13.5	3.0	350-7300 BC	RSZ-7300
325	13.9	3.0	350-7325 BC	RSZ-7325
350	14.2	3.1	350-7350 BC	RSZ-7350
375	14.4	3.2	350-7375 BC	RSZ-7375
400	14.5	3.2	350-7400 BC	RSZ-7400
450	14.9	3.4	350-7450 BC	RSZ-7450
500	15.2	3.6	350-7500 BC	RSZ-7500
550	15.9	3.6	350-7550 BC	RSZ-7550
600	16.5	3.7	350-7600 BC	RSZ-7600
700	17.4	3.9	350-7700 BC	RSZ-7700
800	18.2	4.1	350-7800 BC	RSZ-7800

Round Gel Breast Implants, Cohesive I™

Smooth Round Breast Implants, Moderate Plus Profile, Cohesive I™



Vol cc	Diameter cm	Proj cm	Order Catalog #	Resterilizable Gel Sizer*
125	8.9	2.8	350-1251 BC	RSZ-1251
150	9.5	2.9	350-1501 BC	RSZ-1501
175	10.0	3.1	350-1751 BC	RSZ-1751
200	10.5	3.2	350-2001 BC	RSZ-2001
225	10.9	3.3	350-2251 BC	RSZ-2251
250	11.3	3.4	350-2501 BC	RSZ-2501
275	11.7	3.5	350-2751 BC	RSZ-2751
300	12.0	3.6	350-3001 BC	RSZ-3001
325	12.3	3.8	350-3251 BC	RSZ-3251
350	12.5	3.9	350-3501 BC	RSZ-3501
375	12.8	4.0	350-3751 BC	RSZ-3751
400	13.1	4.0	350-4001 BC	RSZ-4001
450	13.6	4.2	350-4501 BC	RSZ-4501
500	14.1	4.3	350-5001 BC	RSZ-5001
550	14.6	4.5	350-5501 BC	RSZ-5501
600	15.0	4.6	350-6001 BC	RSZ-6001
700	15.8	4.9	350-7001 BC	RSZ-7001
800	16.5	5.1	350-8001 BC	RSZ-8001

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

Round Gel Breast Implants, Cohesive I™

Round Gel Breast Implants, Cohesive I™

Smooth Round Breast Implants, High Profile, Cohesive I™



Vol cc	Diameter cm	Proj cm	Order Catalog #	Resterilizable Gel Sizer*
125	8.3	3.5	350-1254 BC	RSZ-1254
150	8.8	3.7	350-1504 BC	RSZ-1504
175	9.3	3.9	350-1754 BC	RSZ-1754
200	9.7	4.0	350-2004 BC	RSZ-2004
225	10.1	4.2	350-2254 BC	RSZ-2254
250	10.5	4.3	350-2504 BC	RSZ-2504
275	10.8	4.4	350-2754 BC	RSZ-2754
300	11.1	4.5	350-3004 BC	RSZ-3004
325	11.4	4.6	350-3254 BC	RSZ-3254
350	11.7	4.8	350-3504 BC	RSZ-3504
375	12.0	4.8	350-3754 BC	RSZ-3754
400	12.2	5.0	350-4004 BC	RSZ-4004
425	12.5	5.0	350-4254 BC	RSZ-4254
450	12.8	5.1	350-4504 BC	RSZ-4504
500	13.2	5.3	350-5004 BC	RSZ-5004
550	13.6	5.5	350-5504 BC	RSZ-5504
600	14.0	5.6	350-6004 BC	RSZ-6004
650	14.4	5.7	350-6504 BC	RSZ-6504
700	14.8	5.8	350-7004 BC	RSZ-7004
800	15.5	6.0	350-8004 BC	RSZ-8004

Smooth Round Breast Implants, Ultra High Profile, Cohesive I™



Vol cc	Diameter cm	Proj cm	Order Catalog #	Resterilizable Gel Sizer*
135	7.8	4.1	350-5135 BC	RSZ-5135
160	8.2	4.3	350-5160 BC	RSZ-5160
185	8.4	4.4	350-5185 BC	RSZ-5185
215	8.7	4.5	350-5215 BC	RSZ-5215
240	9.0	4.7	350-5240 BC	RSZ-5240
270	9.3	4.8	350-5270 BC	RSZ-5270
295	9.7	5.0	350-5295 BC	RSZ-5295
320	9.8	5.1	350-5320 BC	RSZ-5320
350	10.1	5.2	350-5350 BC	RSZ-5350
375	10.4	5.3	350-5375 BC	RSZ-5375
400	10.6	5.4	350-5400 BC	RSZ-5400
430	10.9	5.6	350-5430 BC	RSZ-5430
455	11.2	5.7	350-5455 BC	RSZ-5455
480	11.4	5.8	350-5480 BC	RSZ-5480
535	12.0	6.1	350-5535 BC	RSZ-5535
590	12.5	6.3	350-5590 BC	RSZ-5590

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

Round Saline Breast Implants

Round Saline Breast Implants

Smooth Round Breast Implants, Moderate Plus Profile



Nominal Fill			Maximum Fill			Order Catalog #	Sterile Saline Sizert
Vol cc	Diameter cm	Proj cm	Vol cc	Diameter cm	Proj cm		
175	9.5	3.5	210	9.3	4.2	350-2175	351-2175SZ
200	10.0	3.7	240	9.7	4.4	350-2200	351-2200SZ
225	10.4	3.8	270	10.1	4.6	350-2225	351-2225SZ
250	10.8	4.0	300	10.5	4.7	350-2250	351-2250SZ
275	11.0	4.1	330	10.8	4.9	350-2275	351-2275SZ
300	11.5	4.3	360	11.2	5.0	350-2300	351-2300SZ
325	11.9	4.4	390	11.5	5.2	350-2325	351-2325SZ
350	12.1	4.5	420	11.7	5.3	350-2350	351-2350SZ
375	12.3	4.6	450	12.0	5.4	350-2375	351-2375SZ
400	12.6	4.7	480	12.3	5.5	350-2400	351-2400SZ
425	12.9	4.8	510	12.5	5.6	350-2425	351-2425SZ
450	13.0	4.9	540	12.8	5.7	350-2450	351-2450SZ
475	13.3	5.0	570	13.0	5.8	350-2475	351-2475SZ
500	13.6	5.1	600	13.2	5.9	350-2500	351-2500SZ
550	14.0	5.3	660	13.7	6.1	350-2550	351-2550SZ
600	14.5	5.5	720	14.1	6.3	350-2600	351-2600SZ
650	14.8	5.6	780	14.4	6.5	350-2650	351-2650SZ
700	15.2	5.8	840	14.8	6.6	350-2700	351-2700SZ
750	15.6	5.9	900	15.1	6.8	350-2750	351-2750SZ
800	16.1	6.0	960	15.5	7.0	350-2800	351-2800SZ
850	16.4	6.3	1020			350-2850	351-2850SZ
900	16.6	6.5	1080			350-2900	351-2900SZ

Smooth Round Breast Implants, High Profile



Nominal Fill			Maximum Fill			Order Catalog #	Sterile Saline Sizert
Vol cc	Diameter cm	Proj cm	Vol cc	Diameter cm	Proj cm		
170	9.0	4.0	200	8.8	4.9	350-3170	
190	9.3	4.1	225	9.1	5.1	350-3190	351-3190SZ
210	9.6	4.2	250	9.5	5.3	350-3210	
230	10.0	4.3	275	9.8	5.5	350-3230	351-3230SZ
250	10.2	4.5	300	10.0	5.6	350-3250	
270	10.4	4.6	325	10.2	5.8	350-3270	351-3270SZ
290	10.8	4.7	350	10.5	5.9	350-3290	351-3290SZ
310	11.0	4.8	375	10.7	6.0	350-3310	351-3310SZ
330	11.3	4.8	400	11.0	6.2	350-3330	351-3330SZ
380	11.7	5.2	450	11.4	6.4	350-3380	351-3380SZ
420	12.0	5.4	500	11.7	6.7	350-3420	351-3420SZ
460	12.4	5.5	550	12.1	6.9	350-3460	
500	12.8	5.6	600	12.4	7.1	350-3500	351-3500SZ
560	13.4	5.7	675	13.1	7.3	350-3560	
630	13.8	5.9	750	13.4	7.4	350-3630	351-3630SZ

tSingle use only.

Note: Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

Round Saline Breast Implants

SILTEX™ Round Saline Breast Implants

Smooth Round Breast Implants, Moderate Profile



Nominal Fill			Maximum Fill			Order Catalog #	Sterile Saline Sizert
Vol cc	Diameter cm	Proj cm	Vol cc	Diameter cm	Proj cm		
125	9.5	3.0	150	9.1	3.4	350-1610	351-125SZ
150	10.0	3.1	175	9.6	3.7	350-1615	
175	10.6	3.3	200	10.0	3.9	350-1620	351-175SZ
200	11.0	3.4	225	10.4	3.9	350-1625	
225	11.5	3.5	250	10.9	4.0	350-1630	351-225SZ
250	11.9	3.6	275	11.2	4.1	350-1635	
275	12.3	3.7	300	11.7	4.2	350-1640	351-275SZ
300	12.6	3.7	325	12.1	4.3	350-1645	
325	13.0	3.8	375	12.3	4.5	350-1650	351-325SZ
350	13.3	3.9	400	12.6	4.7	350-1655	
375	13.6	4.0	425	12.8	4.8	350-1660	351-375SZ
425	14.2	4.1	475	13.5	4.9	350-1670	351-425SZ
475	14.8	4.2	525	14.2	4.9	350-1680	351-475SZ
525	15.0	4.2	575	14.5	5.1	350-1685	
575	15.0	4.5	625	14.8	5.1	350-1690	
625	15.2	4.6	700	14.9	5.6	350-1695	
700	15.6	4.9	775	14.9	5.8	350-1697	

SILTEX™ Round Breast Implants, Moderate Profile



			Maximum Fill	Order Catalog #	Sterile Saline Sizert
Vol cc	Diameter cm	Proj cm	Vol cc		
125	9.5	3.0	150	354-2610	351-125SZ
150	10.0	3.1	175	354-2615	
175	10.6	3.3	200	354-2620	351-175SZ
200	11.0	3.4	225	354-2625	
225	11.5	3.5	250	354-2630	351-225SZ
250	11.9	3.6	275	354-2635	
275	12.3	3.7	300	354-2640	351-275SZ
300	12.6	3.7	325	354-2645	
325	13.0	3.8	375	354-2650	351-325SZ
350	13.3	3.9	400	354-2655	
375	13.6	4.0	425	354-2660	351-375SZ
425	14.2	4.1	475	354-2670	351-425SZ
475	14.8	4.2	525	354-2680	351-475SZ

†Single use only.

Note: Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

SILTEX™ CONTOUR PROFILE™ Saline Breast Implants


SILTEX™ CONTOUR PROFILE™ Breast Implants, Moderate Profile



Vol cc	Width cm	Height cm	Proj cc	Maximum Fill Volume	Order Catalog #
175	10.2	8.5	4.3	200	354-2910
225	11.2	9.4	4.5	250	354-2911
275	12.2	10.3	4.5	300	354-2912
350	13.1	11.0	4.7	400	354-2913
425	13.8	11.6	5.2	475	354-2914
525	14.9	12.5	5.6	575	354-2915

SILTEX™ CONTOUR PROFILE™ Saline Breast Implants

SILTEX™ CONTOUR PROFILE™ Breast Implants, High Profile




Vol cc	Width cm	Height cm	Proj cm	Max Fill Vol cc	Order Catalog #	Sterile Saline Sizert
275	11.5	9.5	5.1	300	354-2711	351-0711SZ
350	12.3	10.5	5.3	400	354-2712	351-0712SZ
450	13.2	11.0	6.1	500	354-2713	351-0713SZ
550	14.0	11.9	6.4	600	354-2714	351-0714SZ
650	15.0	12.7	6.6	725	354-2715	351-0715SZ

†Single use only.
Note: Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

SPECTRUM™ Post-Operatively Adjustable Saline Breast Implants

Smooth Round SPECTRUM™ Implants



Temp Vol cc	Min Vol cc	Nominal Fill			Maximum Fill			Order Catalog #
		Final Vol cc	Diameter cm	Proj cm	Final Vol cc	Diameter cm	Proj cm	
125	105	125	9.4	2.7	150	9.3	3.2	350-1410
175	150	175	10.3	3.1	210	10.2	3.8	350-1420
225	190	225	10.8	3.5	270	10.9	4.3	350-1430
275	230	275	11.6	3.8	330	11.5	4.6	350-1440
325	275	325	12.4	4.0	390	12.1	4.9	350-1450
375	320	375	13.5	4.1	450	13.0	4.9	350-1460
425	360	425	14.0	4.2	510	13.7	4.8	350-1470
475	405	475	14.3	4.4	570	13.8	5.4	350-1480
525	450	525	15.0	4.4	630	14.4	5.6	350-1485
575	490	575	14.7	4.8	690	14.2	6.0	350-1490


SPECTRUM™ Post-Operatively Adjustable Saline Breast Implants

SILTEX™ Round SPECTRUM™ Implants



Temp Vol cc	Min Vol cc	Nominal Fill			Maximum Fill			Order Catalog #
		Final Vol cc	Diameter cm	Proj cm	Final Vol cc	Diameter cm	Proj cm	
125	105	125	9.4	3.0	150	9.2	3.5	354-2410M
175	150	175	10.5	3.1	210	10.4	3.9	354-2420M
225	190	225	11.3	3.3	270	11.1	4.2	354-2430M
275	235	275	12.0	3.6	330	11.4	4.4	354-2440M
325	275	325	12.9	3.7	390	12.6	4.5	354-2450M
375	320	375	13.4	3.9	450	13.1	4.8	354-2460M
425	360	425	14.0	4.2	510	13.6	5.0	354-2470M
475	400	475	14.5	4.5	570	14.2	5.7	354-2480M

SILTEX™ CONTOUR PROFILE™ Expander/Implants




Temp Vol cc	Min Vol cc	Final Vol cc	Nominal Fill			Maximum Fill				Order Catalog #
			Width cm	Height cm	Proj cm	Final Vol cc	Width cm	Height cm	Proj cm	
275	235	275	11.5	9.5	5.1	330	11.1	9.4	6.3	354-2511
350	300	350	12.3	10.5	5.3	420	11.9	10.1	6.9	354-2512
450	380	450	13.2	11.0	6.1	540	12.7	11.0	7.5	354-2513
550	470	550	14.0	11.9	6.4	660	13.7	11.9	7.9	354-2514
650	550	650	15.0	12.7	6.6	780	14.5	12.5	8.2	354-2515

Other Tissue Expanders


Other Tissue Expanders

Smooth Round Tissue Expanders



Device Volume	Diameter	Projection	Order Catalog #
400 cc	11.3 cm	5.9 cm	350-4305M
550 cc	13.4 cm	5.8 cm	350-4307M
700 cc	14.6 cm	6.7 cm	350-4304M
850 cc	15.0 cm	7.2 cm	350-4308M
1000 cc	15.5 cm	7.8 cm	350-4306M

Smooth Crescent Tissue Expanders



Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
50	9.0	5.7	2.7	350-4312M
100	11.0	6.9	3.7	350-4313M
200	14.6	7.5	4.3	350-4314M
300	16.5	10.3	4.7	350-4315M
400	18.2	11.4	5.2	350-4316M
500	19.6	12.4	5.5	350-4317M

Smooth Elliptical Tissue Expanders



Device Vol cc	Length cm	Width cm	Nominal Proj cm	Proj cm	Max Volume cc	Order Catalog #
25	6.0	4.0	2.3	4.5	75	350-5301M
50	8.0	5.0	2.5	5.3	150	350-5302M
75	10.0	6.0	2.1	6.0	225	350-5303M
150	12.0	7.0	3.2	7.2	450	350-5304M
250	14.0	8.0	4.5	9.0	750	350-5305M
125	15.0	5.0	3.3	6.5	375	350-5306M

Smooth Rectangle Tissue Expanders



Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
100	6.6	5.1	5.6	350-4301M
250	10.4	7.3	5.6	350-4302M
400	10.6	9.3	6.6	350-4303M
550	13.7	7.1	7.7	350-4310M
700	15.0	8.0	7.9	350-4311M

Sterile Resterilizable Gel Breast Implant Sizers*

Sterile Resterilizable Gel Breast Implant Sizers*

Round Gel Breast Implant Sizers, Moderate Profile



Volume	Diameter	Projection	Order Catalog #
125 cc	10.0 cm	2.2 cm	RSZ-7125
150 cc	10.3 cm	2.3 cm	RSZ-7150
175 cc	11.2 cm	2.4 cm	RSZ-7175
200 cc	11.7 cm	2.5 cm	RSZ-7200
225 cc	12.2 cm	2.6 cm	RSZ-7225
250 cc	12.3 cm	2.8 cm	RSZ-7250
275 cc	13.2 cm	2.9 cm	RSZ-7275
300 cc	13.5 cm	3.0 cm	RSZ-7300
325 cc	13.9 cm	3.0 cm	RSZ-7325
350 cc	14.2 cm	3.1 cm	RSZ-7350
375 cc	14.4 cm	3.2 cm	RSZ-7375
400 cc	14.5 cm	3.2 cm	RSZ-7400
450 cc	14.9 cm	3.4 cm	RSZ-7450
500 cc	15.2 cm	3.6 cm	RSZ-7500
550 cc	15.9 cm	3.6 cm	RSZ-7550
600 cc	16.5 cm	3.7 cm	RSZ-7600
700 cc	17.4 cm	3.9 cm	RSZ-7700
800 cc	18.2 cm	4.1 cm	RSZ-7800

Round Gel Breast Implant Sizers, Moderate Plus Profile



Volume	Diameter	Projection	Order Catalog #
125 cc	8.9 cm	2.8 cm	RSZ-1251
150 cc	9.5 cm	2.9 cm	RSZ-1501
175 cc	10.0 cm	3.1 cm	RSZ-1751
200 cc	10.5 cm	3.2 cm	RSZ-2001
225 cc	10.9 cm	3.3 cm	RSZ-2251
250 cc	11.3 cm	3.4 cm	RSZ-2501
275 cc	11.7 cm	3.5 cm	RSZ-2751
300 cc	12.0 cm	3.6 cm	RSZ-3001
325 cc	12.3 cm	3.8 cm	RSZ-3251
350 cc	12.5 cm	3.9 cm	RSZ-3501
375 cc	12.8 cm	4.0 cm	RSZ-3751
400 cc	13.1 cm	4.0 cm	RSZ-4001
450 cc	13.6 cm	4.2 cm	RSZ-4501
500 cc	14.1 cm	4.3 cm	RSZ-5001
550 cc	14.6 cm	4.5 cm	RSZ-5501
600 cc	15.0 cm	4.6 cm	RSZ-6001
700 cc	15.8 cm	4.9 cm	RSZ-7001
800 cc	16.5 cm	5.1 cm	RSZ-8001

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

Sterile Resterilizable Gel Breast Implant Sizers*

Sterile Resterilizable Gel Breast Implant Sizers*

Round Gel Breast Implant Sizers, High Profile



Volume	Diameter	Projection	Order Catalog #
125 cc	8.3 cm	3.5 cm	RSZ-1254
150 cc	8.8 cm	3.7 cm	RSZ-1504
175 cc	9.3 cm	3.9 cm	RSZ-1754
200 cc	9.7 cm	4.0 cm	RSZ-2004
225 cc	10.1 cm	4.2 cm	RSZ-2254
250 cc	10.5 cm	4.3 cm	RSZ-2504
275 cc	10.8 cm	4.4 cm	RSZ-2754
300 cc	11.1 cm	4.5 cm	RSZ-3004
325 cc	11.4 cm	4.6 cm	RSZ-3254
350 cc	11.7 cm	4.8 cm	RSZ-3504
375 cc	12.0 cm	4.8 cm	RSZ-3754
400 cc	12.2 cm	5.0 cm	RSZ-4004
450 cc	12.8 cm	5.1 cm	RSZ-4504
500 cc	13.2 cm	5.3 cm	RSZ-5004
550 cc	13.6 cm	5.5 cm	RSZ-5504
600 cc	14.0 cm	5.6 cm	RSZ-6004
650 cc	14.4 cm	5.7 cm	RSZ-6504
700 cc	14.8 cm	5.8 cm	RSZ-7004
800 cc	15.5 cm	6.0 cm	RSZ-8004

Round Gel Breast Implant Sizers, Ultra High Profile



Volume	Diameter	Projection	Order Catalog #
135 cc	7.8 cm	4.1 cm	RSZ-5135
160 cc	8.2 cm	4.3 cm	RSZ-5160
185 cc	8.4 cm	4.4 cm	RSZ-5185
215 cc	8.7 cm	4.5 cm	RSZ-5215
240 cc	9.0 cm	4.7 cm	RSZ-5240
270 cc	9.3 cm	4.8 cm	RSZ-5270
295 cc	9.7 cm	5.0 cm	RSZ-5295
320 cc	9.8 cm	5.1 cm	RSZ-5320
350 cc	10.1 cm	5.2 cm	RSZ-5350
375 cc	10.4 cm	5.3 cm	RSZ-5375
400 cc	10.6 cm	5.4 cm	RSZ-5400
430 cc	10.9 cm	5.6 cm	RSZ-5430
455 cc	11.2 cm	5.7 cm	RSZ-5455
480 cc	11.4 cm	5.8 cm	RSZ-5480
535 cc	12.0 cm	6.1 cm	RSZ-5535
590 cc	12.5 cm	6.3 cm	RSZ-5590

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

CPG™ Implant Sizers 311
Cohesive III™, Low Height,
Moderate Projection


Vol cc	Width cm	Height cm	Proj cm	Lower Pole Arc A-B cm	Order Catalog #
120	9.5	8.5	3.2	6.9	334-0951
140	10.0	8.9	3.4	7.2	334-1001
160	10.5	9.4	3.6	7.5	334-1051
180	11.0	9.8	3.7	7.9	334-1101
210	11.5	10.3	3.9	8.2	334-1151
235	12.0	10.7	4.1	8.5	334-1201
270	12.5	11.2	4.2	8.8	334-1251
300	13.0	11.6	4.4	9.2	334-1301
335	13.5	12.0	4.6	9.5	334-1351
375	14.0	12.5	4.7	9.8	334-1401
415	14.5	12.9	4.9	10.1	334-1451
460	15.0	13.4	5.1	10.5	334-1501
510	15.5	13.8	5.2	10.8	334-1551
560	16.0	14.3	5.4	11.1	334-1601
615	16.5	14.7	5.6	11.4	334-1651

CPG™ Implant Sizers 321
Cohesive III™, Medium Height,
Moderate Projection


Vol cc	Width cm	Height cm	Proj cm	A cm	Arc Length cm	Order Catalog #
120	9.0	8.5	3.3	3.1	4.3	RSZ-0908
135	9.5	8.9	3.5	3.2	4.5	RSZ-0958
155	10.0	9.4	3.7	3.4	4.7	RSZ-1008
180	10.5	9.9	3.8	3.6	5.0	RSZ-1058
215	11.0	10.3	3.9	3.6	5.1	RSZ-1108
245	11.5	10.8	4.0	3.8	5.3	RSZ-1158
280	12.0	11.3	4.2	4.0	5.6	RSZ-1208
315	12.5	11.8	4.4	4.1	5.8	RSZ-1258
355	13.0	12.2	4.6	4.3	6.0	RSZ-1308
395	13.5	12.7	4.7	4.5	6.3	RSZ-1358
440	14.0	13.2	4.9	4.6	6.5	RSZ-1408
480	14.5	13.6	5.0	5.0	6.9	RSZ-1458
530	15.0	14.1	5.2	5.1	7.2	RSZ-1508
640	16.0	15.0	5.6	5.5	7.6	RSZ-1608
775	17.0	16.0	5.9	5.8	8.1	RSZ-1708

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.


Not all units will conform exactly to the dimensions noted above.

CPG™ Implant Sizers 331
Cohesive III™, Tall Height,
Moderate Projection



Vol cc	Width cm	Height cm	Proj cm	A cm	Arc Length cm	Order Catalog #
125	9.0	9.2	3.2	2.4	3.9	RSZ-0903
150	9.5	9.7	3.3	2.6	4.1	RSZ-0953
175	10.0	10.2	3.4	2.7	4.3	RSZ-1003
200	10.5	10.7	3.5	2.8	4.6	RSZ-1053
230	11.0	11.3	3.6	3.0	4.8	RSZ-1103
265	11.5	11.8	3.7	3.1	5.0	RSZ-1153
300	12.0	12.3	3.9	3.2	5.2	RSZ-1203
340	12.5	12.8	4.0	3.4	5.4	RSZ-1253
380	13.0	13.3	4.1	3.5	5.6	RSZ-1303
425	13.5	13.8	4.3	3.6	5.9	RSZ-1353
475	14.0	14.3	4.5	3.8	6.1	RSZ-1403
530	14.5	14.8	4.6	3.9	6.3	RSZ-1453
585	15.0	15.3	4.8	4.0	6.5	RSZ-1503
645	15.5	15.9	5.0	4.2	6.7	RSZ-1553

CPG™ Implant Sizers 312
Cohesive III™, Low Height,
Moderate Plus Projection



Vol cc	Width cm	Height cm	Proj cm	A cm	Arc Length cm	Order Catalog #
125	9.0	8.0	3.8	2.5	4.6	RSZ-0907
145	9.5	8.4	4.0	2.6	4.8	RSZ-0957
170	10.0	8.8	4.2	2.7	5.1	RSZ-1007
195	10.5	9.3	4.4	2.9	5.4	RSZ-1057
225	11.0	9.7	4.7	3.0	5.6	RSZ-1107
255	11.5	10.2	4.9	3.1	5.9	RSZ-1157
290	12.0	10.6	5.1	3.3	6.1	RSZ-1207
330	12.5	11.1	5.3	3.4	6.4	RSZ-1257
370	13.0	11.5	5.5	3.6	6.6	RSZ-1307
415	13.5	12.0	5.7	3.7	6.9	RSZ-1357
465	14.0	12.4	5.9	3.8	7.1	RSZ-1407
515	14.5	12.8	6.1	4.0	7.4	RSZ-1457
570	15.0	13.3	6.3	4.1	7.7	RSZ-1507
690	16.0	14.2	6.8	4.4	8.2	RSZ-1607

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

CPG™ Implant Sizers 322
Cohesive III™, Medium Height,
Moderate Plus Projection



Vol cc	Width cm	Height cm	Proj cm	A cm	Arc Length cm	Order Catalog #
140	9.0	8.5	3.8	2.9	4.7	RSZ-0905
165	9.5	8.9	4.0	3.0	4.9	RSZ-0955
195	10.0	9.4	4.2	3.2	5.2	RSZ-1005
225	10.5	9.9	4.4	3.3	5.5	RSZ-1055
255	11.0	10.3	4.7	3.5	5.7	RSZ-1105
295	11.5	10.8	4.9	3.7	6.0	RSZ-1155
330	12.0	11.3	5.1	3.8	6.3	RSZ-1205
375	12.5	11.8	5.3	4.0	6.5	RSZ-1255
420	13.0	12.2	5.5	4.1	6.8	RSZ-1305
475	13.5	12.7	5.7	4.3	7.0	RSZ-1355
525	14.0	13.2	5.9	4.5	7.3	RSZ-1405
585	14.5	13.6	6.1	4.6	7.6	RSZ-1455
650	15.0	14.1	6.3	4.8	7.8	RSZ-1505

CPG™ Implant Sizers 332
Cohesive III™, Tall Height,
Moderate Plus Projection



Vol cc	Width cm	Height cm	Proj cm	A cm	Arc Length cm	Order Catalog #
145	9.0	9.4	3.8	2.8	4.6	RSZ-0909
175	9.5	9.9	4.0	3.0	4.9	RSZ-0959
205	10.0	10.4	4.2	3.1	5.1	RSZ-1009
235	10.5	10.9	4.4	3.3	5.4	RSZ-1059
270	11.0	11.5	4.7	3.5	5.7	RSZ-1109
305	11.5	12.0	4.9	3.6	5.9	RSZ-1159
350	12.0	12.5	5.1	3.8	6.2	RSZ-1209
395	12.5	13.0	5.3	3.9	6.4	RSZ-1259
445	13.0	13.5	5.5	4.1	6.7	RSZ-1309
495	13.5	14.1	5.7	4.2	7.0	RSZ-1359
555	14.0	14.6	5.9	4.4	7.2	RSZ-1409
615	14.5	15.1	6.1	4.6	7.5	RSZ-1459
680	15.0	15.6	6.3	4.7	7.7	RSZ-1509

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

CPG™ Implant Sizers 313
Cohesive III™, Low Height,
High Projection


Vol cc	Width cm	Height cm	Proj cm	A cm	Arc Length cm	Order Catalog #
130	9.0	8.1	4.4	2.6	5.2	RSZ-0906
155	9.5	8.6	4.5	2.7	5.4	RSZ-0956
180	10.0	9.0	4.6	2.9	5.7	RSZ-1006
210	10.5	9.5	4.8	3.0	6.0	RSZ-1056
240	11.0	9.9	4.9	3.2	6.3	RSZ-1106
270	11.5	10.3	5.0	3.3	6.6	RSZ-1156
310	12.0	10.8	5.2	3.5	6.9	RSZ-1206
350	12.5	11.2	5.4	3.6	7.2	RSZ-1256
395	13.0	11.7	5.6	3.8	7.4	RSZ-1306
440	13.5	12.2	5.8	3.9	7.7	RSZ-1356
490	14.0	12.6	6.1	4.0	8.0	RSZ-1406
545	14.5	13.1	6.3	4.2	8.3	RSZ-1456
605	15.0	13.5	6.6	4.3	8.6	RSZ-1506

CPG™ Implant Sizers 323
Cohesive III™, Medium Height,
High Projection


Vol cc	Width cm	Height cm	Proj cm	A cm	Arc Length cm	Order Catalog #
165	9.0	8.5	4.6	3.0	5.3	RSZ-0902
195	9.5	8.9	4.8	3.2	5.6	RSZ-0952
225	10.0	9.4	5.1	3.4	5.9	RSZ-1002
260	10.5	9.9	5.3	3.5	6.2	RSZ-1052
300	11.0	10.3	5.6	3.7	6.5	RSZ-1102
345	11.5	10.8	5.8	3.9	6.8	RSZ-1152
390	12.0	11.3	6.0	4.0	7.1	RSZ-1202
440	12.5	11.8	6.2	4.2	7.4	RSZ-1252
495	13.0	12.2	6.5	4.4	7.7	RSZ-1302
555	13.5	12.7	6.7	4.5	8.0	RSZ-1352
620	14.0	13.2	6.9	4.7	8.3	RSZ-1402
685	14.5	13.6	7.1	4.9	8.6	RSZ-1452

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

CPG™ Gel Breast Implant Sizers*

Sterile Saline Breast Implant Sizers†

CPG™ Implant Sizers 333 Cohesive III™, Tall Height, High Projection



Vol cc	Width cm	Height cm	Proj cm	A cm	Arc Length cm	Order Catalog #
145	9.0	9.4	3.8	2.8	4.6	RSZ-0909
175	9.5	9.9	4.0	3.0	4.9	RSZ-0959
205	10.0	10.4	4.2	3.1	5.1	RSZ-1009
235	10.5	10.9	4.4	3.3	5.4	RSZ-1059
270	11.0	11.5	4.7	3.5	5.7	RSZ-1109
305	11.5	12.0	4.9	3.6	5.9	RSZ-1159
350	12.0	12.5	5.1	3.8	6.2	RSZ-1209
395	12.5	13.0	5.3	3.9	6.4	RSZ-1259
445	13.0	13.5	5.5	4.1	6.7	RSZ-1309
495	13.5	14.1	5.7	4.2	7.0	RSZ-1359
555	14.0	14.6	5.9	4.4	7.2	RSZ-1409
615	14.5	15.1	6.1	4.6	7.5	RSZ-1459
680	15.0	15.6	6.3	4.7	7.7	RSZ-1509

Round Saline Breast Implant Sizers, Moderate Profile



Nominal Fill			Maximum Fill			Order Catalog #
Vol cc	Diameter cm	Proj cm	Vol cc	Diameter cm	Proj cm	
125	9.5	3.0	150	9.1	3.4	351-125SZ
175	10.6	3.3	200	10.0	3.9	351-175SZ
225	11.5	3.5	250	10.9	4.0	351-225SZ
275	12.3	3.7	300	11.7	4.2	351-275SZ
325	13.0	3.8	375	12.3	4.5	351-325SZ
375	13.6	4.0	425	12.8	4.8	351-375SZ
425	14.2	4.1	475	13.5	4.9	351-425SZ
475	14.8	4.2	525	14.2	4.9	351-475SZ

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

†Single use only.

Note: Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

Sterile Saline Breast Implant Sizers†

Round Saline Breast Implant Sizers, Moderate Plus Profile



Vol cc	Width cm	Height cm	Proj cm	Vol cc	Sterile Saline Sizer	Order Catalog #
175	9.5	3.5	210	9.3	4.2	351-2175SZ
200	10.0	3.7	240	9.7	4.4	351-2200SZ
225	10.4	3.8	270	10.1	4.6	351-2225SZ
250	10.8	4.0	300	10.5	4.7	351-2250SZ
275	11.0	4.1	330	10.8	4.9	351-2275SZ
300	11.5	4.3	360	11.2	5.0	351-2300SZ
325	11.9	4.4	390	11.5	5.2	351-2325SZ
350	12.1	4.5	420	11.7	5.3	351-2350SZ
375	12.3	4.6	450	12.0	5.4	351-2375SZ
400	12.6	4.7	480	12.3	5.5	351-2400SZ
425	12.9	4.8	510	12.5	5.6	351-2425SZ
450	13.0	4.9	540	12.8	5.7	351-2450SZ
475	13.3	5.0	570	13.0	5.8	351-2475SZ
500	13.6	5.1	600	13.2	5.9	351-2500SZ
550	14.0	5.3	660	13.7	6.1	351-2550SZ
600	14.5	5.5	720	14.1	6.3	351-2600SZ
650	14.8	5.6	780	14.4	6.5	351-2650SZ
700	15.2	5.8	840	14.8	6.6	351-2700SZ
750	15.6	5.9	900	15.1	6.8	351-2750SZ
800	16.1	6.0	960	15.5	7.0	351-2800SZ

Sterile Saline Breast Implant Sizers†

Round Saline Breast Implant Sizers, High Profile



Vol cc	Diameter cm	Proj cm	Vol cc	Diameter cm	Proj cm	Order Catalog #
190	9.3	4.1	225	9.1	5.1	351-3190SZ
230	10.0	4.3	275	9.8	5.5	351-3230SZ
270	10.4	4.6	325	10.2	5.8	351-3270SZ
290	10.8	4.7	350	10.5	5.9	351-3290SZ
310	11.0	4.8	375	10.7	6.0	351-3310SZ
330	11.3	4.8	400	11.0	6.2	351-3330SZ
380	11.7	5.2	450	11.4	6.4	351-3380SZ
420	12.0	5.4	500	11.7	6.7	351-3420SZ
500	12.8	5.6	600	12.4	7.1	351-3500SZ
630	13.8	5.9	750	13.4	7.4	351-3630SZ

†Single use only.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above.

Sterile Saline Implant Sizers†

CONTOUR PROFILE™ Saline
Breast Implant Sizers, High Profile



Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
275 + 25	11.5	9.5	5.1	351-0711SZ
350 + 50	12.3	10.5	5.3	351-0712SZ
450 + 50	13.2	11.0	6.1	351-0713SZ
550 + 50	14.0	11.9	6.4	351-0714SZ
650 + 75	15.0	12.7	6.6	351-0715SZ


†Single use only.
Note: Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

BodyLogic™ System Preoperative Worksheet

Patient Name: _____

Age: _____ Weight: _____ Height: _____ PreOp Bra Size: _____


Desired Outcome:



Skin Envelope:

Right: ☐ Tight ☐ Moderate ☐ Loose ☐ with Pits

Left: ☐ Tight ☐ Moderate ☐ Loose ☐ with Pits



Patient Dimensions:

Right	Left	Right	Left
Sternal Notch to Nipple (SN-N): _____	Nipple to IMF (N-IMF): _____		
Clavicle to Nipple (C-N): _____	Nipple* to IMF (stretched): _____		
Breast Base Width (BBW): _____	Compare with Lower Pole Arc measurement when using CPG™ Breast Implants		
Breast Height (BH): _____	Upper Pole Pinch: _____		
Medial Pinch: _____	Areolar Diameter: _____		
Lateral Pinch: _____	Intermammary Distance: _____		

Observations:

Breast Shape: ☐ Tubular ☐ Round ☐ Conical ☐ Wide ☐ Narrow ☐ Other _____

Breast Exam: _____

Estimated Breast Volume: Right: _____ Left: _____ Difference: _____

Asymmetry: Chest Wall: _____ Breast: _____

Nipple Level Discrepancy: _____ IMF Level Discrepancy: _____

Surgery Planning:

Surgery Type: ☐ Augmentation ☐ Revision ☐ Reconstruction

Surgery Date: _____ Location: _____ Inframammary Fold: ☐ Raise ☐ Keep ☐ Lower

Preferred External Size: Volume: _____ Profile: _____ Breast Base Width: _____

Incision Approach: ☐ Inframammary ☐ Periareolar ☐ Axillary ☐ Other _____

Implant Placement: ☐ Subpectoral ☐ Subglandular ☐ Dual Plane _____

Implant Base Diameter: _____

Breast Base Width – 1/2 medial pinch – 1/2 lateral pinch

BodyLogic™ System Preoperative Worksheet

Product Selection:

Patient Name: _____

MENTOR® Breast Implants: ☐ Memory Gel™ Breast Implants ☐ CPG™ Gel Breast Implants ☐ Saline

☐ SPECTRUM™ Post-Operatively Adjustable Saline Breast Implant

☐ Tissue Expander ☐ Low ☐ Medium ☐ Tall

Implant Specs

Right Left

Catalog # _____

Base Width _____

Height _____

Profile _____

Lower Pole Arc _____

(CPG™ Implants only)

Gel Volume _____

Saline Nominal Fill Volume _____

Saline Final _____

Tissue Expander Specs

Right Left

Catalog # _____

Base Width _____

Height _____

Volume _____

Final Fill Volume _____

Please find BodyLogic worksheets at www.mentormedical.co.uk

SILTEX™ Round MemoryGel Breast Implants

Organised by base diameter

Cohesive I™

Diameter cm	Volume cm	Projection cm	Profile	Order Catalog#
8.0	135	4.3	Ultra High	354-5135
8.1	100	2.7	Moderate Plus	354-1001
8.4	125	3.6	High	354-4125
8.4	160	4.5	Ultra High	354-5160
8.6	185	4.6	Ultra High	354-5185
8.8	100	2.5	Moderate	354-1007
8.8	125	2.9	Moderate Plus	354-1251
8.9	150	3.8	High	354-4150
8.9	215	4.7	Ultra High	354-5215
9.2	240	4.9	Ultra High	354-5240
9.3	125	2.8	Moderate	354-1257
9.4	150	3.0	Moderate Plus	354-1501
9.4	175	4.0	High	354-4175
9.5	270	5.0	Ultra High	354-5270
9.9	200	4.1	High	354-4200
9.9	295	5.2	Ultra High	354-5295
10.0	175	3.2	Moderate Plus	354-1751
10.0	320	5.3	Ultra High	354-5420
10.2	150	2.7	Moderate	354-1507
10.2	225	4.3	High	354-4225
10.3	350	5.4	Ultra High	354-5350
10.5	200	3.3	Moderate Plus	354-2001
10.5	250	4.5	High	354-4250
10.5	375	5.5	Ultra High	354-5375
10.7	175	2.8	Moderate	354-1757
10.8	400	5.6	Ultra High	354-5400
10.9	225	3.5	Moderate Plus	354-2251
10.9	275	4.6	High	354-4275

SILTEX™ Round MemoryGel Breast Implants

Organised by base diameter

Cohesive I™

Diameter cm	Volume cm	Projection cm	Profile	Order Catalog#
11.1	300	4.7	High	354-4300
11.1	430	5.8	Ultra High	354-5430
11.2	200	2.8	Moderate	354-2007
11.3	250	3.6	Moderate Plus	354-2501
11.3	455	5.9	Ultra High	354-5455
11.4	225	3.0	Moderate	354-2257
11.5	250	3.3	Moderate	354-2507
11.5	325	4.8	High	354-4325
11.6	480	6.0	Ultra High	354-5480
11.7	275	3.7	Moderate Plus	354-2751
11.7	350	4.9	High	354-4350
12.0	300	3.7	Moderate Plus	354-3001
12.0	375	5.0	High	354-4375
12.1	535	6.3	Ultra High	354-5535
12.3	325	3.8	Moderate Plus	354-3251
12.3	400	5.1	High	354-4400
12.4	275	3.4	Moderate	354-2757
12.5	425	5.2	High	354-4425
12.6	300	3.5	Moderate	354-3007
12.6	350	3.8	Moderate Plus	354-3501
12.6	590	6.5	Ultra High	354-5590
12.7	450	5.2	High	354-4450
12.9	325	3.6	Moderate	354-3257
12.9	375	3.9	Moderate Plus	354-3751
13.2	400	4.0	Moderate Plus	354-4001
13.2	500	5.4	High	354-4500
13.4	350	3.7	Moderate	354-3507
13.4	375	3.8	Moderate	354-3757
13.5	400	3.9	Moderate	354-4007
13.5	550	5.6	High	354-4550

SILTEX™ Round MemoryGel Breast Implants

Organised by base diameter

Cohesive I™

Diameter cm	Volume cm	Projection cm	Profile	Order Catalog#
13.7	450	4.1	Moderate Plus	354-4501
13.7	450	4.1	Moderate Plus	354-4501
13.9	450	4.1	Moderate	354-4507
14.0	600	5.7	High	354-4600
14.1	500	4.2	Moderate Plus	354-5001
14.2	500	4.3	Moderate	354-5007
14.3	650	5.8	High	354-4650
14.4	550	4.4	Moderate Plus	354-5501
14.7	600	4.5	Moderate Plus	354-6001
14.7	700	6.0	High	354-4700
14.8	550	4.4	Moderate	354-5507
15.4	600	4.5	Moderate	354-6007
15.4	800	6.3	High	354-4800
15.7	700	4.8	Moderate Plus	354-7001
16.6	800	5.0	Moderate Plus	354-8001
16.8	700	4.3	Moderate	354-7007
17.2	800	4.6	Moderate	354-8007

SILTEX™ Round MemoryGel Breast Implants

Organised by base diameter

Cohesive II™

Diameter cm	Volume cm	Projection cm	Profile	Order Catalog#
8.4	125	3.6	High	324-4125
8.8	125	2.9	Moderate Plus	324-5125
8.9	150	3.8	High	324-4150
9.1	100	2.6	Moderate	354-1009
9.3	125	2.8	Moderate	354-1259
9.4	150	3.0	Moderate Plus	324-5150
9.4	175	4.0	High	324-4175
9.9	200	4.1	High	324-4200
10.0	175	3.2	Moderate Plus	324-5175
10.1	150	2.4	Moderate	354-1509
10.2	225	4.3	High	324-4225
10.5	175	2.9	Moderate	354-1759
10.5	200	3.3	Moderate Plus	324-5200
10.5	250	4.5	High	324-4250
10.9	225	3.5	Moderate Plus	324-5225
10.9	275	4.6	High	324-4275
11.1	300	4.7	High	324-4300
11.2	200	2.8	Moderate	354-2009
11.3	250	3.6	Moderate Plus	324-5250
11.5	225	3.0	Moderate	354-2259
11.5	325	4.8	High	324-4325
11.6	250	3.5	Moderate	354-2509
11.7	275	3.7	Moderate Plus	324-5275
11.7	350	4.9	High	324-4350
12.0	300	3.7	Moderate Plus	324-5300
12.0	375	5.0	High	324-4375
12.3	325	3.8	Moderate Plus	324-5325
12.3	400	5.1	High	24-4400
12.5	425	5.2	High	324-4425

SILTEX™ Round MemoryGel Breast Implants

Organised by base diameter.

Cohesive II™

Diameter cm	Volume cm	Projection cm	Profile	Order Catalog#
12.6	275	3.1	Moderate	354-2759
12.6	300	3.6	Moderate	354-3009
12.6	350	3.8	Moderate Plus	324-5350
12.7	450	5.2	High	324-4450
12.9	375	3.9	Moderate Plus	324-5375
13.1	325	3.6	Moderate	354-3259
13.2	400	4.0	Moderate Plus	324-5400
13.2	500	5.4	High	324-4500
13.4	350	3.5	Moderate	354-3509
13.4	375	3.9	Moderate	354-3759
13.5	550	5.6	High	324-4550
13.7	400	3.8	Moderate	354-4009
13.7	450	4.6	Moderate	354-4509
13.7	450	4.1	Moderate Plus	324-5450
14.0	600	5.7	High	324-4600
14.1	500	4.2	Moderate Plus	324-5500
14.2	500	4.8	Moderate	354-5009
14.3	650	5.8	High	324-4650
14.4	550	4.4	Moderate Plus	324-5550
14.5	550	4.5	Moderate	354-5509
14.7	600	4.5	Moderate Plus	324-5600
14.7	700	6.0	High	324-4700
15.3	600	4.6	Moderate	354-6009
15.4	800	6.3	High	324-4800
15.7	700	4.8	Moderate Plus	324-5700
16.3	700	5.0	Moderate	354-7009
16.6	800	5.0	Moderate Plus	324-5800
16.7	800	5.4	Moderate	354-8009

MemoryGel™ Breast Implants Low, Medium and Tall Heights

The following fold-out pages are sizing guides for Two-Stage Breast Reconstruction, when moving from MENTOR CONTOUR® PROFILE™ Tissue Expander to a CPG™ Gel Breast Implant.

Each page reads from left to right:

First locate the size of the CPX™ Expander used during Stage 1.

Now move across to the Stage 2 section. This lists a selection of CPG™ Gel Breast Implants to use as a long-term implant after expansion is complete.

Stage 1 Selection of the CPX™4 and CPX™4 with Suture Tabs Expander*					Stage 2 Selection of the CPG™ Gel Breast Implant to use as a long-term implant after expansion is complete.														
Low Height					CPG™ Implant 311 Cohesive iG™ Low Height, Moderate Projection				CPG™ Implant 312 Cohesive iG™ Low Height, Moderate Plus Projection				CPG™ Implant 313 Cohesive iG™ Low Height, High Projection						
Vol. cc	Width cm	Height cm	Proj. cm	Order Catalog #	Vol. cc	Width cm	Height cm	Proj. cm	Order Catalog #	Vol. cc	Width cm	Height cm	Proj. cm	Order Catalog #	Vol. cc	Width cm	Height cm	Proj. cm	Order Catalog #
250	11.4	8.1	8.1	354-9121	280	11.0	9.8	3.7	334-1100	225	11.0	10.7	4.7	334-1107	240	11.0	9.9	4.9	334-1106
					220	11.5	10.3	3.9	334-1153	225	11.5	10.2	4.9	334-1157	270	11.5	10.3	5.0	334-1156
					270	12.5	11.2	4.2	334-1253	330	12.5	11.1	5.3	334-1257	350	12.5	11.2	5.4	334-1256
350	12.7	9.4	9.4	354-9122	300	13.0	11.6	4.4	334-1303	370	13.0	11.5	5.5	334-1307	395	13.0	11.7	5.6	334-1306
450	14.0	10.2	10.2	354-9123	375	14.0	12.5	4.7	334-1403	465	14.0	12.4	5.9	334-1407	490	14.0	12.6	6.1	334-1406
550	15.0	10.9	10.9	354-9124	460	15.0	13.4	5.1	334-1503	570	15.0	13.3	6.3	334-1507	650	15.0	13.5	6.6	334-1506
650	15.7	11.2	11.2	354-9125	530	15.5	13.8	5.2	334-1553										
					560	16.0	14.2	5.4	334-1603	690	16.0	14.2	6.8	334-1607					
750	16.5	11.9	11.9	354-9126	615	16.5	14.7	5.6	334-1653										

Note how the pages are divided into Low Height, Medium Height and Tall Height.

Stage 1

Selection of the CPX™4 and CPX™4 with Suture Tabs Expander*

LOW HEIGHT				
CPX™4 Expander & CPX™4 w/Suture Tabs Expander Styles 8100, 9100				
Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
250	11.4	8.1	8.1	354-9121
350	12.7	9.4	9.4	354-9122
450	14.0	10.2	10.2	354-9123
550	15.0	10.9	10.9	354-9124
650	15.7	11.2	11.2	354-9125
750	16.5	11.9	11.9	354-9126

Stage 2

Selection of the CPG™ Gel Breast Implant to use as a long-term implant after expansion is complete.

CPG™ Implant 311 Cohesive III™ Low Height, Moderate Projection					CPG™ Implant 312 Cohesive III™ Low Height, Moderate Plus Projection					CPG™ Implant 313 Cohesive III™ Low Height, High Projection				
Vol cc	Width cm	Height cm	Proj cm	Order Catalog #	Vol cc	Width cm	Height cm	Proj cm	Order Catalog #	Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
180	11.0	9.8	3.7	334-1101	225	11.0	19.7	4.7	334-1107	240	11.0	9.9	4.9	334-1106
210	11.5	10.3	3.9	334-1151	225	11.5	10.2	4.9	334-1157	270	11.5	10.3	5.0	334-1156
270	12.5	11.2	4.2	334-1251	330	12.5	11.1	5.3	334-1257	350	12.5	11.2	5.4	334-1256
300	13.0	11.6	4.4	334-1301	370	13.0	11.5	5.5	334-1307	395	13.0	11.7	5.6	334-1306
375	14.0	12.5	4.7	334-1401	465	14.0	12.4	5.9	334-1407	490	14.0	12.6	6.1	334-1406
460	15.0	13.4	5.1	334-1501	570	15.0	13.3	6.3	334-1507	650	15.0	13.5	6.6	334-1506
510	15.5	13.8	5.2	334-1551	690	16.0	14.2	6.8	334-1607					
560	16.0	14.2	5.4	334-1601										
615	16.5	14.7	5.6	334-1651										

MEDIUM HEIGHT				
CPX™4 Expander & CPX™4 w/Suture Tabs Expander Styles 8200, 9200				
Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
275	10.7	9.3	6.2	354-9221
350	11.7	10.0	6.6	354-9222
450	12.7	10.8	7.0	354-9223
550	13.5	11.7	7.4	354-9224
650	14.6	12.6	7.6	354-9225
800	15.6	13.3	8.0	354-9226

CPG™ Implant 311 Cohesive III™ Medium Height, Moderate Projection					CPG™ Implant 312 Cohesive III™ Medium Height, Moderate Plus Projection					CPG™ Implant 313 Cohesive III™ Medium Height, High Projection				
Vol cc	Width cm	Height cm	Proj cm	Order Catalog #	Vol cc	Width cm	Height cm	Proj cm	Order Catalog #	Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
160	10.5	9.4	3.6	334-1051	195	11.0	19.7	4.7	334-1107	210	10.5	9.5	4.8	334-1056
180	11.0	9.8	3.7	334-1101	225	11.5	10.2	4.9	334-1157	240	11.0	9.9	4.9	334-1106
210	11.5	10.3	3.9	334-1151	255	12.5	11.1	5.3	334-1257	270	11.5	10.3	5.0	334-1156
235	12.0	10.7	4.1	334-1201	290	13.0	11.5	5.5	334-1307	310	12.0	10.8	5.2	334-1206
270	12.5	11.2	4.2	334-1251	330	14.0	12.4	5.9	334-1407	350	12.5	11.2	5.4	334-1256
300	13.0	11.6	4.4	334-1301	370	15.0	13.3	6.3	334-1507	395	13.0	11.7	5.6	334-1306
335	13.5	12.0	4.6	334-1351	415	15.0	13.3	6.3	334-1507	440	13.5	12.2	5.8	334-1356
415	14.5	12.9	4.9	334-1451	515	15.0	13.3	6.3	334-1507	545	14.5	13.1	6.3	334-1456
460	15.0	13.4	5.1	334-1501	570	15.0	13.3	6.3	334-1507	606	15.0	13.5	6.6	334-1506
510	15.5	13.8	5.2	334-1551	690	15.0	13.3	6.3	334-1507					
560	16.0	14.3	5.4	334-1601										

Stage 1

Selection of the CPX™4 and CPX™4 with Suture Tabs Expander*

MEDIUM HEIGHT

CPX™4 Expander & CPX™4 w/Suture Tabs Expander Styles 8200, 9200

Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
275	10.7	9.3	6.2	354-9221
35	11.7	10.0	6.6	354-9222
450	12.7	10.8	7.0	354-9223
550	13.5	11.7	7.4	354-9224
650	14.6	12.6	7.6	354-9225
800	15.6	13.3	8.0	354-9226

TALL HEIGHT

CPX™4 Expander & CPX™4 w/Suture Tabs Expander Styles 8300, 9300

Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
250	11.4	8.1	8.1	354-9121
350	12.7	9.4	9.4	354-9122
450	14.0	10.2	10.2	354-9123
550	15.0	10.9	10.9	354-9124
650	15.7	11.2	11.2	354-9125
750	16.5	11.9	11.9	354-9126

Stage 2

Selection of the CPG™ Gel Breast Implant to use as a long-term implant after expansion is complete.

CPG™ Implant 321 Cohesive III™ Medium Height, Moderate Projection

Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
160	10.5	9.9	3.8	354-1058
215	11.0	10.3	3.9	354-1108
245	11.5	10.8	4.0	354-1158
280	12.0	11.3	4.2	354-1208
315	12.5	11.8	4.4	354-1308
355	13.0	12.2	4.6	354-1358
395	13.5	12.7	4.7	354-1351
480	14.5	13.6	5.0	354-1458
530	15.0	14.1	5.2	354-1508
640	16.0	15.0	5.6	354-1608

CPG™ Implant 322 Cohesive III™ Medium Height, Moderate Plus Projection

Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
225	10.5	9.9	4.4	334-1055
255	11.5	10.3	4.7	334-1105
295	11.5	10.8	4.9	334-1155
330	12.0	11.3	5.1	334-1205
375	12.5	11.8	5.3	334-1255
420	13.0	12.2	5.5	334-1305
475	13.5	12.7	5.7	334-1355
585	14.5	13.6	6.1	334-1455
650	15.0	14.1	6.3	334-1505

CPG™ Implant 323 Cohesive III™ Medium Height, High Projection

Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
260	10.5	9.9	5.3	334-1052
300	11.0	10.3	5.6	334-1102
345	11.5	10.8	5.8	334-1152
390	12.0	11.3	6.0	334-1202
440	12.5	11.8	6.2	334-1252
495	13.0	12.2	6.5	334-1302
555	13.5	12.7	6.7	334-1352
585	14.5	13.6	7.1	334-1452

CPG™ Implant 331 Cohesive III™ Tall Height, Moderate Projection

Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
180	11.0	9.8	3.7	334-1101
210	11.5	10.3	3.9	334-1151
270	12.5	11.2	4.2	334-1251
300	13.0	11.6	4.4	334-1301
375	14.0	12.5	4.7	334-1401
460	15.0	13.4	5.1	334-1501
510	15.5	13.8	5.2	334-1551
560	16.0	14.2	5.4	334-1601
615	16.5	14.7	5.6	334-1651

CPG™ Implant 332 Cohesive III™ Tall Height, Moderate Plus Projection

Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
225	11.0	19.7	4.7	334-1107
225	11.5	10.2	4.9	334-1157
330	12.5	11.1	5.3	334-1257
370	13.0	11.5	5.5	334-1307
465	14.0	12.4	5.9	334-1407
570	15.0	13.3	6.3	334-1507
690	16.0	14.2	6.8	334-1607

CPG™ Implant 333 Cohesive III™ Tall Height, High Projection

Vol cc	Width cm	Height cm	Proj cm	Order Catalog #
240	11.0	9.9	4.9	334-1106
270	11.5	10.3	5.0	334-1156
350	12.5	11.2	5.4	334-1256
395	13.0	11.7	5.6	334-1306
490	14.0	12.6	6.1	334-1406
605	15.0	13.5	6.6	334-1506

*PLEASE NOTE: The CONTOUR PROFILE™ Tissue Expander Injection Dome is not MRI compatible. For detailed indications, contraindications, warning and precautions associated with the use of this device, please refer to the Product Insert Data Sheet. Individual implant dimensions may vary slightly in products of this type. Not all units conform exactly to the dimensions noted.

Important safety information.

MENTOR® Breast Implants are indicated for breast augmentation, in women who are at least 18 years old, or for breast reconstruction. Breast implant surgery should not be performed in women with active infection anywhere in their body, with existing cancer or pre-cancer of their breast(s) who have not received adequate treatment for those conditions or who are pregnant or nursing.

There are risks associated with breast implant surgery. Breast implants are not lifetime devices and breast implantation is not necessarily a one-time surgery.

The most common complications with MENTOR® MemoryGel™ Breast Implants include re-operation, implant removal, capsular contracture, asymmetry, and breast pain. A lower risk of complication is implant rupture, which is most often silent.

The health consequences of a ruptured silicone gel-filled breast implant have not been fully established. Screenings such as mammography, MRI, or ultrasound are recommended after initial implant surgery to assist in detecting implant rupture.

The most common complications with MENTOR® Saline-Filled Implants include re-operation, implant removal, capsular contracture, wrinkling, deflation, asymmetry, and breast pain.

MENTOR® CONTOUR PROFILE™ Tissue Expanders and MENTOR® Smooth Round Tissue Expanders with Remote Dome are used for breast reconstruction following mastectomy. These expanders are intended for temporary subcutaneous or submuscular implantation and are not intended for use beyond six months. CONTOUR PROFILE™ Tissue Expanders are devices that contain magnetic injection domes and are NOT MRI compatible. Do not use the CONTOUR PROFILE™ Tissue Expander in patients where an MRI may be needed. DO NOT use the CONTOUR PROFILE™ Tissue Expander in patients that have a previously implanted device that could be affected by a magnetic field.

The device could be moved by the MRI causing pain or displacement, potentially resulting in a revision surgery. The incidence of extrusion of the expander has been shown to increase when the expander has been placed in injured areas: scarred, heavily irradiated or burned tissue, crushed bone areas or where severe surgical reduction of the area has previously been performed.

Your patient needs to be informed and understand the risks and benefits of breast implants, and provided with an opportunity to consult with you prior to deciding on surgery. Your patient needs to be informed and understand the risks and benefits of breast implants, and provided with an opportunity to consult with you prior to deciding on surgery.

For detailed indications, contraindications, warning and precautions associated with the use of all MENTOR® Implantable Devices, please refer to the Product Insert Data Sheet provided with each product, or review the Important Safety Information provided at www.mentorwwllc.eu