

SYMLE MOUSE-

Here is a great solution for those frustrated using computers due to limited arm movement - A Head & face gesture mouse control switch via webcam is a perfect answer. Useful for persons with ALS, Carpal Tunnel, Cerebral Palsy, Multiple Sclerosis, Muscular Dystrophy, Spinal Cord Injury. Don't need the Pesky Wearables, Sip/Puff/Bite Switches, Or Stickers-on-Face!




The picture shows a person lying down and accessing a computer

The Smyle Mouse Advantage is its ease of set up requiring no hardware, accessories to purchase or install, no dots, electrical contacts or headgear on the face, and no bite sticks or sip/puff tubes in the mouth.

It includes a Head mouse, Smile & Dwell Clicker, and Adaptive Switch Emulator allowing use both indoors or outdoors and can be used for work or play, art or entertainment, socializing or education. It is designed for Microsoft Windows 7/8/10

<https://smylemouse.com/>

https://youtu.be/GzFZp_WT4aY



The picture shows a screenshot of a video of a person -the lips and the spectacle frame are showing a green outer frame highlighting of smyle mouse being detected.