

## **Tub Transfer Bench**



The tub transfer bench permits a sliding transfer from a wheelchair into a tub. The bench can perform in this manner because it extends out beyond the edge of the tub. It is also useful for patients who have difficulty in stepping over the side of the tub.

### **Operating Instructions**

- Install the bench by placing the bench in the tub in the appropriate position, placing the faucet end with the 2 outer legs resting on the bathroom floor outside the tub.
- Adjust the outside legs so the seat tilts slightly toward the inside of the tub. This position will allow any water falling on the bench to drain into the tub rather than out on the floor.

### **Precautions**

- The rubber or plastic footings on this bench are slip resistant, not slip proof.
- When sitting or raising your self from the bench, always be careful to only exert force straight down. A force in a horizontal direction may cause the bench to slip.