



KANAD VIDYUT
कणाद विद्युत्



HEMA

Digital Hemoglobinometer



LED Source



NO Filters

"HEMA" DIGITAL HAEMOGLOBINOMETER

- ☞ **For Direct Hb Results.**
- ☞ **No Maintenance**
- ☞ **With AUTOZERO.**
- ☞ **With Hb, CAL and Absorbance modes.**
- ☞ **Numeric Entry of Known Hb for CAL mode.**
- ☞ **Can Use Round or Square Cuvettes.**
- ☞ **Drabkin, SLS, AHD or other non poison methods**



HEMA

KANAD VIDYUT
कणाद विद्युत्

Digital Haemoglobinometer

MAIN FEATURES

- "HEMA" Haemoglobinometer can measure either by Drabkin Cyanmeth method, by SLS method or non poisonous.
- "HEMA" Haemoglobinometer has CALIBRATION mode. Hence quality control can be maintained by the user.
- Time Saving: AUTOZERO Switch saves time in setting Blank.
- "HEMA" Haemoglobinometer has a keyboard for entering Numbers and Modes. Standard Concentration Values can be entered directly for calibrating "HEMA" Haemoglobinometer. Then "HEMA" Haemoglobinometer can give direct reading of hemoglobin of samples.
- "HEMA" Haemoglobinometer gives 3 digit ABS values and Hb values in g/dl, as follows:

Line 1:' ABS=0.364

Line 2:' Hb = 15.0 g/dl '

- "HEMA" Haemoglobinometer are Reliable, Stable and Sturdy in daily operation. No Maintenance: "HEMA" Haemoglobinometer do not have bulbs.
- Light weight Non Rusting Aluminum Body with Epoxy Powder Coating keeps "HEMA" Haemoglobinometer 'Clean and New' even after many years of daily use in Laboratory, Hospital or Camps.
- Very low power consumption: "HEMA" Haemoglobinometer can operate on Torch Cells [Dry cells] or motorcycle batteries, for a long time of many hours of use. Hence "HEMA" Haemoglobinometer can be used outdoors or blood donation camps.

Designed & Manufactured By :: **KANAD VIDYUT कणाद विद्युत्**

Office : 10, Deenanath Mangeshkar Hospital Road, Erandawane Pune-411004

Contact No.: (+91)(India)-9850848722, 020-25462809

Email: kanad.pune@gmail.com

Website:: www.kanadvidyut.in, www.kanadvidyut.com

