


## Anaesthesia Cannula

Disposable ophthalmic cannula used for administering local anaesthesia.



### Product Description

Product Code	Description with Gauge	Needle Length	Qty. per box
MIPL/E6/A4	Retrobulbar Cannula (Atkinson), 19G	38 mm	10 pcs. (Sterile)
MIPL/E6/A5	Retrobulbar Cannula (Atkinson), 21G	38 mm	10 pcs. (Sterile)
MIPL/E6/A1	Retrobulbar Cannula (Atkinson), 23G	38 mm	10 pcs. (Sterile)
MIPL/E6/A2	Retrobulbar Cannula (Atkinson), 25G	38 mm	10 pcs. (Sterile)
MIPL/E6/A3	Retrobulbar Cannula (Atkinson), 27G	38 mm	10 pcs. (Sterile)
MIPL/E6/B4	Peribulbar Cannula (Atkinson), 19G	25 mm	10 pcs. (Sterile)
MIPL/E6/B5	Peribulbar Cannula (Atkinson), 21G	25 mm	10 pcs. (Sterile)
MIPL/E6/B1	Peribulbar Cannula (Atkinson), 23G	25 mm	10 pcs. (Sterile)
MIPL/E6/B2	Peribulbar Cannula (Atkinson), 25G	25 mm	10 pcs. (Sterile)
MIPL/E6/B3	Peribulbar Cannula (Atkinson), 27G	25 mm	10 pcs. (Sterile)
MIPL/E6/C1	Sub Tenon Cannula, 19G	25 mm	10 pcs. (Sterile)



## Features

- Thin wall cannula used for maximum flow
- CE certified

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

## Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)