



Comfort Mobility Inc.


Home Health Care

2707 Temple Drive, Windsor ON N8W 5E5



519-988-1234

Joey Lift

VSL-4000



When it comes to ease of operation, quality engineering and unsurpassed safety -- Bruno's Joey interior scooter and power chair lift sets the standard. Just drive your mobility device onto the platform, press a button and Bruno's Joey lifts your mobility device and tucks it gently into a mini-van or full-sized van. Lifting and stowing your scooter or power chair doesn't get any easier than with a Joey.





Comfort Mobility Inc.
Home Health Care
2707 Temple Drive, Windsor ON N8W 5E5
519-988-1234

Joey Lift VSL-4000



Options

Safety Barrier

Prevents mobility device from entering passenger area during sudden stops



Highlights

- Fully-powered, one-button operation
- Compact design allows for second row seating in most applications
- Self-leveling platform for slightly uneven surfaces
- Easy drive on/off platform with traction coating
- Securement straps and obstruction safety sensing
- Manual backup for peace of mind
- No scooter or power chair modifications
- Optional industry exclusive safety barrier prevents mobility device from entering passenger area during sudden stops
- Potential to reinstall in another applicable vehicle
- Eligible for auto manufacturer potential rebates

All sales, installation and service provided by Comfort Mobility Inc.

Retractable Self-Tensioning Securement Belts

Easy to use, one-hand operation



Back-Off Device

Easily removes back of scooter/powerchair to fit into shorter vehicle openings without having to remove the entire seat

Specifications

- Lift capacity: 350 lbs/159 kg
- Vehicle compatibility: Contact Comfort Mobility Inc.
- Maximum mobility device width: 28 in/71 cm
- Maximum mobility device wheelbase: 40 in/102 cm
- Installed weight: 175 lb/80 kg
- Vertical travel: 27 in/68 cm
- Installation depth: 44 in/112 cm
- 3-year limited warranty

