VANCOUVER OFFICE

5601 E 18th St. #102 Vancouver, WA 98661 PH: (360) 726-6928 FAX: (360) 828-5769



SPOKANE OFFICE

220 W Indiana Ave Spokane, WA 99205 PH: (509) 326-7615 FAX: (509) 464-6266

Western Medical Specialties

Reorder No.	Dimensions in Inches	Pouches per Box/ Boxes per Case	Cost per Box/ Cost per Case
SS03	2¼ x 2¾	200 per box, 20 boxes per case	\$6.50/Box \$130/Case
SS04	2¼ x 4	200 per box, 20 boxes per case	\$6.50/Box \$130/Case
4461	3.5 x 5.25	200 per box	\$5.95/Box
SS08	2¼ x 8	200 per box, 10 boxes per case	\$7.50/Box \$75/Case
SS06	3½ x 5¼	200 per box, 10 boxes per case	\$7.50/Box \$75/Case
SS09	3½ x 9	200 per box, 10 boxes per case	\$7.50/Box \$75/Case
SS10	5 x 10	200 per box, 10 boxes per case	\$17.50/Box \$175/Case
SS14*	7½ x 14	100 per box, 10 boxes per case	\$13.50/Box
SS15*	10 x 15	100 per box, 10 boxes per case	\$16.50/Box
SS19*	12 x 15	100 per box, 10 boxes per case	\$17.50/Box

^{*}Limited Stock on Reorder Numbers SS14, SS15, SS19



ADENNA® SELF-SEALING STERILIZATION POUCHES


ADENNA® Self-Sealing Sterilization Pouches are designed for use in steam (vapor) and ethylene oxide (ETO) sterilizers. They are manufactured in an ISO certified facility under strict quality control. The pouches feature both internal and external indicators to confirm sterility of the contents, meeting the CDC's recommendation.

"The internal indicator is an important feature."

"Pouches are easy to seal and open."







FEATURES AND BENEFITS

- Process Indicators Internal and external indicators that change colors to confirm sterilization is complete.
- 2) Triple Seal, Triple Heat Process The edges of the pouch are sealed by three separately heat-sealed rails, which is stronger than single sealed rail, to prevent separation or tearing by instruments.
- Blue Tinted Transparent Film Easy viewing of instruments and identifying any punctures or tears should they occur.
- Medical Grade Paper Durable medical grade paper, wet strength treated for high performance.
- 5) Wide Self-Seal Adhesive Strip Easily and securely seal pouch.
- 6) Perforated Fold Neatly fold the adhesive strip to form a tight seal.
- Corner Tack Seals Prevent corners of pouch from curling after sterilization.
- 8) Thumb Notch Easy opening of pouch.

