



Clinical Slit Lamps



- 1. Lamp Cap Plug
- 2. LED / Halogen Cover
- 3. Fixation Knob of Lamp Cap
- 4. Aperture Window / Slit Height
- 5. Filter Selection Lever / Display Mark
- 6. Slit Height Control Knob
- 7. 12.5x Eyepieces
- 8. Magnification Select Dial
 - 9. Prism Box
 - 10. Center Knob
 - 11. Inclination Lever
- 12. Slit Width Control Knob
- 13. JOystiCk (Camera Freeze & Capture)
- 14. Joystick
- 15. Slit Width Control Knob
- 16. Base Locking Screw
 - 17. Aperture & Filter Selection Dial
 - 18. Reflective Mirror
- 19. Chinrest
- 20. Fixation Target
- 21. The Mark of Microscope Arm
- 22. The Fixation Knob of Chinrest Paper
- 23. Chinrest Elevation Adjustment Knob
- 24. The Mark of Microscope Arm
- 25. Microscope Arm Locking Knob
- 26. Patient's Handle
- 27. Wheel Shields
- 28. Work Table





Features:

- The ergonomic design of CSL-7U provides a clear and detailed illumination and observations for patients eye.
- The Slit Lamp CSL-7U has a 14mm aperature that can help you observe the whole cornea and iris easily and makes it is convenient to examine the foreign object.
- CSL-7U offers a high eye point and open field eyepieces designed and guarantee a clear image.
- CSL-7U has a powerful LED or halogen bulb for optimum color rendition and illumination.





LED / Halogen Clinal Slit Lamps

CSL-7U

Product Introduction



CSL-7U a classic styling and outstanding converging stereoscope optics, complete with moisture and mildew proof lenses and a three step drum magnification system.



CSL-7U halogen bulb (6V/30W) provides a steady lighting for observation, and optimum color rendition which allows maximum illumination during the observation.



CSL-7U offers a convenient beam splitter, digital camera adapter and CCD adapter. The modular design is an excellent optical system to exam all patients eye.



CSL-7U slit lamps examinations provides a stereoscopic magnified view of the eye structures in details with built in four different filters that can be inserted.



A key feature of the CSL 7U is the maneouvred cross slide joystick offered with a high precision focus control and repeatability measurements for multiple fixation points.



CSL-7U features an operation illumination unit, rotating slit image from 0mm to 14mm, can change the aperature & slight height, filter selection and diffusion lens.



Clinical Slit Lamps

Continuously variable from 0 mm to 14mm (at 14mm, slit becomes a circle)

Continuous from 1 mm to 14mm Ø9mm, Ø8mm, Ø5mm, Ø3mm,

0° to 180° continuously adjustable

Thermal safety, UV, red-free, and blue

Ø2mm, Ø1mm,Ø0.2mm

from vertical to horizontal

Four Steps 5°, 10°, 15°, 20°

6V / 30W halogen bulb

Microscope

Туре	Galilean-Type
Model of Magnifying	Revolving Drum
Eyepiece	12.5x
Total Magnification Ratio	10x 16x 25x (standard configuration)
Field of View	Ø5.5mm Ø8.5mm 13.5mm Ø 22mm Ø34.7mm
Angle Between Eyepieces	13°
Range of P. D. Adjustment	55mm to 78mm
Diopter adjustment	±6D



Input Voltage Input Frequency Input Power Output Voltage Flectrical Safe St 80 mm Green LED

90 mm

100 mm

15mm

30 mm

Power Source

Product Specifications

Illumination

Slit Width

Slit Height

Slit Angle

Slit Inclination

Illumination Bulb

Movement Base

Fine Movement

Vertical Movement

Chin-Rest Parts

Vertical Movement

Fixation Target

Fore and Back Movement

Left and Right Movement

Filter Piece

Aperture Diameter

Input Voltage	100 / 220V ~±10%
Input Frequency	50 / 60 Hz ± 1Hz
Input Power	30 VA
Output Voltage	Illumination bulb 6V Fixation target 3V
Electrical Safe Standard	Conform to Standard IEC601-1, Class I Type B

Dimensions and Weight

Packing Box	740mm x 450mm x 500mm
Total Weight	25 Kg
Net Weight	24 Kg

US. CORPORATE OFFICE 8380 NW 64th ST MIAMI, FL 33169

TEL: 1.305.406.3915 FAX: 1.305.938.5012 CUSTOMERSERVICE@VIEWLIGHTUSA.CON

