



Video Otoscope

Store pictures in NOAH, retrieve results easily for follow-up visits showing the patients their progress

The remarkably clear resolution of the MedRx Video Otoscope is particularly helpful when comparing a clean ear canal with normal eardrum to a canal that has excessive cerumen, foreign bodies, external otitis, PE Tube status, cholesteatoma, TM perforations or other abnormalities. The Complete Package: Otoscope Probe, High Resolution Color Camera, Light Source and Software

Features

- Capture, View and Store Digital Images
- Perform an Otoscopic Examination and Place the Probe Microphone Tube in the Ear Canal
- Conveniently Shows a List of Captured Pictures, Comments and Previews of the Selected Images
- View a Picture then Zoom In and Enhance it
- NOAH™, TIMS®, Blueprint OMS and Sycle Compatible
- An Integrated Video Module is Standard in all MedRx Software

Analog Camera with Probe

Plug the Video Otoscope straight into a monitor or connect to your computer via a USB Capture Dongle

The LED Light Handle

The 360° rotatable Video Otoscope Handle produces virtually no heat and comes with two sets of rechargeable batteries



USB Capture Dongle

MedRx®
www.medrx-usa.com

Call Today! 888-392-1234