

Timpac Freezpen Made in Belgium 



cryosurgery
made simple



Timpac Freezpen is a powerful, precise cryosurgery instrument for the removal of benign skin lesions.

- **Simply powerful**
High power (55-bar) microjet achieves in-tissue temperature of -89°C, for complete tissue necrosis in seconds.
- **Simply precise**
Ultrafine sharp jet for total precision and no damage to healthy tissue. Range of applicators.
- **Simple procedure**
Assess the lesion, choose the right applicator, point, freeze, wipe, freeze again. Done!
- **Simply painless**
No pain, no cuts or bleeding, no anesthesia - even for children. Slight sting signals when treatment is complete.
- **Simply perfect**
Perfect results with no scarring. Healthy tissue is unaffected; frozen tissue disappears in a few weeks.
- **Simply practical**
Easy handling and transport. Handy disposable cartridges; store at room temperature.
- **Its simply good business!**
Each treatment takes just seconds. Return on investment in 10 weeks at just six treatments/week.

What can be treated with Timpac Freezpen?

Acne (scars)	Lichen planus
Acrochordon	Lichen sclerosus
Actinic cheilitis	Lichen simplex
Actinic keratosis	Melanotic macules of the lips
Angiocheratoma of Mibelli	Milii
Angiocheratoma of the scrotum and vulva	Milium colloid
Angiofibromas	Mucoid cyst
Angioma	Nevoid hyperkeratosis of the nipple and the mammary areola
Annular granuloma	Nevus sebaceous
Ano-genital warts	Nodular chondrodermatitis of the helix
Callosity	Nodular prurigo
Capillary hemangioma	Papillary hyperplasia of the palate
Cutaneous leishmaniasis	Perforating dermatoses
Cutaneous sarcoidosis	Piercing elastosis
Dermatofibroma	Plantar warts
Dermatoses due to mechanical traumas	Porokeratosis linear
Dermatosis papulosa nigra	Pruritus anirhinophyma
Facial eosinophilic granuloma	Pyogenic granuloma
Fibrous papule of the face	Rosacea
Filiform and typed warts	Seborrheic keratosis
Fixed discoid erythematosus lupus	Stuccocheratosi
Flat warts	Syringoma
Granuloma fissuratum	Tilomi
Herpesvirus infections	Viral lesions in HIV positive and negative patients
Horns	Vulgar warts
Idiopathic guttata hypomelanosis	Warty carcinoma
Infectious mollusk	Warty epidermal nevus
Keloid and hypertrophic scars	Xantomi Kaposi
Keratoacanthoma	
Lentigo solar	
Leukoplakia	

This list was established based on the results of an extensive study performed by Dr. Renata STRUMIA surgeon, specialist in dermatology and venereology as described in her book La Crioterapia in Dermatologia (p. 261-262).


Freezpen can remove most common skin pathologies in just a few seconds.



Pigmented Spot



Hemangioma



Warts



Molluscum Contagiosum



Acrochordon



Skin Tag



Pigmented Spot



Seborrheic Keratosis



Wart



Skin Tag



Wart




Solar lentigo



Before treatment



After treatment





8g Freezen set:

Freezen 8g instrument set - 2 applicators
Blue and white code applicators
User manual and tweezers for filter placement.

16g Freezen set:

Freezen 16g instrument set - 3 applicators
Blue, white and green code applicators
User manual and tweezers for filter placement.



Box with 15 N₂O
cartridges of 8g each
carton with 12 boxes.



Box with 6 N₂O
cartridges of 16g each
carton with 12 boxes.



- 1-3mm applicator ----- |
- 2-6mm applicator ----- ||
- 4-8mm applicator ----- |||



Timpac Healthcare Private Limited

An ISO 9001:2008 Company

C-168, 1st Floor, Naraina Industrial
Area Phase-I, New Delhi-110028 India

Tel: +91-11-45673300, 25794065

E-mail: wecare@timpac.in

Timpac Freezen Made in Belgium