


## Irrigating & Aspirating Cannula

Disposable ophthalmic cannula used for maintaining or forming anterior chamber and irrigating cornea during surgery. It may be used for injection or removal of air, fluids & viscoelastic substances during an ophthalmic surgery.



### Product Description

Product Code	Description with Gauge	CE Certified	Qty. per box
MIPL/E11/A6	I & A Cannula (Simcoe Regular) - 21G	No	10 Pcs. (Sterile)
MIPL/E11/A1	I & A Cannula (Simcoe Regular) - 22G	Yes	10 Pcs. (Sterile)
MIPL/E11/A2	I & A Cannula (Simcoe Regular) - 23G	Yes	10 Pcs. (Sterile)



## Features

- Thin wall cannula used for maximum flow
- CE certified

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

## Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)