





JHHS Supply Chain Coronavirus Supply Request

March 21, 2020

Testing Materials

MPN	Manuf.	Description	Specifications	Image
220531	BD	Universal Viral Transport System	<ul style="list-style-type: none"> Designed to transport viruses, chlamydiae, mycoplasmas and ureaplasmas at room temperature May be stored and transported at 2°C to 25°C 	
220246	BD	BD Eswab Flocked Swab	<ul style="list-style-type: none"> Collect clinical specimens containing aerobic, anaerobic and fastidious bacteria from the collection site, and transport them to the testing laboratory Flocked swab ensures optimal elution of the specimen into the transport medium 	
R723480	Remel	Eswab (White Cap-Regular)	<ul style="list-style-type: none"> Uses: Wounds 	
R723481	Remel	Eswab (Green Cap – Mini Tip)	<ul style="list-style-type: none"> Uses: Pediatric Samples, Eyes, Ears 	





Additional Swabs (cotton swabs do not work):

- E-Swab by Copan (Catalogue # 481C and 482C) with regular or flex minitip applicator
- Opti-Swab by Puritan (Catalogue # LA-117), swab included in kit (Catalog#3317-H).
- Puritan: 25-3317-H, 25-3316-U, 25-3316-H, 25-3317-U, 25-3318-U, 25-3318-H, 25-3319-H, 25-3320-U, 25-3320-U EMB 100MM, 25-3320-U EMB 80MM, 25-3320-H and 25-3320-H EMB 80MM
- Copan: 501CS01, 503CS01, 516CS01, 518CS01 and 534CS01

Transport Medium:

- Vial Transport from Puritan (Catalogue UT-100 and UT-300) for 1ML and 3ML sizes
- Hanks Balanced Salt Solution (HBSS), with calcium and magnesium, no phenol red, 500 ml or 1 liter bottles



N95 Masks

MPN	Manuf.	Description	Specifications	Fit Tested?	Image
1860	3M Health	MASK PARTICULATE RESP N95 REG	<ul style="list-style-type: none"> • BFE: 99% • NIOSH approved N95 • FDA cleared for use as a surgical mask • Fluid resistant according to ASTM F1862 	Yes	
46767	Halyard Health	MASK RESP N95 REG W/SAFETY	<ul style="list-style-type: none"> • BFE: $\geq 99\%$ • PFE: $\geq 99\%$ • Splash Resistance: 160 • ASTM Level of Protection: 3 	TBD	
1870+	3M Health	MASK RESPIRATOR N95 1870+	<ul style="list-style-type: none"> • BFE: $\geq 99\%$ • FDA cleared for use as a surgical mask • Fluid Resistance: 160 mmHg 	Yes	
8210	3M Health	MASK PARTICULATE RESPIRATOR	<ul style="list-style-type: none"> • BFE: $\geq 95\%$ • FDA cleared for use as a filtration mask • Fluid Resistance: None 	NO	



Surgical Masks

MPN	Manuf.	Description	Specifications	Image
AT7509	Cardinal Health	MASK SURGICAL FILTER DISP	<ul style="list-style-type: none"> • BFE: ≥95% • PFE: ≥95% • Delta P: <4.0 • Fluid Resistance: 80 mm Hg • Flame spread: Class 1 	
AT73035	Cardinal Health	MASK SENSITIVE SKIN WH LF		
AT71235	Cardinal Health	MASK FOG FREE SURG TIE BLUE LF		
1001124	Halyard Health	MASK SURGICAL W/TIES SUBMICRON	<ul style="list-style-type: none"> • BFE: ≥99% • PFE: ≥99% • Fluid Resistance: None 	





Procedure Masks

MPN	Manuf.	Description	Specifications	Image
AT70021	Cardinal Health	<u>MASK PROCEDURE YELLOW W/EAR LOOP</u>	<ul style="list-style-type: none"> • BFE: $\geq 95\%$ • PFE: $\geq 95\%$ • Delta P: < 4.0 • Fluid Resistance: 80 mm Hg • Flame spread: Class 1 	
AT74531	Cardinal Health	<u>MASK FOG FREE IMPERV EAR LP LE</u>	<ul style="list-style-type: none"> • BFE: ≥ 98 • PFE: $\geq 98\%$ • Delta P: < 5.0 • Fluid Resistance: 160 mm Hg • Flame spread: Class 1 	





Surgical Gowns

MPN	Manuf.	Description	Specifications	Image
9041	CARDINAL HEALTH	Gown Surg AAMI 4	<ul style="list-style-type: none"> SmartSleeve™ Poly-Reinforced Breathable Sleeve 	
9571	CARDINAL HEALTH	GOWN REINFORCED - AAMI 3	<ul style="list-style-type: none"> Fabric-Reinforced AAMI Level 3 Surgical Gowns 	

Isolation Gowns

MPN	Manuf.	Description	Specifications	Image
NONTH150	MEDLINE INDUSTRIES	GOWN ISO POLY THUMBS UP REG BLUE	<ul style="list-style-type: none"> • Polyethylene Thumb Loop Style • Protection Level: ASTM F1670 	
CR15010	MEDLINE INDUSTRIES	GOWN ISO THMBLP POLYETHYLENE BLUE	<ul style="list-style-type: none"> • Polyethylene thumb loop • Over-the-head style neck with dual side ties 	
47003	HALYARD HEALTH	GOWN ISO MEDWT COATED OVERHEAD REG YEL	<ul style="list-style-type: none"> • Polyethylene coated spunbond fabric, Tie waist, Thumb hooks 	
NON26278	MEDLINE INDUSTRIES	GOWN ISO POLYCTD OVR-HD OPN BK KNT CF LG	<ul style="list-style-type: none"> • Overhead gown with thumb loops • Material passes the ASTM 1670, 1671 test 	

PAPR Hoods

MPN	Manuf.	Description	Specifications	Image
BE1020	3M	PAPR HOOD ENHANCED PRECAUTIONS W/DRAPE R	<ul style="list-style-type: none"> • APF: 1000 • Double bib hood 	
BE10L3	3M	PAPR HOOD ENHANCED PRECAUTIONS W/DRAPE R		
S-433L-5	3M	HOOD VERSAFLO M/L WHITE	<ul style="list-style-type: none"> • APF: 1000 • Material: Polypropylene coated non-woven • Operating temp range: -10°C to 55°C 	
S-103L-20	3M	HEADCOVER VERSAFLO M/L	<ul style="list-style-type: none"> • APF: 25 • CFM Capacity per Minute: 100 CFM, 50 CFM • Material: Polypropylene coated non-woven • Operating temp range: -10°C to 55°C 	
S-103S-20	3M	HEADCOVER VERSAFLO STANDARD S/M		