


INVACARE®
TDX® SP
POWER WHEELCHAIR



The **Invacare® TDX® SP Power Wheelchair** is the first model in the next generation of TDX (Total Driving eXperience) power wheelchairs. The TDX SP wheelchair takes all the best features of the TDX series and improves upon them with enhanced performance, superior ride quality, quieter chassis and an elegantly simple design.

Features

The available Invacare® Formula™ CG Powered Seating System which features up to 7" of weight shift tilt Rehab (ASBA) seat that is growable in width and depth Semi-recline captain's seat is available in 16", 18", 20" and 22" widths

Optional heavy-duty package

Transport Ready Option (TRRO) available*

Invacare Core Technologies


Invacare® SureStep® Suspension with Stability Lock provides a smooth and stable transition

Invacare® G-Trac™ Technology provides precise tracking for imprecise environments

Invacare® MK6i™ Electronics are the benchmark in innovative technology

Motion Concepts power positioning systems can be interfaced with the Invacare TDX SP power base. Contact Motion Concepts for more information at 888-433-6818

HCPCS Code: K0848, K0849, K0850, K0851, K0856, K0858, K0861, K0862



1. Enhanced Invacare® SureStep® Suspension with Stability Lock
2. Invacare® G-Trac™ Technology
3. Invacare® MK6i™ Electronics
4. Motion Concepts Powered Positioning System Compatibility

*TRRO includes four factory-installed brackets and wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form No. 11-099
Rev. 09/12
120335



Model #	Overall Dimensions**	Seat Dimensions	Seat-to-Floor Height	Drive Wheels	Ground Clearance	Turning Radius	Incline
TDXSP, TDXSP-CG, TDXSP-MCG, TDXSP-CG-GT	41.5" L 27.6" W 35.5" H	12"-24" x 12"-22"	16.5"-21"	14"	3"	20.5"	9°
	Total Weight	Weight Capacity	Motor	Max Speed	Battery	Battery Range	Limited Warranty
	Up to 399 lb.	300 lb. 400 lb. HD	4-Pole	5.8 mph 5 mph HD	Group 24, 22NF	Up To 20 Miles	5 Years Base Frame 3 Years Seat Frame 1 Year Electronics/ Motor/Gearbox

**Rehab seat and 90° swingaway footrest
Do not operate on roads, streets or highways.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.