

Clinical Slit Lamp



C64

0



U



Features:

- The ergonomic design of CSL-11U provides a clear and detailed illumination and observations for patients eye.
- The Slit Lamp CSL-11U is a basic model offering specialized optics and increased magnification levels.
- This basic design allows the operator to easily maneuver between using a two-position of magnification. (10x & 16x)
 CSL-11U slit lamps utilized a Galilean Optical system.
- CSL-11U examinations provides a stereoscopic magnified view of the eye structures in details.



Magnification Ratio 10x 16x



Maneouvred Cross Slide Joystick



Halogen Clinal Slit Lamps





CSL-11U A classic style and outstanding Galilean optics, complete with moisture and mildew proof lenses.



CSL-11U slit lamps examinations provides a stereoscopic magnified view of the eye structures in details with built in filters.



CSL-11U halogen bulb (12V 30W) provides a steady lighting for observation, optimum color rendition which allows maximum illumination during the observation.




CSL-11U features a illumination tower allows examination in horizontal, vertical, and oblique section including option to offset the slit for corneal trans-illumination.



A key feature of the CSL-11U is the short reflex mirror which is recommended for posterior segment examinations. The illumination beam must be in a 3-5 degree angle.



Viewlight's CSL-11U Tonometer can be mounted and installed as an accessory used for measuring ocular pressure and applanation of the corneal surface.



- 12. Magnification Changer Lever
- 13. Center Knob
- 14. Inclination Lever
- 15. Slit Width Control Knob
- 16. Patient's Handle
- 17. Angle Mark Ring
- 18. Joystick (Capture
- 19. Base Locking Screw
- 20. Wheel Shields
- 21. Single Device Table
- 22. Actuator
- 23. Up and Down Electrical Swithces
- 24. External / Power Switch
- 25. Wheel Locks



21



Clinical Slit Lamp

Product Specifications

		No. of the second se
Microscope	Microscope Type	Converging Seterocope Microscope
	Magnification Selection	Two Steps Magnification
	Eyepiece	10x
	Magnification Ration (Field of View)	10x (18mm), 16x (14.5mm)
	Pupillary Adjustment	50mm - 82mm
	Diopter Adjustment	± 5D
Base	Side Shift	100mm
	Depth Shift	90mm
	Height Shift	30mm
Chinrest	Height Shift	80mm
	Fixation Target	Red LED
Illumination System	Slith Width	0-9mm Continuous (at 9mm, slit becomes a circle)
	Slit Height	1-8mm Continuous
	Slit Angle	0 - 180 adjustable
	Slit Inclination	0 -23
	Diameter of Light Spot	9mm, 8mm, 5mm, 3mm,2mm, 1mm, 0,2mm
	Filter	Heat Absorption, Grey, Red-Free, Cobalt Blue
Power	Illumination Bulb	12v 50w halogen lamp
	Input Voltage	110V / 220V AC
	Frequency	50/60Hz
	Power Consumption	68VA
Weight & Size	Gross Weight	18kg
	Dimension	690mm (L) x 440mm (W) x 420mm (H)

US. CORPORATE OFFICE 8380 NW 64th ST MIAMI, FL 33169 TEL: 1.305.406.3915 FAX: 1.305.938.5012 CUSTOMERSERVICE@VIEWLIGHTUSA.CO

