

Hygienists, meet your new assistant...

The Mirror Suction

Imagine having your mouth mirror, aspiration, and retraction all in one instrument.

Water flow can be better controlled during ultra-sonic scaling. Sealants can be applied without assistance.

That's why the Mirror Suction provides strong suction power by connecting to your HVE valve.

- **Suit Your Preference** with double and single-sided mirrors
- **Easy-to-Clean** because all components are autoclavable
- **Reduces Fogging** through airflow in the mirror head
- **For best results** use in tandem with the saliva ejector

Order from your dealer using the below part numbers.



Double-Sided (#5 mirror head)



Single-Sided (#4 mirror head)

Mirror Suction - Double-Sided

| | |
|--------|---|
| 254069 | Mirror Suction Double-Sided "Introductory" Kit: (5) Handles with #5 double-sided mirrors, (2) Hoses, (2) 11 mm adapters and (1) Cleaning brush |
| 254070 | Mirror Suction Double-Sided - "Mini" Kit - (2) Handles with (2) #5 double-sided mirrors, (1) Hose, (1) 11 mm adapter & (1) Cleaning brush all packaged in a sterilizable cassette |
| 254077 | Mirror Suction Double-Sided #5 Mirror Head - Pack of 6 heads |

Mirror Suction - Single-Sided

| | |
|--------|--|
| 254050 | Mirror Suction "Introductory" Kit: (5) Handles with #4 mirrors, (2) Hoses, (2) 11 mm adapters and (1) Cleaning brush |
| 254068 | Mirror Suction "Mini" Kit: (2) Handles with #4 mirrors, (1) Hose, (1) 11 mm adapter and (1) Cleaning brush |
| 254058 | Mirror Suction "Single-Sided" #4 Mirror Head - Pack of 6 heads |



Mirror Suction Introductory Kit



Mirror Suction Mini Kit



(800) 328-2335
www.hagerworldwide.com



BioGuard[®]

by **BULLFROG**

Biological Vacuum Cleaning System

- » Only Needs to be Used Twice a Week
- » Amalgam Separator Friendly - pH 7.9
- » Provides Continuous Residual Action for the Entire Evacuation Line System
- » Non-Foaming Formulation
- » Eliminates Offensive Odors - Leaves a Fresh Scent
- » Environmentally Friendly - Biodegradable

Maintain Maximum Evacuation System Suction



Before
BioGuard



After
BioGuard

Save Time
Only needs to
be used
twice a week!

*"One operatory had 2 HVE lines with one barely functioning.
It now has great suction."*

Deborah
RDH, Florida

Request Your **FREE BioGuard Evaluation Kit
by calling (800) 328-2335**



(800) 328-2335 • www.hagerworldwide.com