



5614 N. Union Blvd.  
Colorado Springs, CO 80918  
719.471.3200  
[www.lowvisionsoco.com](http://www.lowvisionsoco.com)

## Using your LED Stand Magnifier

- Hold the base of the stand magnifier directly on the page and keep it there
- Adjust your distance from the magnifier until you see it as clearly as possible. The closer your eye is to the magnifier, the more area you will see on the page.
- To read, slide the magnifier slowly across the page without picking it up
- Since your magnifier has its own illumination, avoid overhead lights and windows to decrease glare and reflections
- Reading may be tiring at first. Start with large print and read only a few minutes at a time. Gradually decrease print size and increase reading time.
- Reading stands may allow improved postural comfort while reading
- Clean your magnifier with a mild soap (non-alcohol) and water and dry with a soft clean cloth (no Kleenex).
- Your LED bulbs are "lifetime" bulbs and do not need to be replaced. The LED handle uses 3 AA batteries, which will last about 30-40 hours.
- If you have any questions, call our office at **719-471-3200**

**Remember: Patience & Practice**