



"THIS CHAIR HAS MET AND EXCEEDED MY EXPECTATIONS. I'VE GONE MANY MILES AND HOURS ... IN THE ALMOST FOUR YEARS I'VE HAD IT. IT HAS NEVER BROKEN DOWN ON ME OR LET ME DOWN. IT'S ALWAYS GOT ME HOME, THAT'S IMPORTANT. I HAVE CONFIDENCE WHEN I LEAVE ON MY DAILY EXCURSIONS."

Gary Hays - Texas, USA

### FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.



### REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be custom designed to a fixed position. It's whatever you need.

### UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



### TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.

# MAGIC MOBILITY FRONTIER V6 ALL TERRAIN

AUSTRALIAN MADE  
-  
BUILT TO LAST  
-  
CUSTOMIZED FOR YOUR LIFESTYLE



201 Growth Parkway  
Angola, IN 46703  
United States

Tel 1-260-665-2769

Email [iim@mobility-usa.com](mailto:iim@mobility-usa.com)

Web [mobility-usa.com](http://mobility-usa.com)



# MAGIC MOBILITY FRONTIER V6 ALL TERRAIN



For ultimate indoor-outdoor freedom. This is it, our most popular wheelchair worldwide due to its versatility. An engineering beauty that adapts to any environment. The Frontier V6 All Terrain is the ideal outdoor wheelchair, yet its size makes it a capable indoor performer. Its suspension and outdoor performance lets you tackle the bumpy terrain life throws at you.

#### Basic specifications:

- Length: 41½"
- Width: 28"
- Seat to floor height: 17"
- Turning radius 23"
- Weight limit: 400lbs  
340lbs with lift and tilt
- Top speed: 6.2mph
- Motors: 800W Gear in Line
- Wheels: Off-road 14"  
(145/70-6) black
- Castors: 2.50-4 black
- Batteries: Group 24, 75 Ah  
LA/AGM and Gel (Group 34,  
60 Ah Gel option)
- Penny and Giles controller  
with LCD display (R-Net)
- Centre post legrest with  
one piece footplate
- Sling backrest and rehab  
seat pan (from 12"-23" wide)\*
- WC19 crash test compliant  
tie down points.

#### Rehab positioning:

- Power options: 12" seat lift,  
tilt (45 degrees), power recline,  
lie down, power elevating  
legrests, transfer tilt (8 degrees  
and up to 2½" floor to seat  
height reduction)
- Seating and Backrests: MPS,  
rehab backrest and adjustable  
seat pans and custom seating
- Headrests and mounts: MPS,  
rehab
- Legrests: Centre post, side mount,  
swingaway, power elevating
- Footplates: 1 or 2 piece fixed/  
flip up/angle adjustable, any  
custom requirement
- A full range of armpads,  
straps, pads, harness, restraints,  
supports are all available from  
our extensive catalogue.

#### Options and accessories:

- 24V LED Lighting kit
- USB charger
- Jack
- Luggage rack, bag hooks, gun  
mounts, bags, pockets, keep safe  
secret boxes, tables, drinks  
holder, oxy and vent carriers
- Special controls: Attendant,  
head, chin, environmental  
controls all available
- Hybrid indoor wheel kit: narrower,  
grey and black drive wheels and  
castors, when off-road capabilities  
are not needed.  
This also reduces the size of  
the chair further.

#### The right wheelchair for:

- Adventurous off-road recreation
- Everyday off-road, country living  
and working
- Urban outdoor living and working
- Roomy indoor spaces.

#### Base options:



V6 HYBRID - BLACK WHEELS

"WE'VE BEEN  
CREATING CUSTOM  
MADE VEHICLES  
FOR 22 YEARS,  
SO THAT YOU CAN  
LIVE THE LIFE  
YOU CHOOSE"

Ashley Daff  
Director




#### ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.

#### BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



#### ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

#### PROUDLY AUSTRALIAN OWNED AND DESIGNED. ASSEMBLED IN THE USA

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

\*Specifications subject to change.  
Check online for current version  
\*Various model options excluded

**Magic**  
MOBILITY  
BECAUSE THE WORLD'S NOT FLAT

Visit us online to see how we change lives every day:  
[magicmobility.com.au](http://magicmobility.com.au)

**Postal Address**  
Magic Mobility Pty Ltd  
2/112 Browns Road  
Noble Park VIC 3174 Australia

**Tel** +61 3 8791 5600  
**Fax** +61 3 9701 0144  
**Email** [sales@magicmobility.com.au](mailto:sales@magicmobility.com.au)



Join the Magic Mobility social conversations!

