

Contact Lens Wearing Schedule

Manufacturer	Lens	Wearing Schedule	Should I Sleep in Them?	Replace Every	Dk Value
Alcon	Air Optix w/Hydraglyde	Daily wear; remove nightly	No	30 Days	110
Alcon	Air Optix Focus Dailies	One day – throw away	No	Day	26
Alcon	Aqua Comfort Dailies	One day – throw away	No	Day	26
Alcon	Dailies Total 1	One day – throw away	No	Day	140
Alcon	Dailies Total 1 for Astigmatism	One day – throw away	No	Day	140
Alcon	Air Optix Astigmatism	Daily wear or extended wear up to 6 nights	Yes – if wearing as extended	30 Days	110
Alcon	Air Optix Night & Day	30-day continuous wear	Yes	30 Days	140
Alcon	Air Optix w/HG MFs	Daily wear or extended wear up to 6 nights	Yes – if wearing as extended	30 Days	110
Alcon	Precision Dailies	One day – throw away	No	Day	90
Alcon	Air Optix Colors	Daily wear; remove nightly	No	30 Days	110
Alcon	Freshlook Dimensions	Daily wear; remove nightly	No	2 Weeks	16
Alcon	Freshlook Color Blend	Daily wear; remove nightly	No	2 Weeks	16
Alcon	Focus Dailies	One day – throw away	No	Day	26
Alcon	Precision 1/Astig.	One day – throw away	No	Day	90
Alcon	Total 30/Astig.	Daily Wear	No	30 Days	123
Bausch and Lomb	Bio-True 1-Day Astig/Presbyopia	One day – throw away	No	Day	42
Bausch and Lomb	Eureka	One day – throw away	No	Day	42
Bausch and Lomb	Infuse/Infuse Asti.	One day – throw away	No	Day	107

Bausch and Lomb	Optima 38	Daily wear	No	30 Days	8.5
Bausch and Lomb	PureVision 2 (All)	Extended wear	Yes	30 days	91
Bausch and Lomb	Ultra Monthly	Daily wear or extended wear up to 6 nights	Yes - if wearing as extended	30 Days	114
Clerio	H2OExtreme Daily	Daily wear; remove nightly	No	30 Days	21
Clerio	H2O Extreme Monthly	Daily wear; or extended wear up to 6 nights	Yes – if wearing as extended	30 Days	21
Clerio	H2O Extreme 54%	Daily wear; or extended wear up to 6 nights	Yes – if wearing as extended	30 Days	21
Clerio	H2O Extreme 54% Toric	One day – throw away	No	Day	21
Clerio	H2O 59% XTRA	One day – throw away	No	Day	21
CooperVision	MyDay 1-Day	One day – throw away	No	Day	80
CooperVision	MyDay Energys Multifocal/Toric	One day – throw away	Np	Day	80
CooperVision	Avaira Vitality/Toric	Daily wear	No	2 Weeks	90
CooperVision	Biofinity (All)	Extended wear up to 6 nights	Yes	30 Days	128
CooperVision	Clariti (All)	One day – throw away	No	Day	60
CooperVision	MiSight (*)	Myopia management – 1 Day-throw away	No	Day	27
CooperVision	ProClear (All)	Daily wear	No	30 days	21
Johnson & Johnson	Acuvue 1-Day Define	One day – throw away	No	Day	21
Johnson & Johnson	Acuvue 1-Day Moist (Also Astigmatism)	One day – throw away	No	Day	21
Johnson & Johnson	Acuvue 2	1 week extended wear; 2 weeks daily wear	Yes – if wearing as extended	2 Weeks if DW 1 Week if EW	28
Johnson & Johnson	Acuvue Oasys (All)	Daily wear; or extended wear up to 6 nights	Yes – if wearing as extended	2 Weeks	103
Johnson & Johnson	Acuvue Oasys for Astigmatism	Daily wear; or extended wear up to 6 nights	Yes – if wearing as extended	2 Weeks	103
Johnson & Johnson	Acuvue Oasys Multifocals	Daily wear; or extended wear up to 6 nights	Yes – if wearing as extended	2 Weeks	103

Johnson & Johnson	Acuvue Oasys Transitions	Daily wear; up to two weeks	No	2 weeks	103
Johnson & Johnson	Acuvue Oasys 1-Day	One day – throw away	No	Day	103
Johnson & Johnson	Acuvue Vita	Daily wear; remove nightly; replace monthly	No	30 Days	103
Johnson & Johnson	Acuvue Vita for Astigmatism	Daily wear; remove nightly; replace monthly	No	30 Days	103

(*) Recommended only for children 8-12 years old for myopia control.

(Compiled August 2023)