

SOLARIS 3-axis high load, low profile, high precision motorized positioner

Eric Acome, Eric Van Every, Alex Deyhim, and Marcin Zajac

Citation: [AIP Conference Proceedings](#) **1741**, 030021 (2016); doi: 10.1063/1.4952844

View online: <http://dx.doi.org/10.1063/1.4952844>

View Table of Contents: <http://scitation.aip.org/content/aip/proceeding/aipcp/1741?ver=pdfcov>

Published by the [AIP Publishing](#)

Articles you may be interested in

[Development of compact high precision two degree of freedom X Y piezoelectric stepping positioner](#)

Rev. Sci. Instrum. **79**, 026110 (2008); 10.1063/1.2841807

[Development of compact high precision linear piezoelectric stepping positioner with nanometer accuracy and large travel range](#)

Rev. Sci. Instrum. **78**, 075112 (2007); 10.1063/1.2756627

[Design method of three-dimensional flux magnetic leakage profiles for high precise motor drive](#)

J. Appl. Phys. **91**, 6970 (2002); 10.1063/1.1452668

[Simple, high precision linear-motor-driven XYΘ positioner \(Walker\)](#)

Rev. Sci. Instrum. **62**, 2492 (1991); 10.1063/1.1142220

[Precision low cost V block for heavy loads](#)

Rev. Sci. Instrum. **51**, 393 (1980); 10.1063/1.1136193

SOLARIS 3-Axis High Load, Low Profile, High Precision Motorized Positioner

Eric Acome¹, Eric Van Every¹, Alex Deyhim^{1, a)} and Marcin Zajac²

¹ADC USA Inc. 126 Ridge Road Lansing NY, 14882 USA, ²National Synchrotron Radiation Centre Solaris Jagiellonian University ul. Czerwone Maki 98/p.3.03 0-392 Krakow, Poland ^{a)}adc@adc9001.com

Abstract. A 3-axis optical table, shown in Figure 1, was designed, fabricated, and assembled for the SOLARIS synchrotron facility at the Jagiellonian University in Krakow, Poland. To accommodate the facility, the table was designed to be very low profile, as seen in Figure 2, and bear a high load. The platform has degrees of freedom in the vertical (Z) direction as well as horizontal transversal (X and Y) directions. The table is intended to sustain loads as large as 1500 kg which will be sufficient to support a variety of equipment to measure and facilitate synchrotron radiation. After assembly, the table was tested and calibrated to find its position error in the vertical direction. ADC has extensive experience designing and building custom complex high precision motion systems [1,2].




FIGURE 1. (a) Computer assisted design model of the 3-axis platform (b) actual 3-axis platform



FIGURE 2. Overall height of motorized positioner (low profile)

Specifications

Each of the motorized stages has a normally closed limit switch. For vertical and horizontal transversal directions of motion, a linear encoder and zero position markers must be provided. The three axes of motion are supported on THK rails. A NEMA 23 stepper motor in combination with a planetary gearbox serves to drive the table in each direction of motion. All three motors drive a ball screw and are equipped with adjustable limit switches that can change the motion within the maximum range of travel. The specifications of the system including range of motion, resolution, weight, and load data are presented in Table 1.

TABLE 1. Key specifications of the 3-axis platform

Description	X-Axis (horizontal)	Y-Axis (horizontal)	Z-Axis (vertical)
Range of Motion	+387mm/-387mm	+20mm/-20mm	+12.5mm/-12.5mm
Resolution (unit/step)	5 micron/full step	1.1 micron/full step	0.1 micron/full step
Minimum dynamic load capacity	1500 kg		
Estimated System Mass	993 kg		
Encoder Manuf.	Renishaw		

Following design, fabrication, and assembly of the system, a test was performed to assess the accuracy of the platform's vertical motion control. The system was loaded appropriately, and the Keyence high accuracy optical micrometer was controlled with a stepper motor controller and used to take the measurements associated with the system's mechanical repeatability. This data is tabulated and presented graphically in Table 2.

TABLE 2. Mechanical repeatability data for the 3-axis platform

Mechanical Repeatability Test Record									
FORM-00009 Rev1.6, 2015-01-29									
Measurments taken using the Keyence high accuracy optical micrometer, controlled using a stepper motor controller.									
					Device Name			3 Axis Platform	
Microstepping	64	ustep/step			System Serial Number:			150624-1	
Run Velocity	1000	(steps/sec)			Job number:			14-044	
Accel. & Decl.	1000	(steps/sec*sec)			Test Date			6/24/2015	
Steps Made per Move	30000	steps/position			Open or Closed Loop?			Open	
Motor Current Rating	4.2	Amps							
Run Current	1.5	Amps			Description Text in Yellow				
Controller % Hold Current	24	%			Entered data in Green				
Hold Current	0.36	Amps			Calculated values in Blue				
					Grey is unused cells				
Z Axis -									
Position	Actual measured move (mm)						Uni-Directional error	Bi-Directional error	
	+ direction 1'st pass	- direction 2'nd pass	+ direction 3'rd pass	- direction 4'th pass	+ direction 5'th pass	- direction 6'th pass			
1		29.7548		29.7559		29.7567	0.0019		
2	26.7640	26.7603	26.7647	26.7612	26.7641	26.7614	0.0011	0.0044	
3	23.7781	23.7747	23.7788	23.7764	23.7792	23.7764	0.0017	0.0045	
4	20.8107	20.8085	20.8111	20.8089	20.8117	20.8093	0.0010	0.0032	
5	17.8452	17.8435	17.8456	17.8437	17.8461	17.8443	0.0009	0.0026	
6	14.8943	14.8926	14.8958	14.8928	14.8955	14.8930	0.0015	0.0032	
7	12.7080		12.7086		12.7082		0.0006		
Positioning Error (RMS)							0.0013	0.0037	

Finite Element Analysis (FEA)

Before actual assembly, the system was tested virtually with FEA techniques. The main platforms were tested under a load of 1500 kg to obtain a realistic model of the maximum load case for the specifications for this table. All deformation values from the FEA testing are within tens of microns order of magnitude. The theoretical deformation spread for the various positioning plates are presented graphically below in figures 3, 4, and 5.




FIGURE 3. Finite Element Analysis of Top Plate




FIGURE 4. Finite Element Analysis of Second Plate




FIGURE 5. Finite Element Analysis of Third Plate

Discussion of Data

The 3-axis optical platform was completed and tested. The Keyence high accuracy optical micrometer was used to measure positioning of the table over 6 passes while the table was under loading. This test, results shown in Figure 6, included 3 passes in both the positive and negative direction (i.e. upwards and downwards). The root mean square error of the system after traveling in only one direction was 1.3 microns whereas that of the system after traveling one pass in each direction was 3.7 microns. This error was sufficiently small that the table was well within normal operating limits and determined to be fully functional. The shortness of the table also fulfilled the low profile requirements of the SOLARIS lab, and the table was tested under the desired loading requirements to ensure reliability of the platform. Photos of the Assembled 3-axis platform under zero load and under 1500 Kg load, are shown in Figures 7 and 8.




FIGURE 6. Position error of 3-axis platform in vertical axis




FIGURE 7. Photograph of the assembled 3-axis platform under zero load



FIGURE 8. Photograph of the assembled 3-axis platform under 1500 Kg load

CONCLUSIONS

This design proved to be very reliable and meet the project requirements bot for low profile and high load as well as reputability.

REFERENCES

1. J. Kulesza, D. Waterman, A. Deyhim, E. Van Every, "Development of Six Sample Positioning Stages for X-ray diffraction, X-ray Scattering, X-ray Imaging and Tomography Experiments" *The 10th International Conference on Synchrotron Radiation Instrumentation* (Melbourne, Australia, 2009)
2. A. Deyhim, E. Acome, E. Van Every, Y. Wang, T. Yu, "Development of a Nine Axis Custom Design High Precision Motion System for Advance Photon Source" *7th International Conference on Mechanical Engineering Design of Synchrotron Radiation Equipment and Instrumentation* (Shanghai, China, 2012)