### **stryker**



# OPERON

**Accessories catalog** 

Like the OPERON line of surgical tables, our accessories are characterized by robust construction, ergonomic functionality, and a wide range of applications designed to meet the needs of nearly every surgical discipline. Our extensive portfolio range is designed for safe patient positioning and is arranged in this catalog by specialty and application to help ensure you have the best options possible.

#### **Dimensions:**

Accessory dimensions are given in standard measurements. US standard operating room table side rails – 3/8" W x 1  $^{1/8}$ " H. All dimensions are L x W x H or L x W. Standard pads are 2  $^{3/8}$ " thick.

#### **Weight limitations:**

Accessories have weight limitations based on the patient's total weight proportional to the body part being supported. When multiple accessories are used, the patient weight should be limited to the lowest capacity accessory. Known weight limitations are listed with the accessories throughout this catalog.

# Table of contents

#### Standard accessories Clamps 4 1.1 6 1.2 Storage systems Anesthesia aids 6 Restraint straps 8 8 1.5 General **Patient supports** Armboards & hand/arm tables 12 Body supports 14 Leg supports 16 Tabletop extensions 18 **Accessories for surgical applications** 22 3.1 Bariatric 3.2 ENT/Ophthalmology 24 Gynecology/Urology 26 Neuro/Spine 30 3.4 Ortho/Trauma 34 3.5 Vascular/Cardiovascular 38 **Positioning solutions** Patient positioning aids 40

52

Pad accessory reference

#### 1.1 Clamps



#### Radial Setting Clamp OZ 1001041

Accommodates post diameters up to .79'' as well as square profiles up to .79'' square. Fine radial adjustments in  $6^{\circ}$  increments.

Clamp fits standard EU and US side rails.



### Radial Setting Clamp, Simultaneous Locking OZ 1001075

Accommodates post diameters up to .79" as well as square profiles up to .79" square. Fine radial adjustments in 6° increments.

Clamp fits standard EU and US side rails.



### Quick Clamp, Round Post OZ 7500775

Snap locks to the side rail. Fits round rods up to .63" in diameter.

Fits US standard table side rails.



### Rail Clamp, Flat Blade OZ 7500875

Attaches at OR table rail ends or at rail notches. Allows horizontal or vertical installation of accessories. Accommodates flat blade accessories up to  $1" \times .25"$ .

Fits US standard table side rails.



### Quick Clamp, Flat Blade OZ 7516175

Snap lock on mechanism to lock onto side rail at any position. Accommodates flat blade accessories up to  $1'' \times .25''$ .

Fits US standard table side rails.



#### Imaging Extension Rail Clamp OZ 7522575

This clamp fits the profile of our imaging extensions and provides a 6" length of side rail to be used for mounting an armboard or side rail mounted accessory in combination with the respective clamp. Fits OZ 7520175 and OZ 7521275.

With US standard side rails.



#### Inside Rail Clamp OZ 7527675

Allows the use of standard width footboard and other flat blade accessories that are designed for a standard 20" wide table top. Accommodates flat blade accessories up to 1"x .25".

Fits US standard table side rails.

#### 1.2

### Storage systems



### Accessory Stand OZ 1205075

Mobile storage cart. One storage tray, three storage baskets, six side rails, two sleeves for knee crutch storage. No equipment included.

Dimensions: 22.2" x 19.7" x 58"



#### Accessory Cart, Split Legs OZ 7532175

Cart provides storage for split leg sections, single leg plates, clamps, and stirrups.

Dimensions: 21.3" x 21.3" x 54.3"



### Shoulder Chair Cart OZ 7529475

Storage cart for integrated shoulder chair accessory. Lower surface provides additional storage for accessories.

## 1.3

#### Anesthesia aids



### Anesthesia Screen OZ 1001101

Vertical height adjustment of 25.6", .71" post diameter.

Recommended clamp: OZ 1001075



### Flexible Anesthesia Tube Holder OZ 1001200

Single locking knob secures 2 joints simultaneously providing flexible positioning for a variety of anesthesia tube diameters. Accommodates tube sizes of .57"/.87"/.93".

Clamp fits standard EU and US side rails.

## Restraint straps



### Body Restraint Strap, Padded OZ 7500975

Secures patient with two integrated locking buckles. Padded surface at patient contact area.

Dimensions: Strap: 80" x 2" Pad: 18" x 5"

Patient weight capacity: 500lbs

1.5 General



### Side Rail Extender OZ 7503975

Attaches to side rail with integrated clamp to provide additional rail length.

Dimensions: L: 8.3"



#### X-Ray Cassette Tray OZ 7518475

Cassette tray securely holds film and can be inserted from either end of the table. Tray will pass over column for complete body imaging. Integrated handle stiffens during use and is flexible for convenient storage.



# Patient Transfer Boards OZ 7565475 OZ 7565575 OZ 7565675

Designed to be better for patients, easier for nurses, and more convenient for patient transfer. Foldable for easy storage, includes carry-on strap. Short length, especially recommended for orthopaedics, gynecology and urology.

OZ 7565475 long and wide version. Dimensions: 71.5" x 20.5"

OZ 7565575 short and wide version. Dimensions: 32" x 20.5"

OZ 7565675 long and narrow version. Dimensions:  $40'' \times 16.5''$ 

Patient weight capacity: 500lbs

1.5 General



### Foot Control OZ 7518675

Controls height, tilt and Trendelenburg.



## Filler Plates OZ 7519975 OZ 7522675

OZ 7519975 (Perineal Filler) includes removable standard pad. Closes perineal cut-out at seat and back section, when applicable, of standard table top. Radiolucent. Attaches with hook and loop tape.

OZ 7522675 (Imaging Extension Filler) includes removable standard pad.

Compatible with OZ 7520175, OZ 7521275 and OZ 7520375 when attached at the foot end of the table.

Patient weight capacity: 500lbs



#### Kidney Support OZ 7520075

Slides onto kidney elevator plate for mounting. Anatomically curved for comfortable body support. Can be used for lateral patient positioning.

Compatible with all D series tables.



### Footboard/Foot Extension OZ 7523275

Attaches to side rails. Adjustable width to also fit most other tables.

Dimensions: 12" x 24"

Recommended clamps: 2 x OZ 7500875 Patient weight capacity: 300lbs



### X-Ray Grid Encasement OZ 7537475

 $14^{\prime\prime}$  x  $17^{\prime\prime}$  (drop on) grid encasement customized to fit the integrated x-ray channel of the OPERON D series tables.

Works with 14" x 17" film and CR cassettes.

#### 2.1

#### Armboards & hand/arm tables



### Armboard, Ball Joint, One Hand OZ 1251875

Height adjustable armboard with one handed manipulation of tilt and rotation. Integrated clamp.

Dimensions: 23.6" x 6.3" Height adjustment: 3.94"

Pads: OZ 1251800 (concave) OZ 1251810 (flat)

Restraint strap: OZ 1201830 Patient weight capacity: 440lbs



#### Armboard, Multitask OZ 7500475

Ball joint permits a wide range of positioning flexibility. Includes pad and restraint strap. Clamp sold separately.

Dimensions: 20" x 5.5"
Height adjustment: 16.5"
Longitudinal adjustment: 7"
Recommended clamp: OZ 7500775
Patient weight capacity: 500lbs



### Trigger Armboard OZ 7504175

Integrated clamp snap locks onto rail at any position. One hand adjustable. Includes standard pad and restraint strap. Radiolucent.

Dimensions: 26" x 6"

Patient weight capacity: 400lbs



### Universal Arm/Hand Table OZ 7515775

Attaches to side rail via round post clamps. Includes removable support leg and standard pad. Radiolucent.

Dimensions:  $30'' \times 15''$ 

Recommended clamps: 2 x OZ 7500775 Patient weight capacity: 500lbs



### Hand/Arm Table, Hourglass Shape OZ 7516875

Attaches to the side rail. Contoured for better access. Includes removable support leg and standard pad.

Dimensions: 31.1" x 22.8"

Width at table center: 15" Recommended clamps: 2 x OZ 7500775

Patient weight capacity: 500lbs



### Armboard, Bariatric OZ 7519175

Integrated clamp snap locks onto rail at any position. One hand adjustable. Includes standard gel pad and restraint strap. Radiolucent.

Dimensions:  $26'' \times 8''$ 

Patient weight capacity: 1,000lbs

#### 2.1

#### Armboards & hand/arm tables



### Prone Arm Support OZ 7531375

Allows for effective position of the arms to help prevent potential patient injury. The malleable hand rest is designed to enable the surgeon to place the patient's hand in any desired orientation. Clamps sold separately.

Recommended clamp: OZ 7500775 Patient weight capacity: 500lbs



### Trigger Armboard, Carbon Fiber OZ 7547675

Integrated radiolucent carbon fiber armboard with integrated clamp easily attaches to side rail at any location on the table. Standard pad included.

Dimensions: 26" x 6.3" 2" Pad: OZ 7547775

Patient weight capacity: 400lbs

## 2.2 Body supports



### Footboard/Foot Extension OZ 7523275

Attaches to side rails. Adjustable width to also fit most other tables.

Dimensions: 12" x 24"

Recommended clamps: 2 x OZ 7500875 Patient weight capacity: 500lbs



#### 18" Kneeboard OZ 7540075

18" Kneeboard for patient support in kneeling position. Clamps included.

Dimensions: 21" x 18"

Mounts onto table side rails.

Support Roll: 2 x OZ 1003310 (optional)

Recommended clamps per Support Roll: OZ 1001075 or OZ 1001041

Patient weight capacity: 500lbs



### Shoulder Support TrenStop OZ 7531875

Designed to provide support for extreme Trendelenburg positioning. Flexible adjustment for optimal positioning to provide protection against nerve brachial plexus nerve damage. Clamps sold separately.

Recommended clamps: 2 x OZ 7500875 Patient weight capacity: 450lbs



### Lateral Positioning Set OZ 7537375

Stabilizes patients for procedures requiring optimal lateral positioning. Attaches to side rail. Vertically, horizontally, and ball joint adjustable.

Recommended clamps: OZ 7500875 or OZ 7516175

Lateral Supports: OZ 7501075

Dimensions: 8.7" x 3.4"

Lateral Padding: OZ 1609104, round dimensions: Ø: 3.9"

Patient weight capacity: 400lbs

## Leg supports



### Knee Crutches OZ 7501575

Includes hook and loop tape straps for security and foam padding. Rotation range  $360^{\circ}$ . Clamps sold separately.

Recommended clamps: OZ 7500775 Patient weight capacity: 500lbs



### Premium Leg Stirrups, Pair OZ 7501775

Lift assist mechanism provides easy one-hand adjustment for lithotomy positioning. Soft pad and contoured boot comfortably secures patient's foot and lower leg.

Lithotomy movement range +84° to -35°, abduction range +90° to -35°.

Recommended clamps: 2 x OZ 7500875 or 2 x OZ 7516175 Accessory Stand: OZ 7532175

Replacement Pad Set: OZ 7555475 Patient weight capacity: 350lbs



#### Arthroscopy Leg Holder OZ 7511975

Attaches to the side rail with clamp. Integrated ratchet adjustment allows stable leg positioning. Includes two disposable pads.

Recommended clamp: OZ 7500775

Disposable pad sets: OZ 7512075 (12 single-use pads units per set)

Patient weight capacity: 500lbs



#### Split Legs, Pair OZ 7521375

Attaches to OPERON standard table connector. Double pivot allows dual articulation, includes standard pad.

Dimensions: L: 30"

Recommended footboards: OZ 7521875 Patient weight capacity: 600lbs



### Footboard for Split Leg Plates, Pair OZ 7521875

Attaches to side rail with clamps. Provides patient support for extreme Reverse Trendelenburg positioning. Includes removable standard pad.

Dimensions: 10.2" square

Recommended clamps: 2 x OZ 7500775
Patient weight capacity: 600lbs



### Leg Transfer Board OZ 7522475

Locks into table using x-ray cassette channel. Used to support patient's legs before and after lithotomy positioning. Lightweight and rigid. Includes removable standard pad.

Dimensions: 35.4" x 21"

Patient weight capacity: 300lbs



### Leg Transfer, Sling Style OZ 7560175

Lightweight fabric style. Used to support the patient's legs before and after lithotomy positioning. Supports user ergonomics for easy attach and release of the section. Connects to the x-ray channel of the table.

Replacement cover: OZ 7560275 Patient weight capacity: 600lbs

## Leg supports



### Leg Transfer Board, Adjustable OZ 7529075

Panel position slides to allow use on both leg and back side of the seat section for normal and reverse orientation. Used to support the patient's legs before and after lithotomy positioning. Includes standard pad.

Dimensions: 36" x 21"
Attaches to OPERON's standard section connectors.

Patient weight capacity: 750lbs standard table position 350lbs reverse table postion



#### Leg Holder, J Shape OZ 7540275

Attaches to side rail with clamp. Stainless steel tube construction.

Vertical adjustment: 28"- 41" Recommended clamps: OZ 7500775 Patient weight capacity: 300lbs

## 2.4Tabletop extensions



## Filler Plates OZ 7519975 OZ 7522675

OZ 7519975 (Perineal Filler) includes removable standard pad. Closes perineal cut-out at seat and back section, when applicable, of standard table top. Radiolucent. Attaches with hook and loop tape.

OZ 7522675 (Imaging Extension Filler) includes removable standard pad.

Compatible with OZ 7520175, OZ 7521275, and OZ 7520375 when attached at the foot end of the table.

Patient weight capacity: 500lbs



### Imaging Extension Plate, 36" Short, Carbon Fiber OZ 7520175

Attaches in place of back plate or leg plate to allow  $360^{\circ}$  of radiolucent imaging / 3D imaging.

Dimensions: 36"x 21"

Use Filler Plate OZ 7522675 if attached in place of leg section.

Patient weight capacity: 400lbs

## 2.4 Tabletop extensions



### Imaging Extension Plate, 44" Long, Carbon Fiber OZ 7521275

Attaches in place of back plate or leg plate to allow  $360^\circ$  of radiolucent imaging. Leg support included.

Dimensions: 44" x 21"

Use Filler Plate OZ 7522675 if attached in place of leg section. Use Filler Plate OZ 7519975 if attached to the seat section in place of the back section.

Patient weight capacity: 750lbs with support leg



#### Split Legs, Pair OZ 7521375

Attaches to OPERON standard table connector. Double pivot allows dual articulation, includes standard pad.

Dimensions: L: 30"

Recommended footboards: OZ 7521875 Patient weight capacity: 600lbs



### Footboard for Split Leg Plates, Pair OZ 7521875

Attaches to side rail with clamps. Provides patient support for extreme Reverse Trendelenburg positioning. Includes removable standard pad.

Dimensions: 10.2" square

Recommended clamps: 2 x OZ 7500775 Patient weight capacity: 600lbs



### Pediatric Table, Carbon Fiber OZ 7520375

Attaches to table connectors in place of back section. Integrated side rails. Scalloped shape to permit close access to patient. Use Filler Plate OZ 7522675 if attached in place of leg section. Includes standard pad.

Dimensions:

36" x 14.25" (rail to rail) 36" x 12.25" (without rails)

Patient weight capacity: 75lbs

#### 3.1 Bariatric



### Tabletop Width Extender OZ 7501275

Extends table surface for larger patients. Clips to the side rail with integrated clamp. Includes removable pad. Integrated side rail accommodates attachment of additional accessories.

Dimensions: 4.1" x 20" Patient weight capacity: 400lbs



#### Bariatric Leg Stirrups, Pair OZ 7518975

Lift assist mechanism provides easy one-hand adjustment for lithotomy positioning. Soft pad and contoured boot designed for larger/ heavier patients. Comfortably secures patient's foot and lower leg.

Lithotomy movement range  $+84^{\circ}$  to  $-33^{\circ}$ Abduction range  $+90^{\circ}$  to  $-35^{\circ}$ 

Clamps sold separately.
Replacement pad set: OZ 7555575
Recommended clamps: OZ 7500875
Patient weight capacity: 720lbs



#### Leg & Body Strap, Bariatric OZ 7519075

Secures patient with four integrated locking buckles. Safe and convenient fixation through additional extra wide and long padding.

Dimensions: Strap: 80" x 2" Pad: 30" x 11.8"

Patient weight capacity: 1,000lbs



### Armboard, Bariatric OZ 7519175

Integrated clamp snap locks onto rail at any position. One hand adjustable. Includes standard gel pad and restraint strap. Radiolucent.

Dimensions: 26" x 8"
Patient weight capacity: 1,000lbs




#### Nissen Straps, Bariatric, Pair OZ 7520575

Bariatric Nissen straps offer the support larger patients need during steep Reverse Trendelenburg positioning. Attaches via integrated clamps to the side rail. Oversized pad offers comfortable support.

Dimensions: 72.3" x 2"

Pad: 33" x 7"

Patient weight capacity: 1,000lbs



### Bariatric Width Extenders set, Short Leg, with $2^{3/8}$ " pad OZ 7556975

Extends table surface to accommodate larger patients. Mounts onto side rail and secures via integrated clamps. Includes standard removable pad (2  $^{3/8}$ "). Fully radiolucent carbon fiber also includes standard side rails to attach accessories.

Patient weight capacity: 1,000lbs



### Bariatric Width Extenders set, Long Leg, with 2<sup>3/8</sup>" pad OZ 7557175

Extends table surface to accommodate larger patients. Mounts onto side rail and secures via integrated clamps. Includes standard removable pad (2  $^{3/8}$ "). Fully radiolucent carbon fiber also includes standard side rails to attach accessories.

Patient weight capacity: 1,000lbs


3.1 Bariatric



### Bariatric Width Extenders set, Short Leg, without pad OZ 7557075

Extends table surface to accommodate larger patients. Mounts onto side rail and secures via integrated clamps. Does not include standard pads. Fully radiolucent carbon fiber also includes standard side rails to attach accessories.

Patient weight capacity: 1,000lbs



### Bariatric Width Extenders set, Long Leg, without pad OZ 7557275

Extends table surface to accommodate larger patients. Mounts onto side rail and secures via integrated clamps. Does not include standard pads. Fully radiolucent carbon fiber also includes standard side rails to attach accessories.

Patient weight capacity: 1,000lbs

## 3.2 ENT/Ophthalmology



#### Horseshoe Headrest OZ 1001351

Height and width adjustable. Ball joint facilitates head positioning for ENT, oral and plastic surgery. Use in combination with wrist support OZ 1001352.

Mounts using neuro adapter part numbers OZ 7515575 or OZ 7520675.



### ENT Horseshoe Headrest, Adjustable OZ 1001400

Dual joint adjustability for ENT, plastic, and ophthalmic applications. Single locking lever for both joints provides ease of positioning.

Mounts using neuro adapters OZ 7515575 or OZ 7520675.



### Wrist Support for Horseshoe Headrest OZ 1001352

Provides hand resting surface and support used in combination with the Horseshoe Headrest OZ 1001351.

Mounts to Horseshoe Headrest OZ 1001351.



### Headrest, Multi Adjustable OZ 1001410

Multi adjustable head support for ENT, plastic and ophthalmic surgery with dual joint operation each secured by a single lever.

Mounts using neuro adapters OZ 7515575 or OZ 7520675.

3.2 ENT/Ophthalmology



### Headrest, Narrow OZ 7517375

Adaptable with appropriate neuro adapter (not included.) Non-articulating headrest contoured to provide optimum access for ENT and Ophthalmology. Includes standard removable pad.

Length: 15.8"

Recommended neuro adapter: OZ 7515575

Patient weight capacity: 500lbs


## 3.3 Gynecology/Urology



### Elbow Support, Piece OZ 1252675

Elbow support for the surgeons hand/forearm during procedures. Multiple adjustments to accommodate positions. Integrated clamp.

Interface: Mounts onto table side rails.



#### Attachment Strap for Drainage Bags OZ 7501875

Flexible band provides excellent access while providing superior fluid containment via drainage bags. Integrated clamp.

Use in combination with OZ 7502075.



#### Drainage Bags OZ 7502075

Non-sterile fine mesh drainage bags for gynecology/urology procedures (10 per box).

Compatible with all common suction systems for liquid collection. Use in combination with OZ 7501875.



### Instrument Surface OZ 7504075

Attaches to the side rail with clamp. Provides surface for instruments and supplies for gynecology/urology procedures.

Dimensions: 10" x 24"

Recommended clamp: OZ 7500775

3.3
Gynecology/Urology



### Kidney Support OZ 7520075

Slides onto kidney elevator plate for mounting. Anatomically curved for comfortable body support. Can be used for lateral patient positioning.



#### Leg Transfer Board OZ 7522475

Locks into table using x-ray cassette channel. Used to support patient's legs before and after lithotomy positioning. Lightweight and rigid. Includes removable standard pad.

Dimensions: 35.4" x 21"

Patient weight capacity: 300lbs



### Leg Transfer, Sling Style OZ 7560175

Lightweight fabric style. Used to support the patient's legs before and after lithotomy positioning. Supports user ergonomics for easy attach and release of the section. Connects to the x-ray channel of the table.

Replacement cover: OZ 7560275 Patient weight capacity: 600lbs



### Leg Transfer Board, Adjustable OZ 7529075

Panel position slides to allow use on both leg and back side of the seat section for normal and reverse orientation. Used to support the patient's legs before and after lithotomy positioning. Includes standard pad.

Dimensions: 36" x 21"

Attaches to OPERON's standard section connectors.

Patient weight capacity: 750lbs standard table position 350lbs reverse table postion



#### Leg Holder, J Shape OZ 7540275

Attaches to side rail with clamp. Stainless steel tube construction.

Vertical adjustment: 28" - 41" Recommended clamps: OZ 7500775 Patient weight capacity: 300lbs



### **TUR Module OZ 7519675**

Provides patient support for extreme Reverse Trendelenburg and extended imaging access for gynecology and urology cases. Equipped with radiolucent support plate. Cannot be used with split leg plates. Includes standard gel pad.

Dimensions: L: 11.9"

Attaches to OPERON's standard section connectors in place of leg plate.

Patient weight capacity: 800lbs

## 3.4 Neuro/Spine



### Knee-Elbow Positioning Device OZ 1005275

For knee-elbow position during spinal surgery or proctology. 7.9" height adjustment via hand crank. Attaches to OPERON's standard section connectors.

Support Roll:  $2 \times OZ$  1003310 (optional) Recommended clamps per support roll: OZ 1001075 or OZ 1001041 Patient weight capacity: 440lbs



#### Neuro Adapter OZ 7515575

Attaches to side rail. Accepts commonly used Mayfield, Horseshoe Headrest and specialty headrests. Adjusts laterally to fit nearly all tables.

Attaches to OPERON's standard section connectors. Ideal for mounting Mayfield and Horseshoe Headrest. Patient weight capacity: 500lbs



### Pediatric Table, Carbon Fiber OZ 7520375

Attaches to table connectors in place of back section. Integrated side rails. Scalloped shape to permit close access to patient. Use Filler Plate OZ 7522675 if attached in place of leg section. Includes standard pad.

Dimensions: 36" x 14.25" (rail to rail) 36" x 12.25" (without rails)

Patient weight capacity: 75lbs



### Neuro Adapter, Low Profile OZ 7520675

Attaches to OPERON's standard section connectors. Ideal for mounting Mayfield, Horseshoe Headrest and special headrests.

Patient weight capacity: 1,000lbs



### Spine Frame OZ 7529375

Attaches easily to any OPERON surgical table for versatile prone, lateral or supine positioning to facilitate spinal surgery. Allows intraoperative flex eye and airway visibility and access. Unrestricted C-Arm access Radiolucent spars. Small storage footprint. Accessories store on the frame. Supine top set option to convert the open spine frame to a solid image surface. Wingsets allow abdomen to hang free.

Dimensions: 48" x 17"

Patient weight capacity: 500lbs



### Prone Arm Support OZ 7531375

Allows for effective position of the arms to help reduce potential patient injury. The malleable hand rest enables the surgeon to place the patient's hand in any desired orientation. Clamps sold separately.

Recommended clamp: OZ 7500775 Patient weight capacity: 500lbs

3.4
Neuro/Spine



#### 18" Kneeboard OZ 7540075

18" Kneeboard for patient support in kneeling position. Clamps included.

Dimensions: 21" x 18" Mounts onto table side rails.

Support Roll: 2 x OZ 1003310 (optional)

Recommended clamps per Support Roll: OZ 1001075 or OZ 1001041 Patient weight capacity: 500lbs



### Arm Cradle, Disposable Pad OZ 7547275

Disposable arm supports. 6 pairs per case.



#### Spine Frame (Wilson Style) OZ 7547175

Provides a stable and convenient platform for spine surgeries such as laminectomies, decompressions and disc surgery. Intraoperatively adjustable to the desired lordosis to get optimal surgical site exposure. Tapered, adjustable viscoelastic pads provide support and comfort. C-arm compatible.

Comes complete with transport storage dolly and crank handle, disposable headrest, and two arm cradles.



### Headrest, Disposable Pad OZ 7547375

Disposable slotted headrests. 6 per case.



Pad Set, Spine Frame OZ 7547475

Wilson style spine frame replacement pad set.



Mask, Disposable OZ 7547575

For use with integrated shoulder chairs OZ 7519475 and OZ 7523075.  $12\ masks$  per box.

## 3.5 Ortho/Trauma



### **Supine Tops OZ 7530475**

Radiolucent 4 piece set that fits spine frame OZ 7539375.

Length: 12"



### Prone Leg Pad for Spine Frame System OZ 7530675

Provides comfortable, anatomically correct positioning of the legs while prone patient is engaged with spine frame.

Dimensions: 42" x 5"



### Extension for Upper Arm (Weinberger) OZ 1001450

Upper arm extension for intramedullary nailing of the upper arm/elbow surgery. Set consisting of: Upper armboard with traction unit, fine linear adjustment for the patient's arm length and applied traction forces. Upper arm countertraction post adjustable via clamps. Includes pad.

Radiolucent Weinberger handcuff with integrated clamp. Attaches to the upper armboard with traction unit. Radiolucent. Attaches to side rail via clamps.

Recommended clamps: 2 x OZ 7500775



#### Arthroscopy Leg Holder OZ 1002540 OZ 1002560

OZ 1002540 normal version OZ 1002560 XL version. Adjustable horizontally with hand crank and radially with clamp for secure positioning of the leg. Includes three removable pads.

Width adjustment 5.5" to 10.5" XL version also fits legs of larger/heavier patients.

#### Dimensions:

OZ 1002540: (2 x pads) 7.5" x 3.3" (1 x pad) 7.1" x 5.5" OZ 1002560: (2 x pads) 10.55"x 2.4" (1 x pad) 7.1" x 5.5" Recommended clamps: OZ 1001075 or OZ 1001041



#### Armboard, Multitask OZ 7500475

Ball joint permits a wide range of positioning flexibility. Includes pad and restraint strap. Clamp sold separately.

Dimensions: 20" x 5.5" Height adjustment: 16.5" Longitudinal adjustment: 7" Recommended clamp: OZ 7500775 Patient weight capacity: 500lbs



### Arthroscopy Leg Holder OZ 7511975

Attaches to the side rail with clamp. Integrated ratchet adjustment allows stable leg positioning. Includes two disposable pads.

Recommended clamp: OZ 7500775

Disposable pad sets: OZ 7512075 (12 single-use pads units per set)

Patient weight capacity: 500lbs

## 3.5 Ortho/Trauma



### Chin Straps Single-Use for Shoulder Chair OZ 7518375

Single-use soft padded strap with additional chin layer securely restrains patient's head. 12 sets per case.

Compatible with shoulder chairs OZ 7523075 and OZ 7519475.




### Integrated Shoulder Chair, D8XX OZ 7519475

Mounts to the head end of the seat section on D model tables. 3-section back plate and ball joint. Adjustable head holder. Two removable side sections provide surgical access. Includes lateral brace to provide support. Rear-mounted side rail for the mounting of shoulder support devices. Attaches to OPERON's standard seat section connectors.

Chin straps: OZ 7518375

Disposable mask: OZ 7547575 (12 sets per box) Replacement Lateral Brace: OZ 7541675 Patient weight capacity: 500lbs



#### Integrated Shoulder Chair, D7XX OZ 7523075

Mounts on the leg-end of the seat section on D7XX model tables. 3-section back plate and ball joint. Adjustable head holder. Two removable side sections provide surgical access. Includes lateral brace to provide support. Rear-mounted side rail for the mounting of shoulder support devices. Attaches to OPERON's standard seat section connectors.

Disposable mask: OZ 7547575 (12 sets per box)

Chin straps: OZ 7518375

Replacement Lateral Brace: OZ 7541675
Patient weight capacity: 500lbs



### Orthopaedic Extension Unit, Carbon Fiber OZ 7532975

Converts your OPERON surgical table from a general purpose table into an orthopaedic specialty table. The extension unit is designed to be used for hip and lower extremity procedures such as: hip pinning, ender pinning, intramedullary nailing of femur, tibia and fibula surgery.

The extension is a simple and secure attachment to the table with the help of the included cart which enables the user to remove and store the extension and accessories when not in use. Abductor bars and most accessories are constructed of high strength, lightweight, and radiolucent carbon fiber material.

Patient weight capacity: 400lbs with support leg 350lbs without leg support

3.6

#### Vascular/Cardiovascular



## Filler Plates OZ 7519975 OZ 7522675

OZ 7519975 (Perineal Filler) includes removable standard pad. Closes perineal cut-out at seat and back section, when applicable, of standard table top. Radiolucent. Attaches with hook and loop tape.

OZ 7522675 (Imaging Extension Filler) includes removable standard pad.

Compatible with OZ 7520175, OZ 7521275 and OZ 7520375 when attached at the foot end of the table.

Patient weight capacity: 500lbs



### Imaging Extension Plate, 36" Short, Carbon Fiber OZ 7520175

Attaches in place of back plate or leg plate to allow 360° of radiolucent imaging / 3D imaging.

Dimension: 36"x 21"

Use Filler Plate OZ 7522675 if attached in place of leg section.

Patient weight capacity: 400lbs



### Imaging Extension Plate, 44" Long, Carbon Fiber OZ 7521275

Attaches in place of back plate or leg plate to allow  $360^\circ$  of radiolucent imaging. Leg support included.

Dimensions: 44" x 21"

Use Filler Plate OZ 7522675 if attached in place of leg section.

Use Filler Plate OZ 7519975 if attached to the seat section in place of the back section.

Patient weight capacity: 750lbs with support leg.

#### 4.1

#### Patient positioning aids



### Support Roll, Horizontal OZ 1001691

Multiuse positioner adjusts vertically and horizontally with additional joint for pad rotation. Multiple uses such as upper arm positioning for fractures, or thoracic support in lateral positions. Includes removable pad.

Dimensions: L x  $\emptyset$ : 10.8" x 5.3" Attaches via clamp to the side rail. Recommended clamp: OZ 1001075



#### Tibia Support / Knee Rest OZ 1003310

Used to support the patient's leg during knee/ tibia applications. Clamp sold separately.

Dimensions: L x Ø: 10.8" x 5.3" Attaches via clamp to the side rail. Recommended clamps: OZ 1001075 or OZ 1001041



### Knee-Elbow Positioning Device OZ 1005275

For knee-elbow position during spinal surgery or proctology. 7.9" height adjustment via hand crank. Attaches to OPERON's standard section connectors.

Support Roll:  $2 \times OZ$  1003310 (optional) Recommended clamps per support roll: OZ 1001075 or OZ 1001041 Patient weight capacity: 440lbs



### Anesthesia Cushion OZ 2200010

Stabilizes the position of the head in supine position. Equipped with electrically conductive top cover.

Dimensions: 16.53" x 18.9" x 3.93



#### Half-Moon Cushion OZ 2200020

Half-moon/semi-circular cushion. Suitable for breast, pelvis and leg support/elevation. Equipped with electrically conductive top cover.

Dimensions: 19.7" x 8.7" x 4.3"



### Leg Rest Cushion OZ 2200030

Designed to stabilize leg in a comfortable position. Cushion for one leg. Equipped with electrically conductive top cover.

Dimensions: 23.6" x 8.6" x 7.1"



### Tunnel Cushion OZ 2200040

Designed to protect the legs from pressure damage in lateral position. Equipped with electrically conductive top cover.

Dimensions: 23.6" x 15.3" x 7.1"

4.1

#### Patient positioning aids



### Lateral Positioning Cushion OZ 2200050

Wedge-shaped cushion with cut-out designed to reduce pressure of the arm beneath the patient. Mainly suitable for entire lateral patient positions. Equipped with electrically conductive top cover.

Dimensions: 27.6" x 19.7" x 6.3"



### Spinal Cushion OZ 2200060

Foam cushion with cut-out for spinal interventions/prone positioning. X-ray- translucent. Equipped with electrically conductive top cover.

Dimensions: 33" x 18" x 3.93" / 14" x 7.87" (cut-out)



### Vacuum Pump with Foot Stand OZ 2200110

Manual foot pump for air suction of vacuum mattress OZ 22001X0.

Dimensions: L: 23.6"



#### Vacuum Repair Kit OZ 2200120

Repair kit to close small holes of the vacuum mattress surface.



#### Non-Slip Sheet Roll Tabletop OZ 2200130

Designed to prevent vacuum mattress from sliding on the tabletop in extreme positions.

Single roll 15.7" x 35.4" to tailor on desired length.



#### Universal Vacuum Mattress OZ 2200140

Single hand operable valve. X-ray permeable and slip-resistant. Universal size. Contours to the shape of the patient when suction is applied. Utilizes vacuum pump OZ 2200110.

Dimensions: 78.7" x 29.2"

Use with non-slip foil for extreme angle positioning: OZ 2200130



### Torso Section Vacuum Mattress with Cutout OZ 2200150

Single hand operable valve. X-ray permeable and slip-resistant. Torso section. Contours to the shape of the patient when suction is applied. Ideal for gynecology, urology, cholecystectomies, and robotic surgical procedures.

Utilizes vacuum pump OZ 2200110.

Dimensions: 39.4" x 29.2"

Non-slip foil for extreme positioning: OZ 2200130

#### 4.1

#### Patient positioning aids



### Universal Split Leg Vacuum Mattress OZ 2200160

Single hand operable valve. X-ray permeable and slip-resistant. Use with split leg configured tables. Contours to the shape of the patient when suction is applied.

Utilizes vacuum pump OZ 2200110. Dimensions/Interface: 78.7" x 29.5"

Non-slip foil for extreme positioning: OZ 2200130



#### Lateral Support, Rectangular OZ 7501075

Vertically, horizontally and ball-joint adjustable. Attaches to side rail via flat blade clamp, pad included.

Recommended clamps: OZ 7500875 or OZ 7516175



### Shoulder Support, Pair OZ 7517075

Attaches to side rail with clamps. Provides patient support for extreme Trendelenburg positioning. Flexible adjustment for optimal positioning. Includes removable pads.

Dimensions: 7.1" x 5.5"

Recommended clamps: 2 x OZ 7500875 Patient Weight Capacity: 500lbs



### Shoulder Support TrenStop OZ 7531875

Designed to provide support for extreme Trendelenburg positioning. Flexible adjustment for optimal positioning to provide protection against nerve brachial plexus nerve damage. Clamps sold separately.

Recommended clamps: 2 x OZ 7500875 Patient weight capacity: 450lbs



### Lateral Positioning Set OZ 7537375

Stabilizes patients for procedures requiring optimal lateral positioning. Attaches to side rail. Vertically, horizontally, and ball joint adjustable.

Recommended clamps:
OZ 7500875 or OZ 7516175
OZ 7501075 Lateral Supports Dimensions: 8.7" x 3.4"
OZ 1609104 Lateral Padding, Round Dimensions: Ø: 3.9"

Patient weight capacity: 400lbs



#### Gel Head Donut OZ 7538575 OZ 7538675

OZ 7538575, for adults OZ 7538675, for adolescents/children Protects and stabilizes the head for supine positioning.

OZ 7538575
Dimensions:
External Ø x Internal Ø x H:
Adults: 8" x 3" x 2"
Adolescents/children: 6" x 2" x 1"

4.1

#### Patient positioning aids



#### Gel Horseshoe Head Donut OZ 7538775 OZ 7546375

Used for various positioning to protect and stabilize the head. Cut-out allows anesthetist access to tubes during operation.

OZ 7538775 for adults OZ 7546375 for adolescents/children

Dimensions: External  $\emptyset$  x Internal  $\emptyset$  x H: 8" x 3" x 2"



#### Gel Heel Pad OZ 7538875


This pad was designed to effectively protect bones on the heels and elbows.



#### Gel Heel Elevators OZ 7538975

This pad was designed to effectively protect bones on the heels and elbows.

Dimensions: 6" x 2" x 2"



#### Gel Round End Armboard Pad OZ 7539075

Universal gel pad for standard size armboards.

Dimensions: 6" x 26" x .24"



### Gel Oval Ulnar Nerve Protector OZ 7539175

Designed to help prevent pressure damage to the ulnar nerve. Compact design ensures full access for anesthesia.

Dimensions:  $13'' \times 6'' \times .5''$ 



#### **Gel Flat Bottom Positioning Rolls**

OZ 7539275 OZ 7539375

OZ 7539475 OZ 7539575

For lifting, stabilizing and positioning various body parts such as the breast, pelvis, knee, heel, and ankle.

#### Dimensions:


OZ 7539275: 13" x 5" x 4"

OZ 7539375: 17" x 5" x 4"

OZ 7539475: 21" x 5" x 4"

OZ 7539575: 21" x 7" x 6"

4.1
Patient positioning aids



#### Gel Standard Positioners OZ 7539675 OZ 7539775

Universal positioning aid. Used for lateral position, or as a support under the heels, knees or shoulder.

Dimensions:

OZ 7539675 10" x 5" x 2" OZ 7539775: 15" x 5" x 2.2"



#### Gel Contoured Armboard Pads OZ 7539875 OZ 7540175

Protects biceps, elbows and forearm. The concave shape reduces potential for the arm rolling down. Ideal for supine and lateral.

Dimensions:

OZ 7539875: 16" x 5.8" x 2.2" OZ 7540175: 20" x 5.8" x 2.2"



### Gel Horseshoe Head Donut, Pediatric OZ 7546375

Used for various positioning to protect and stabilize the head. Cut-out allows excellent anesthetist access to tubes during operation.

Dimensions:

External Ø x Internal Ø x H:

6" x 2.25" x 1



#### Gel Prone Headrest OZ 7546475

Headrest designed for anatomically correct positioning and stabilization of the patient's head during prone positioning. Allows tube access from both sides.

Dimensions: 11" x 9" x 5.9"



#### Gel Knee Crutch Pad OZ 7546575

Comfort and protection for leg positioning (e.g. Goepel leg supports). Mainly for general surgery, gynecology, and urology.

Dimensions: 9" x 11" x .38"



#### Gel Arm Protector OZ 7546675

Universal gel support for armboards to guard against nerve damage. Gentle positioning of the patient's arm.

Dimensions: 15" x 5.9" x .25"

4.1

#### Patient positioning aids



#### Gel Operating Table Pad OZ 7546775 OZ 7546875

Comfort and protection over the entire length of the OR table.
Usable in multiple ways, e.g. folded or rolled. Designed for pressure distribution for prone, supine and lithotomy as well as steep
Trendelenburg positions.

OZ 7546875 with cut-out for lithotomy and steep Trendelenburg positions. Ideal for gynecology and urology.

Dimensions: OZ 7546775: 46" x 20" x .5" OZ 7546875: 40" x 20" x .5"



### Gel Perineal Post Cushion OZ 7546975

Additional padding of the perineal post, especially recommended for cases where high traction forces are executed, e.g. arthroscopies. Universal adaptable on perineal post.

Dimensions:  $7'' \times 5.1'' .5''$ 



### **Gel Insert OZ 7547075**

Gel insert for knee arthroscopy device.

Dimensions: 37" x 3.94" x .38"

#### Pad accessory

## Pad accessory reference

Description	Accessory with pad	Accessory without pad	3" pad (hybridgel)	2 <sup>3/8</sup> " pad (foam)	4" pad (foam)	2 <sup>3/8</sup> " (gel pad)
Trigger armboard	OZ 7504175	None	OZ 7548275	OZ 7509375	OZ 7533175	OZ 7519275
Trigger armboard, carbon fiber	OZ 7547675	None	None	OZ 7547775	None	None
Mulit-task armboard	OZ 7500475	None	None	None	None	None
TUR module	OZ 7519675 (with gel pad)	OZ 1705675	OZ 7548475	OZ 7534675	OZ 7534075	None
Pediatric table extension	OZ 7520375	OZ 7520275 + rails OZ 7546175 side + OZ 7546275 front	OZ 7548675	OZ 7515175	OZ 7527475	None
Leg transfer board adjustable	OZ 7529075	OZ 7522175	OZ 7548875	OZ 7520775	OZ 7526775	None
Imaging extension (21" x 36")	OZ 7520175	OZ 7520875	OZ 7548875	OZ 7520775	OZ 7526775	None
Imaging extension (21" x 44")	OZ 7521275	OZ 7521075	OZ 7549075	OZ 7520975	OZ 7526875	None
Perineal filler	OZ 7519975	OZ 7519775	OZ 7549275	OZ 7519875	OZ 7526475	None
Imaging extension filler	OZ 7522675	OZ 7522875	OZ 7549475	OZ 7522775	OZ 7526675	None
Table width extender	OZ 7501275	None	OZ 7549675	OZ 7544075	OZ 7533075	None
Table width extender, carbon fiber, short leg	OZ 7556975	OZ 7557075	None	OZ 7557875	OZ 7569475	None
Table width extender, carbon fiber, long leg	OZ 7557175	OZ 7557275	None	OZ 7557975	OZ 7569575	
Right split leg	OZ 7521375, pair, with pads	OZ 7521575	OZ 7549875	OZ 7521775	OZ 7543975	None
Left split leg	OZ 7521375, pair, with pads	OZ 7521475	OZ 7550075	OZ 7521675	OZ 7543875	None
Leg transfer board	OZ 7522475	None	OZ 7550275	OZ 7543575	OZ 7533975	None
Light duty transfer, short version, 27"	OZ 7555375	None	None	None	None	None
Bariatric armboard	OZ 7519175 (with gel pad)	OZ 7519375	OZ 7552675	None	OZ 7534175	OZ 7519275
18" Knee board	OZ 7540075	None	None	None	None	

Pad accessory

## Pad accessory reference

Description	Accessory with pad	Accessory without pad	3" pad (hybridgel)	2 <sup>3/8</sup> " pad (foam)	4" pad (foam)	2 <sup>3/8</sup> " (gel pad)
Foot board, split leg, pair	OZ 7521875	None	OZ 7552875	OZ 7522075	None	None
Adj foot extention	OZ 7523275	None	OZ 7553075	OZ 7542175	OZ 7543775	None
Narrow headrest	OZ 7517375	OZ 7516975	OZ 7554275	OZ 7517275	OZ 7527375	None
Arm board	None	OZ 7547675	None	OZ 7547775	None	None
Transfer board, sling, x-ray channel	OZ 7560175, with cover	None	Cover only, OZ 7560275	None	None	None
Vascular arm board, rotational	None	OZ 7548075	None	OZ 7548175	None	None
Leg stirrups, adult	OZ 7501775	None	Boot only, OZ 7555475	None	None	None
Universal arm/hand table	OZ 7515775	OZ 7512875	None	OZ 7515275	None	None
Hour glass arm/hand table	OZ 7516875	OZ 7516675	None	OZ 7516775	None	None
Shoulder support, trenstop	OZ 7531875	None	None	None	None	None

#### Pad accessory

### Pad accessory **reference**

#### OPERON table sections D850/D860/D760/D770

Description	Accessory with pad	3" pad (hybridgel)	2 <sup>3/8</sup> " pad (foam)	4" pad (foam)	2 <sup>3/8</sup> " pad (gel pad)
Short leg section	OZ 1705775(CFK), OZ 7505875 (Phenolic)	OZ 7550675	OZ 7518175	OZ 7528175	OZ 7528275
Long leg section	OZ 1704275(CFK), OZ 1704375 (Phenolic)	OZ 7551675	OZ 7524475	OZ 7525175	None
Seat with cut outs	87120	OZ 7550875	OZ 7524275 (with cut), OZ 7524375 (without cut)	OZ 7524975 (with cut), OZ 7525075 (without cut)	OZ 7528975 (without cut)
Split back-pair	OZ 7554775 (Phenolic), OZ 7554875 (CFK)	OZ 7555775	OZ 7558975	OZ 7559875	None
Back	OZ 1705075 (Phenolic), OZ 1705575 (CFK)	OZ 7551075	OZ 7524175	OZ 7524875	None
Head	OZ 1706275 (CFK), OZ 1706575 (Phenolic)	OZ 7551275	OZ 7524075	OZ 7524775	None
Short split leg, left	OZ 7554975 (Phenolic), OZ 7555075 (CFK)	OZ 7559375	OZ 7559175	OZ 7559775	None
Short split leg, right	OZ 7554975 (Phenolic), OZ 7555075 (CFK)	OZ 7559275	OZ 7559075	OZ 7559675	None

#### Shoulder chair

Description	Accessory with pad	Accessory without pad	3" pad (hybridgel)	2 <sup>3/8</sup> " pad (foam)	4" pad (foam)
Shoulder chair complete set	OZ 7519475, D8XX OZ 7523075,D7XX	OZ 7562675, D8XX OZ 7564375, D7XX	Pad set, OZ 7551875	Pad set, OZ 7561475	Pad set, OZ 7534575
Shoulder chair, left	N/A	N/A	OZ 7552075	OZ 7561675	OZ 7534275
Shoulder chair, right	N/A	N/A	OZ 7552275	OZ 7561875	OZ 7534375
Shoulder chair, center	N/A	N/A	OZ 7552475	OZ 7562075	OZ 7534475

#### Table width extenders

Description	Accessory with pad	Accessory without pad	2 <sup>3/8</sup> " pad (foam)	4" pad (foam)
Table width extenders, short leg section	OZ 7556975	OZ 7557075	Pad set, OZ 7557875	Pad set, OZ 7569475
Table width extenders, long leg section	OZ 7557175	OZ 7557275	Pad set, OZ 7557975	Pad set, OZ 7569575

<sup>\*\*</sup>Contact your Stryker Sales Representative or Customer Service for individual pad part numbers.

#### Table pad sets for D7XX and D8XX (unless otherwise noted)

2 <sup>3/8"</sup> pad set, short leg section, without perineal cut, D850/D760 OZ 7528475	2 <sup>3/8</sup> " pad set, long leg section, without perineal cut, D850/D760 OZ 7523975	2 <sup>3/8</sup> " pad set, long leg section, with perineal cut OZ 7523875	2 <sup>3/8"</sup> pad set, short leg section, with perineal cut OZ 7528375	2 <sup>3/8</sup> " gel pad, without cut, D850/ D760 OZ 7528875	2 <sup>3/8</sup> " gel pad with cut, D850/ D760 OZ 7528775
4" pad set, short leg section, with perineal cut OZ 7528575	4" pad set, short leg section, without perineal cut OZ 7528675	4" pad set, long leg section, with perineal cut	4" pad set, long leg section, without perineal cut OZ 7524675		
B810 2 <sup>3/8</sup> " pad set with cut OZ 7525275	B810 2 <sup>3/8</sup> " pad set without cut OZ 7525375	B810 4" pad set with cut OZ 7525875	B810 4" pad set without cut OZ 7525975	B810 2 <sup>3/8</sup> " pad set without cut OZ 7525375	

#### Index

0

OZ 1001041	Radial Setting Clamp	4	OZ 7501775	Stirrups, Premium, Pair	16
OZ 1001075	Radial Setting Clamp, Simultaneous Locking	4	OZ 7501875	Drainage Bag, Attachment Strap	26
OZ 1001101	Anesthesia Screen	6	OZ 7502075	Drainage Bags	27
OZ 1001200	Anesthesia Tube Holder, Flexible	7	OZ 7503975	Rail Extender (Side)	8
OZ 1001351	Headrest, Horseshoe	24	OZ 7504075	Instrument Surface	27
OZ 1001352	Wrist Support for Horseshoe Headrest	24	OZ 7504175	Trigger Armboard	12
OZ 1001400	Headrest, Horseshoe, ENT, Adjustable	25	OZ 7511975	Leg Holder, Arthroscopy	16, 35
OZ 1001410	Headrest, Multi, Adjustable	25	OZ 7515575	Neuro Adapter	30
OZ 1001450	Extension for Upper Arm (Weinberger)	34	OZ 7515775	Universal Arm/Hand Table	13
OZ 1001691	Support Roll, Horizontal	40	OZ 7516175	Ouick Clamp, Flat Blade	5
OZ 1002540	Leg Holder, Arthroscopy	35	OZ 7516875	Hand/Arm Table, Hourglass Shape	13
OZ 1002560	Leg Holder, Arthroscopy	35	OZ 7517075	Shoulder Support, Pair	44
OZ 1003310	Tibia Support/Knee Rest	40	OZ 7517375	Headrest, Narrow	26
OZ 1005275	Knee-Elbow Positioning Device	30, 40	OZ 7518375	Chin Straps Single-Use for Shoulder Chair	36
OZ 1205075	Accessory Stand	6	OZ 7518475	X-Ray Cassette Tray	8
OZ 1251875	Armboard, Ball Joint, One Hand	12	OZ 7518675	Foot Control	10
OZ 1252675	Elbow Support, Piece	26	OZ 7518975	Leg Stirrups, Bariatric, Pair	22
OZ 7519675	TUR Module	29	OZ 7519075	Strap, Leg & Body, Bariatric	22
OZ 2200010	Cushion, Anesthesia	41	OZ 7519175	Armboard, Bariatric	13, 23
OZ 2200020	Cushion, Half Moon	41	OZ 7519475	Integrated Shoulder Chair, D 8XX	36
OZ 2200030	Cushion, Leg Rest	41	OZ 7519975	Filler Plate	10, 18, 38
OZ 2200040	Cushion, Tunnel	41	OZ 7520075	Kidney Support	11, 28
OZ 2200050	Cushion, Lateral Positioning	42	OZ 7520175	Extension Plate, Imaging, Short, Carbon Fiber	19, 38
OZ 2200060	Cushion, Spinal	42	OZ 7520575	Nissen Straps, Bariatric, Pair	23
OZ 2200110	Vacuum Pump with Foot Stand	42	OZ 7520675	Neuro Adapter, Low Profile	31
OZ 2200120	Vacuum Repair Kit	42	OZ 7521275	Extension Plate, Imaging, 44" Long Carbon Fiber	20, 39
OZ 2200130	Non-Slip Sheet Roll Tabletop	43	OZ 7521375	Split Legs, Pair	17, 20
OZ 2200140	Vacuum Mattress, Universal	43	OZ 7521875	Footboard for Split Leg Plates, Pair	17, 20
OZ 2200150	Vacuum Mattress, Torso Section with Cutout	43	OZ 7522475	Transfer Board, Leg	17, 28
OZ 2200160	Vacuum Mattress, Universal, Split Leg	44	OZ 7522575	Rail Clamp, Imaging Extension	5
OZ 7500475	Armboard, Multitask	12, 35	OZ 7522675	Filler Plate	10, 18, 38
OZ 7500775	Quick Clamp, Round Post	4	OZ 7560175	Transfer, Leg, Sling Style	17, 28
OZ 7500875	Rail Clamp, Flat Blade	4	OZ 7523075	Shoulder Chair, D7XX	36
OZ 7500975	Restraint Strap, Body, Padded	8	OZ 7523275	Footboard/Foot Extension	11, 14
OZ 7501075	Lateral Support, Rectangular	44	OZ 7527675	Rail Clamp, Inside	5
OZ 7501275	Width Extender, Tabletop	22	OZ 7529075	Transfer Board, Leg, Adjustable	18, 29
OZ 7501575	Knee Crutches	16	OZ 7529375	Spine Frame	31

#### OPERON

#### Index

#### 0

OZ 7529475	Shoulder Chair Cart	7
OZ 7530475	Supine Tops	34
OZ 7530675	Leg Pad, Prone for Spine Frame System	34
OZ 7531875	Shoulder Support TrenStop	14, 45
OZ 7532175	Accessory Cart, Split Legs	14, 45
OZ 7532975	Orthopedic Extension Unit - Carbon Fiber	37
OZ 7533175		
	Prone Arm Support	14, 31
OZ 7537375	Lateral Positioning Set	15, 45
OZ 7538575	Head Donut, Gel	45
OZ 7538675	Head Donut, Gel	45
OZ 7538775	Head Donut, Gel, Horseshoe	46
OZ 7538875	Heel Pad, Gel	46
OZ 7538975	Heel Elevators	46
OZ 7539075	Armboard Pad, Gel, Round End	47
OZ 7539175	Ulnar Nerve Protector, Gel, Oval	47
OZ 7539275	Positioning Roll, Gel, Flat Bottom	47
OZ 7539375	Positioning Roll, Gel, Flat Bottom	47
OZ 7539475	Positioning Roll, Gel, Flat Bottom	47
OZ 7539575	Positioning Roll, Gel, Flat Bottom	47
OZ 7539675	Positioners, Gel, Standard	48
OZ 7539775	Positioners, Gel, Standard	48
OZ 7539875	Armboard Pads, Gel, Contoured	48
OZ 7540075	Knee Board, 18"	15, 32
OZ 7540175	Armboard Pads, Gel, Contoured	48
OZ 7540275	Leg Holder, J Shape	19, 29
OZ 7546375	Head Donut, Gel, Horseshoe, Pediatric	48
OZ 7546475	Headrest, Gel, Prone	49
OZ 7546575	Knee Crutch Pad, Gel	49
OZ 7546675	Arm Protector, Gel	49
OZ 7546775	Operating Table Pad, Gel	50
OZ 7546875	Operating Table Pad, Gel	50
OZ 7546975	Perineal Post Cushion, Gel	50
OZ 7547075	Knee Arthroscopy Insert, Gel	51
OZ 7547175	Spine Frame (Wilson Style)	32
OZ 7547275	Disposable Pad, Arm Cradle	32
OZ 7547375	Disposable Pad, Headrest	33
OZ 7547475	Pad Set, Spine Frame	33

OZ 7547575	Disposable Mask	3
OZ 7565475	Patient Transfer Boards	
OZ 7565575	Patient Transfer Boards	
OZ 7565675	Patient Transfer Boards	

#### Index

#### A

Accessory Cart, Split Legs	OZ 7532175	6
Accessory Stand	OZ 1205075	6
Anesthesia Screen	OZ 1001101	6
Anesthesia Tube Holder, Flexible	OZ 1001200	7
Armboard, Ball Joint, One Hand	OZ 1251875	12
Armboard, Bariatric Armboard,	OZ 7519175	13, 23
Multitask Armboard Pad, Gel,	OZ 7500475	12,35
Round End Armboard Pads, Gel,	OZ 7539075	47
Contoured	OZ 7539875	48
	OZ 7540175	
Arm Protector, Gel	OZ 7546675	49

Chin Straps Single-Use for Shoulder Chair	OZ 7518375	36
Cushion, Anesthesia	OZ 2200010	41
Cushion, Half Moon	OZ 2200020	41
Cushion, Lateral Positioning	OZ 2200050	42
Cushion, Leg Rest	OZ 2200030	41
Cushion, Spinal	OZ 2200060	42
Cushion, Tunnel	OZ 2200040	41

#### D

Disposable Mask	OZ 7547575	33
Disposable Pad, Arm Cradle	OZ 7547275	32
Disposable Pad, Headrest	OZ 7547375	33
Drainage Bag, Attachment Strap	OZ 7501875	26
Drainage Bags	OZ 7502075	27

#### E

Elbow Support, Piece	OZ 1252675	26
Extension for Upper Arm (Weinberger)	OZ 1001450	34
Extension Plate, Imaging, 44" Long Carbon Fiber	OZ 7521275	20, 39
Extension Plate, Imaging, Short, Carbon Fiber	OZ 7520175	19, 38

Filler Plate	OZ 7519975 OZ 7522675	10, 18 38
Footboard/Foot Extension	OZ 7523275	11, 14
Footboard for Split Leg Plates, Pair	OZ 7521875	17, 20
Foot Control	OZ 7518675	10

#### Н

Hand/Arm Table, Hourglass Shape	OZ 7516875	13
Head Donut, Gel	OZ 7538575	45
	OZ 7538675	
Head Donut, Gel, Horseshoe	OZ 7538775	46
Head Donut, Gel, Horseshoe, Pediatric	OZ 7546375	48
Headrest, Gel, Prone	OZ 7546475	49
Headrest, Horseshoe	OZ 1001351	24
Headrest, Horseshoe, ENT, Adjustable	OZ 1001400	25
Headrest, Multi, Adjustable	OZ 1001410	25
Headrest, Narrow	OZ 7517375	26
Heel Elevators	OZ 7538975	46
Heel Pad, Gel	OZ 7538875	46

#### **OPERON**

#### Index

Instrument Surface	OZ 7504075	27
Integrated Shoulder Chair, D7XX	OZ 7523075	36
Integrated Shoulder Chair, D8XX	OZ 7519475	36
K		
Kidney Support	OZ 7520075	11, 28
Knee Arthroscopy Insert, Gel	OZ 7547075	51
Knee Board, 18"	OZ 7540075	15, 32
Knee Crutches	OZ 7501575	16
Knee Crutch Pad, Gel	OZ 7546575	49
Knee-Elbow Positioning Device	OZ 1005275	30, 40
Lateral Positioning Set	OZ 7537375	
	OZ 7501075	,
Lateral Support, Rectangular	OZ 7501075 OZ 1002540 OZ 1002560	44
Lateral Support, Rectangular Leg Holder, Arthroscopy	OZ 7501075 OZ 1002540 OZ 1002560 OZ 7511975	44 16, 35
Lateral Support, Rectangular Leg Holder, Arthroscopy Leg Holder, J Shape	OZ 7501075 OZ 1002540 OZ 1002560 OZ 7511975 OZ 7540275	44 16, 35 19, 29
Lateral Support, Rectangular  Leg Holder, Arthroscopy  Leg Holder, J Shape  Leg Pad, Prone for Spine Frame System	OZ 7501075 OZ 1002540 OZ 1002560 OZ 7511975 OZ 7540275 OZ 7530675	44 16, 35 19, 29 34
Lateral Support, Rectangular Leg Holder, Arthroscopy Leg Holder, J Shape	OZ 7501075 OZ 1002540 OZ 1002560 OZ 7511975 OZ 7540275	44 16, 35 19, 29
Lateral Support, Rectangular  Leg Holder, Arthroscopy  Leg Holder, J Shape  Leg Pad, Prone for Spine Frame System	OZ 7501075 OZ 1002540 OZ 1002560 OZ 7511975 OZ 7540275 OZ 7530675	
Lateral Support, Rectangular  Leg Holder, Arthroscopy  Leg Holder, J Shape  Leg Pad, Prone for Spine Frame System  Leg Stirrups, Bariatric, Pair	OZ 7501075 OZ 1002540 OZ 1002560 OZ 7511975 OZ 7540275 OZ 7530675	44 16, 35 19, 29 34
Lateral Support, Rectangular  Leg Holder, Arthroscopy  Leg Holder, J Shape  Leg Pad, Prone for Spine Frame System  Leg Stirrups, Bariatric, Pair	OZ 7501075 OZ 1002540 OZ 1002560 OZ 7511975 OZ 7540275 OZ 7530675 OZ 7518975	44 16, 35 19, 29 34 22
Lateral Support, Rectangular  Leg Holder, Arthroscopy  Leg Holder, J Shape  Leg Pad, Prone for Spine Frame System	OZ 7501075 OZ 1002540 OZ 1002560 OZ 7511975 OZ 7540275 OZ 7530675 OZ 7518975	44 16, 35 19, 29 34 22

0		
Operating Table Pad, Gel	OZ 7546775 OZ 7546875	50
Orthopedic Extension Unit - Carbon Fiber	OZ 7532975	37
P		
Pad Set, Spine Frame	OZ 7547475	33
Patient Transfer Boards	OZ 7565475 OZ 7565575 OZ 7565675	9
Pediatric Table, Carbon Fiber	OZ 7520375	21
Perineal Post Cushion, Gel	OZ 7546975	50
Positioners, Gel, Standard	OZ 7539675 OZ 7539775	48
Positioning Roll, Gel, Flat Bottom	OZ 7539275 OZ 7539375 OZ 7539475 OZ 7539575	47
Prone Arm Support	OZ 7533175	14, 31
Q		
Quick Clamp, Flat Blade	OZ 7516175	5
Quick Clamp, Round Post	OZ 7500775	4
R		
Radial Setting Clamp	OZ 1001041	4
Radial Setting Clamp, Simultaneous Locking	OZ 1001075	4
Rail Clamp, Flat Blade	OZ 7500875	4
Rail Clamp, Imaging Extension	OZ 7522575	5
Rail Clamp, Inside	OZ 7527675	5

OZ 7503975

OZ 7500975

64

Rail Extender, Side

Restraint Strap, Body, Padded

#### **OPERON**

#### Index

Shoulder Chair Cart	OZ 7529475	5
Shoulder Support, Pair	OZ 7517075	44
Shoulder Support TrenStop	OZ 7531875	14, 45
Spine Frame	OZ 7529375	31
Spine Frame (Wilson Style)	OZ 7547175	32
Split Legs, Pair	OZ 7521375	17, 20
Stirrups, Premium, Pair	OZ 7501775	16
Strap, Leg & Body, Bariatric	OZ 7519075	22
Supine Tops	OZ 7530475	34
Support Roll, Horizontal	OZ 1001691	40
Tibia Support/Knee Rest Transfer Board, Leg	OZ 1003310 OZ 7522475	40 17, 28
Transfer Board, Leg, Adjustable	OZ 7529075	18, 29
Transfer, Leg, Sling Style	OZ 7560175	17, 28
Trigger Armboard	OZ 7504175	12
Trigger Armboard, Carbon Fiber	OZ 7547675	15
TUR Module	OZ 7519675	29
U		
Ulnar Nerve Protector, Gel, Oval	OZ 7539175	4
Universal Arm/Hand Table	OZ 7515775	13



W		
Vacuum Repair Kit	OZ 2200120	42
Vacuum Pump with Foot Stand	OZ 2200110	42
Vacuum Mattress, Universal, Split Leg	OZ 2200160	44
Vacuum Mattress, Universal	OZ 2200140	43
Vacuum Mattress, Torso Section with Cutout	OZ 2200150	43



Width Extenders set, Bariatric, Short Leg with 2 3/8" pad	OZ 7556975	23
Width Extenders set, Bariatric, Long Leg with 2 3/8" pad	OZ 7557175	23
Width Extenders set, Bariatric, Short Leg without 2 3/8" pad	OZ 7557075	24
Width Extenders set, Bariatric, Long Leg with 2 3/8" pad	OZ 7557275	25
Width Extender, Tabletop	OZ 7501275	22
Wrist Support for Horseshoe Headrest	OZ 1001352	24



X-Ray Cassette Tray	OZ 7518475	8
X-Ray Grid Encasement	OZ 7537 <b>4</b> 75	14

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.	Stryker Communications 1410 Lakeside Parkway #100 Flower Mound, TX 75028 USA 866 841 5663   www.stryker.com Manual No. 1000902045 Rev A
Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, OPERON.  All other trademarks are trademarks of their respective owners or holders.	Copyright © 2016 Stryker