Product Specification Sheet









Specifications:

Vendor Parts Number: 2020

Function: Placed around the wrist or ankle it is used to exercise the upper & lower extremities.

Product Dimension: The total length of the strap is 17" long. The foam pad is 8" long by 2" wide.

Shipping: 9"x 6"x 1" envelope. 0.5 lbs. UPS or

FedEx



Specifications:

Vendor Part Number: 2021

Function: Placed around the upper arm and thigh it is used to exercise the upper and lower extremities.

Product Dimension: The total length of the strap is 26". The foam pad is 18" long by 2" wide.

Shipping: 6"x 6"x 6" box, 1 lb. UPS or FedEx.







Specifications:

Vendor Part Number: 2022

Function: This long belt is used for trunk exercises in standing, sitting and the recumbent position.

Product Dimension: The total length of the strap is 40". The foam pad is 32" long by 3" wide.

Shipping: 6"x 6" x 6" box, 1 lb. UPS or FedEx.

RehabPro, Inc. Phone: 602-317-7207. E-mail: sales@rehabpropulleys.com
Designed by Physical Therapists. Manufactured in the USA.