


Medaysis

Hematologic Cancer Markers




By 2020, nearly 412,000 people will be diagnosed with some form of leukemia and nearly 537,000 will be diagnosed with lymphoma worldwide. Leukemia can be categorized based on which type of white blood cell is involved (lymphocytes or myeloid cells) as well as whether the illness is chronic, developing slowly and persisting for a long time, or acute, developing very quickly. Lymphoma is a type of cancer that originates in the part of the body's immune system known as the lymphatic system, which includes the lymph nodes and lymphoid tissues.


GLOBOCAN & Memorial Sloan Kettering Cancer Center




FFPE human tonsil stained with anti-Bcl-2




FFPE human seminoma stained with anti-CD117




FFPE human angiosarcoma stained with CD235a [GYPA/280]



FFPE human cervical ca stained with anti-c-Myc [MYC275+909]



FFPE human Hodgkin's lymphoma stained with anti-Fascin [FSCN1/417]



FFPE human tonsil stained with anti-Ki67 [MIB-1]

Hematologic Cancer Markers

Name	Cat. No.	Application
ALK/p80 [MD73R]	RM0461	Specific marker for anaplastic lymphomas
α -1-Antichymotrypsin [AAT/1451]	MC0141	Useful in the identification of neoplastic myeloid cells within extramedullary tissues such as acute myeloid leukemia
α -1-Antitrypsin [AAT/1378]	MC0142	Useful in identification of benign and malignant hepatic tumors and yolk sac carcinomas
Annexin 1 [ANXA1/3566]	MC0388	Specific to hairy cell leukemia that can be used to differentiate this disease from other B-cell lymphomas
Annexin 1 [EH17a]	MC0481	
Bcl-10 [SPM520]	MC0612	Useful tool for the sub-classification of lymphomas
Bcl-2 [8C8]	MC0614	Used to identify lymphoma, Burkitt's lymphoma, Follicular lymphoma, Acute Myeloid Leukemia Adult
Bcl-6 [PG-B6P]	MC0619	For classification of diffuse large B-cell lymphomas, follicular lymphomas and Burkitt's lymphoma/leukemia
Bcl-x [2H12]	MC0620	A high level found in Reed Sternberg cells in Hodgkin's disease. Overexpression observed in primary central nervous system lymphomas that occur in immunosuppressed patients
BOB.1 [BOB1/2424]	MC0239	Used as an aid in the identification of B-cells in germinal centers, mantle cells, and differentiation of lymphomas
BOB.1 [EPR17685]	RM0010	
CD10 [56C6]	MC0277	Used to identify Burkitt's lymphoma/leukemia, follicular lymphoma (except grade III), precursor B-cell lymphoblastic leukemia/ lymphoma and renal cell carcinoma of the clear cell type
CD10 [MME/1870]	MC0272	
CD10 [SP67]	RM0224	
CD103 [EPR22590-27]	RM0019	Present on all enteropathy-type T-cell lymphomas. Additionally, CD103 has been a useful marker for hairy cell leukemia
CD103 [ITGAE/2063]	MC0260	
CD117/c-Kit [C117/370]	MC0233	A marker for Mast cell and gastrointestinal stroma tumor
CD117/c-Kit [EPR22566]	RM0067	
CD11c [ITGAX/1284]	MC0150	May aid in identification of lesions with histiocytic origin. It may also been used as a marker for hairy cell leukemia
CD123 [IL3RA/2947R]	RM0220	Useful in diagnosis of Kikuchi disease, hyaline vascular type of Castleman disease, lupus, and primary cutaneous marginal zone B-cell lymphoma
CD13 [2D8D11]	MC0379	Used as a myeloid marker and identify AML subtype M0 acute lymphoid leukemia (ALL); a sensitive but not entirely specific marker for ALK+ anaplastic large cell lymphomas (ALCLs)
CD138 [B-A38]	MC0510	Used to identify differentiated plasma cells, multiple myeloma, and subclassification of diffuse large B-cell
CD138 [EPR6454]	RM0023	
CD138 [MD76R]	RM0463	

Hematologic Cancer Markers

Name	Cat. No.	Application
CD14 [MD85R]	RM0024	Used as a differentiation marker for monocytes/macrophages, labels Langerhans' cells and dendritic cells
CD15 [FUT4&815]	MC0646	Used to identify Chronic Myelogenous Leukemia, mononuclear Hodgkin's cells and Reed-Sternberg cells in Hodgkin's lymphoma
CD16 [DJ130C]	MC0331	Useful to differentiate hepatosplenic gamma delta T-cell lymphoma and gamma delta T-cell large granular lymphocyte leukemia from other peripheral T-cell lymphomas, such as mucosal and cutaneous gamma delta T-cell lymphoma
CD163 [EPR14336]	RM0027	Might be used for identifying tumors of monocytic origin
CD163 [M130/1210]	MC0367	
CD19 [MD77]	MC0394	A pan B cell marker expressed in a wide range of maturational stages including pre-B cells, found in the majority of B cell-derived malignancies
CD1a [O10]	MC0278	Used for identification of Langerhans Cell Histiocytosis, thymoma, precursor T-lymphoblastic leukemia/lymphoma
CD2 [MD150R]	RM0283	A pan T-cell marker, Useful for identification of precursor and mature T-cell lymphomas. Aberrant loss of CD2 in T-cell lymphomas may help to distinguish them from reactive T-cell proliferations
CD20 (B Cell) [L26]	MC0056	Used to identify B lymphocyte
CD21 [CR2/3124R]	RM0311	Used to identify follicular lymphoma, follicular dendritic cells, mature B cells and angioimmunoblastic T-cell lymphoma
CD21 [EP64]	RM0032	
CD22 [BLCAM/1795]	MC0252	Expressed in B lymphocytes and hairy cell leukemia
CD23 [HD50]	MC0158	Used to identify both normal B cells and malignant lymphomas such as B-cell chronic lymphocytic leukemia/small lymphocytic lymphoma
CD23 [MD78R]	RM0033	
CD235a [EPR8200]	RM0101	strongly labels normal erythroid cells at all stages of differentiation from the erythroblast to the mature red cell and does not react with glycophorin B
CD235a [GYPA/280]	MC0651	
CD25 [EP218]	RM0034	Used to distinguish neoplastic mast cell aggregates from reactive proliferations. It is a minor criterion for the diagnosis of bone marrow involvement in systemic mastocytosis (SM)
CD25 [IL2RA/2395]	MC0395	
CD3 [PC3/188A]	MC0273	Used to identify T- and NK-cell lymphomas and mycosis fungoides
CD30 [Ber-H2]	MC0314	Used to identify Reed-Sternberg cells in Hodgkin's lymphoma and anaplastic large cell lymphoma
CD30 [EP154]	RM0036	
CD31 (PECAM-1) [JC70A]	MC0514	Used for identification of vascular and lymphatic vessels and tumors of vascular origin

Hematologic Cancer Markers

Name	Cat. No.	Application
CD33 [HIM3-4]	MC0205	Useful for acute myeloid leukemia, minimally differentiated (AML-M0) and acute monocytic leukemia (AML-M5), myeloid sarcoma
CD33 [PWS44]	MC0515	
CD34 [QBEnd/10]	MC0066	Used to identify vascular and lymphatic vessels and tumors of vascular origin
CD35 [E11]	MC0663	Abnormal expression associated with follicular dendritic cell tumor, sarcoidosis, gallbladder carcinomas, mesangiocapillary glomerulonephritis, and systemic lupus erythematosus
CD35 [EP197]	RM0039	
CD38 [AT1]	MC0592	Helpful in the identification of plasma cells and tumors with plasmablastic differentiation
CD38 [MD71R]	RM0459	
CD4 [MD143]	MC0527	Used to identify thymocytes, T lymphocyte, anaplastic large cell lymphoma and mycosis fungoides
CD4 [MD142R]	RM0042	
CD4 [MT310]	MC0237	
CD43 [DF-T1]	MC0516	Used for identification of low-grade B-cell lymphoma, peripheral T-cell lymphoma and myeloid disorders
CD43 [MD133R]	RM0434	
CD44 [156-3C11]	MC0666	Expressed on approximately 90% of lymphocytes, monocytes, granulocytes, and, in lower amounts on thymocytes, fibroblasts, and erythrocytes
CD44 [MD132R]	RM0044	
CD45 [2B11&PD7/26]	MC0316	Used to identify tumors of lymphoid origin
CD45RA [111-1C5]	MC0985	Expressed on 40-50% of peripheral CD4+ T-cells, 50% of peripheral CD8+ T-cells, B-cells, and leukemic B-cell lines
CD45RA [4KB5]	MC0517	
CD45RB [BRA-11G]	MC0678	Used to identify tumors of lymphoid origin
CD45RO [UHL-1]	MC0078	Useful marker for reactive T cells and T-cell neoplasms
CD5 [C5/473&54/F6]	MC0681	Used to identify B-cell chronic lymphocytic leukemia/ small lymphocytic lymphoma, mantle cell lymphoma, and T-cell lymphoma and leukemia
CD5 [EP77]	RM0047	
CD56/NCAM [123C3]	MC0290	Used to identify natural killer (NK) cells, NK-like T cells, small cell carcinoma of the lung and carcinoid tumor
CD57/B3GAT1 [NK-1]	MC0519	Used to identify lymphocyte-predominant HL, T-cell large granular lymphocyte disorder, oligodendroglioma and neuroendocrine tumor
CD61 [2f2]	MC0520	Used for diagnosis of megakaryoblastic leukemia and absence of GpIIb/IIIa in Glanzmann's thrombastenia
CD61 [Y2/51]	MC0137	
CD68 [KP1]	MC0084	Used to identify macrophages, M4 and M5 types of acute myeloid leukemia (+), histiocytic sarcoma (+), M1 acute myeloid leukemia (-)
CD7 [EP132]	RM0050	Used for identification of T lymphocyte, Acute Lymphoblastic Leukemia
CD7 [MD146]	MC0274	
CD71 [EPR20584]	RM0284	Useful in identifying erythroid precursors
CD71 [TFRC/1817]	MC0257	

Hematologic Cancer Markers

Name	Cat. No.	Application
CD74 [LN-2]	MC0701	A marker for distinguishing atypical fibroxanthoma from malignant fibrous histiocytoma and distinguishing leiomyosarcoma from leiomyoma
CD74 [MD83R]	RM0051	
CD79a [EP82]	RM0052	Used to identify B lymphocyte and B cells during all stages of differentiation
CD79a [JCB117]	MC0319	
CD8 [C8/468+144B]	MC0706	Used to identify cutaneous T-Cell lymphoma, mycosis fungoides
CD8 [CD8/4391R]	RM0409	
CK20 [KRT20/1993]	MC0174	A well-established marker for colon epithelium, may be used for detecting circulating tumor cells (CTCs) in patients with metastatic colorectal cancer
CK20 [Ks20.8]	MC0114	
CK7 [OV-TL12/30]	MC0754	A marker for lung, breast and urothelial cancer, may be used for detecting circulating tumor cells (CTCs) in patients with metastasis
c-Myc [EP121]	RM0070	A marker for Burkitt lymphoma, and diffused large B-cell lymphoma
c-Myc [MYC275+909]	MC0134	
Cyclin D1 [DCS-6]	MC0732	Used to identify mantle cell lymphoma (+), and B-cell chronic lymphocytic leukemia/small lymphocytic lymphoma (-)
Cyclin D1 [EP12]	RM0071	
Factor VIII Polyclonal	RC0325	Reacts specifically with the endothelial cells of normal, reactive, and neoplastic blood and lymphatic vessels, endocardium, platelets, and megakaryocytes
Factor VIII RA [VWF/1859R]	RM0412	a sensitive marker of benign blood vessels, and aid for classification of tumors derived from megakaryocytes
Fascin [FSCN/417]	MC0787	Highly sensitive for staining Reed-Sternberg cells making it an excellent marker for classic Hodgkin's lymphoma
FOXP1 [EPR4113]	RM0096	Overexpression in a subset of diffuse large B-cell lymphoma (DLBCL), helpful in the subclassification of DLBCL
Galectin-3 [B2C10]	MC0132	Used to differentiate between benign and malignant thyroid neoplasms in both histologic sections and fine needle aspiration biopsy material, and identifying anaplastic large cell lymphoma
Glycophorin A [EPR8200]	RM0101	Useful in identifying cells of the erythroid lineage
Glycophorin A [GYPA/280]	MC0651	
Granzyme B [GZMB/3014]	MC0258	Useful in the identification of natural killer/T-cell lymphoma, and anaplastic large cell lymphoma
Hairy Cell Leukemia [DBA.44]	MC0135	Combining with TRAcP, it is a good tool to detect HCL

Hematologic Cancer Markers

Name	Cat. No.	Application
hGAL (GCET2) [EPR14333]	RM0373	A useful marker for diffuse large B-cell lymphoma (DLBCL) and several germinal center (GC)-like lymphoma cell lines, normal GC lymphocytes and GC-derived malignancies
IgA [IA761]	MC0811	Used to identify plasma cells and related lymphoid cells containing IgA, and for the classification of B-cell neoplasias
IgD [EPR6146]	RM0111	Used to identify mantle cell lymphoma, splenic marginal zone lymphoma, B-cell chronic lymphocytic/leukemia/small lymphocytic lymphoma and rare subsets of multiple myeloma
IgG [B33/20]	MC0812	Used to identify plasma cells, lymphoid cells and IgG plasma cell neoplasia containing IgG
IgM [DA4-4]	MC0815	Used to identify plasma cells and lymphoid cells containing IgM, classification of IgM subtype for B-cell neoplasia
Kappa [HP6053]	MC0338	Used to identify plasma cells and lymphoid cells containing kappa light chains, and classification of kappa light chain restriction for monoclonal gammopathies and amyloidosis
Kappa [MD137R]	RM0117	
Ki67 [MDKI67]	MC0196	
Ki67 [MIB-1]	MC0185	Used for identification of Burkitt's lymphoma/leukemia
Ki67 [SP6]	RM0255	
Lambda [HP6054]	MC0341	Used to identify plasma cells and lymphoid cells containing lambda light chains, and for classification of lambda light chain restriction for monoclonal gammopathies and amyloidosis
Lambda [MD138R]	MC0341	
Langerin/CD207 [12D6]	MC0544	A panel of langerin, CD1a, CD21, CD23, CD35, S100 may distinct Langerhans cell histiocytosis, histiocytic sarcoma, interdigitating dendritic cell sarcoma, follicular dendritic cell sarcoma, disseminated juvenile xanthogranuloma, and Rosai-Dorfman disease
Langerin/CD207 [H4]	MC0325	
LEF1 [B6]	MC0374	May be an aid for differentiation of chronic lymphocytic leukemia/small lymphocytic lymphoma from other small B cell lymphomas
LEF1 [EP310]	RM0381	
LMO2 [MD112]	MC0412	A tool in the identification of lymphomas of B-cell origin
LMO2 [SP51]	RM0256	
Lysozyme [EPR2994(2)]	RM0121	Used to identify histiocytic neoplasias and myeloid leukemias
Macrophage L1 Pro [MAC387]	MC0073	A marker for neutrophils, monocytes, macrophages, and squamous mucosal epithelia to identify macrophages in tissue sections
Macrophage Spec [LN-5]	MC0834	An important tool for the study of histiocytic lesions
Mesothelial C [HBME-1]	MC0343	Used to identify epithelial mesotheliomas
Mesothelin [MSLN/2131]	MC0207	

Hematologic Cancer Markers

Name	Cat. No.	Application
MUM1 [EP190]	RM0135	Used to identify epithelial mesotheliomas
MPO [EPR17996]	RM0137	Useful in identifying myeloid leukemias and granulocytic sarcoma
MPO Polyclonal	RC0291	
OCT2 [EPR542]	RM0147	A marker for various B cell lymphomas with strong expression in germinal center-derived lymphomas, lymphocyte predominance HLs
OCT2 [Oct2/2136]	MC0297	
PAX5 [A11]	MC0203	Used to identify epithelial mesotheliomas
PAX5 [EPR3730(2)]	RM0156	
PAX5 [PAX5/3735]	MC0387	
PD-1 [NAT105]	MC0547	A useful marker for angioimmunoblastic T-cell lymphoma
PD-1 [PDCD1/922]	MC0909	
Perforin [E-5]	MC0980	A specific marker of functionally active cytotoxic T-lymphocytes (CTLs) and natural killer (NK) cells
Perforin [PRF1/2470]	MC0301	
SOX11 [MD70R]	RM0411	Expressed in most mantle cell lymphomas (MCL) and, in particular, in all cyclin D1 negative MCLs
SOX11 [MRQ-58]	MC0572	
SOX11 Polyclonal	RC0300	
Spectrin [SPTA1/2939R]	RM0018	A marker to glycophorin for the identification of erythroid precursors
T-bet/TBX21 [MD40R]	RM0343	A useful new marker for Hodgkin's lymphoma and identification of hairy cell leukemia (HCL).
TCL1 [MD117R]	RM0182	Useful in identifying B-cell lymphomas including follicular lymphoma and Burkitt lymphoma
TCL1 [TCL1/2079]	MC0302	
TdT [MD154R]	RM0305	Used to identify precursor acute lymphoblastic leukemia and thymoma
TdT [TDT/1393]	MC0276	
Top2A [Ki-S1]	MC0575	Used to identify proliferating cells in Hodgkin's Lymphoma
Top2A [TOP2A/1362]	MC0246	
TRAcP [9C5]	MC0385	Utility in providing adjunctive information for the identification of hairy cell leukemia
TRAcP [ACP5/1070]	MC0958	
TF PU.1 [EP18]	RM0190	A marker for B lymphocytes in germinal center and mantle B cells, labels many types of B-cell lymphomas including mantle cell lymphoma, but is not expressed in classical Hodgkin lymphoma (cHL)
TF PU.1 [PU1/2146]	MC0421	
ZAP70 [2F3.2]	MC0372	Used to identify ZAP-70 in a subset of B-cell chronic lymphocytic leukemias/small lymphocytic lymphoma
ZAP70 [EP52]	RM0197	

Research Use Only