


# *Schwartz Custom AFO Orthometry Form - Model SAFO-C*



Date \_\_\_\_\_ Date Needed \_\_\_\_\_  
 Practitioner \_\_\_\_\_  
 P.O. # \_\_\_\_\_  
 Company Name \_\_\_\_\_  
 Ship to: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Patient Name: \_\_\_\_\_  
 Gender: Male / Female / Right / Left  
 Height \_\_\_\_\_ Weight \_\_\_\_\_  
 Diagnosis \_\_\_\_\_

Shoe Type \_\_\_\_\_  
 Rocker Bottom: Yes / No  
 AFO Upper Plastic: Co-poly / Poly-Pro / Carbon  
 Plastic Color: \_\_\_\_\_  
 Valgus/Varus attachment: Yes / No

### INSTRUCTIONS

1. Send cast or take measurements of the ankle and tibia section of the patient
2. Send shoe with cast. Make sure the shoe used has a thick, solid sole. We will be splitting the outer sole and installing a carbon fiber shank the entire length of the sole. A rocker bottom will also be added.
3. If you do not want the rocker bottom, please indicate that on the order form. If a rocker bottom is needed, please send both shoes.
4. Indicate the color of the upper portion of the AFO. Black is standard.
5. Please indicate if you would like the optional Valgus/Varus attachment on the order form.
6. Standard turn-around time is 5 days but rush service is available for an additional charge.



### SPECIAL INSTRUCTIONS

### Components

- |          |  |
|----------|--|
| SAFO-001 | Carbon Foot Shank w/hardware                 |
| SAFO-002 | AFO Custom Upper Section (Copoly or PolyPro) |
| SAFO-003 | Universal "T" Attachment w/strap             |
| SAFO-004 | AFO pad                                      |
| SAFO-005 | AFO strap                                    |
| SAFO-006 | Rocker bottom                                |
| SAFO-007 | Complete Prefab SAFO kit                     |
| SAFO-008 | Custom SAFO (Includes Installation)          |
| SAFO-009 | Install rocker bottom on opp. Shoe           |
| SAFO-010 | Custom Carbon AFO Upper                      |
| SAFO-011 | AFO Prefab Upper Section (Copoly or PolyPro) |

\* kit includes shank, rocker pad, hardware and AFO Prefab upper



1971 Old Covington Road  
 Conyers, GA 30013  
 800-922-5155 / 800-813-8139 Fax