

# Liquid Nitrogen Freezer Report



## General Freezer Information

Company:

Brand:

Report Date:

Serial Number:

Model Number:

Controller Type:

## Liquid Nitrogen Level Tests and Calibration

Level Calibration:

Liquid Level Alarm Test:

Controller Level (inches):

Low Level Alarm (inches):

Low Level Fill Point (inches):

Manual Level (inches):

Hi Level Alarm (inches):

Hi Level Fill Point (inches):

## Temperature Probe Tests and Calibration

Temp A Calibration:

Temp A Reading: 0

Temp A Level:

Temp B Calibration:

Temp B Reading: 0

Temp B Level:

High Temp Alarm Test:

## General Freezer Function Tests

Fill Test:

Manual Fill Amount (inches):

Manual Fill Time (minutes):

Hot Gas Bypass Function:

Battery Backup Function:

Data Download:

## General Freezer Inspections and Parts Replaced

Lid Hinge Inspection:

Liquid Level Tube:

In-Line Filter Change:

Folding Step Inspection:

Thawed Cork:

Solenoid Change:

Wiring:

Lid Gasket Change:

Relief Valve Change:

Notes

Performed By: \_\_\_\_\_

Signature: \_\_\_\_\_

Paramount Gas Products has performed the liquid nitrogen freezer maintenance checks requested by the customer in this report. This report states that the settings and calibrations requested by the customer were set at the time the maintenance was performed. This report in no way guarantees the future performance of the freezer. Paramount Gas Products is not liable for any losses due to freezer performance or malfunction.