

## Kitty Medic Box



Anyone who owns a cat knows there is a real need for a way to medicate the cat or other small animal at home and save a trip to the vet to have the pet medicated. Pet owners know how difficult it is to get cats or other small animals to stay still long enough to give vital shots or give medicine orally. Stubborn cats often struggle to get away, sometimes scratching and injuring the owner. The Kitty Medic™ Box is the perfect solution to make the process easy and safe for both pet owner and animal. This is the perfect solution for clipping nails, checking temperature, or giving shots without taking the animal out of the box.



DanielInvents.com  
Pet-Wel.com

P.O. Box 423  
Grapeland, TX 75844  
Cell: 936-546-4500  
Phone: 936-687-2666  
Fax: 936-687-5416

Email:  
wildwoodrefuge@valornet.com



## Kitty Medic Box

DanielInvents.com/ Pet-Wel.com

The Kitty Medic™ Box is a compact, lightweight 14" x 11" x 8" canvas case that keeps your animal safe and secure while administering medications. It features a special panel on the side for giving shots, and oral medications can be safely given while the cat's head remains in place through the opening.



***Pet medication carrier***

Phone: 936-546-4500  
wildwoodrefuge@valornet.com

## Kitty Medic Box



### The perfect solution

Even though Kitty seems to have **nine lives**, anyone who owns a feline pet knows there is a real need for a way to provide routine care for kitty as well as other small animals at home and, in some cases, saving a traumatic trip to the vet. Pet owners know how difficult it is to get their feline pet or other small animals to stay still long enough to give vital care. A frightened feline will surely struggle to get away, sometimes scratching and injuring the owner. The Pet-Wel Kitty Medic Soft-Sided Carrier is the perfect solution to make the process easy, safe, and less stressful for both pet owner and animal.



As a concerned pet owner, you can now handle routine procedures such as clipping nails, checking temperature, giving injections, and TLC, all without taking the animal out of the carrier, thanks to the patented design of the Pet-Wel Kitty Medic. To use, simply place the pet inside the Pet-Wel Kitty Medic Carrier. Make sure that the animal's head is through the opening. Then adjust the built-in collar device to comfortably constrain kitty. Because kitty can't move around much, kitty stays calm and comfortable inside. FYI: Experience using the carrier



Use the Kitty Medic to transport your pet.

has provided an unexpected benefit: Some cats actually like the Kitty Medic Carrier. Panels are designed to give access to all areas of the animal's anatomy. Now you can easily administer injections quickly and easily without the animal seeing the needle by simply lifting the side panel. If you need to give oral medications, approach the kitty from the front and open its mouth to give drops or other medications. The carrier also works great for other small pets such as Chihuahua dogs and can be special ordered in a variety of sizes for other animals. For larger dogs, wheels can be installed for ease of handling. The Pet-Wel Kitty Medic is the perfect carrier to take your pet to the veterinarian's office or administer care at home.

The Kitty Medic is hand washable.

*Watch our video at [pet-wel.com](http://pet-wel.com).  
Also visit [danielinvents.com](http://danielinvents.com) for  
more information about our  
company.*

*Or call Daniel at 936-546-4500*