

Advanced 3D Material Makes the Difference in Burn Management

TheraBond 3D looks different, feels different and is different from previous generation dressings used in burn management. While offering rapid, sustained antimicrobial activity, TheraBond 3D provides unmatched fluid and exudate management capabilities, conformability and comfort.

Active Fluid and Exudate Management

Unlike previous generation dressings, TheraBond 3D is a single piece dressing designed with struts within the fabric that act as capillaries, moving excess fluid and exudate away from the wound to an inexpensive, outer absorbent dressing. This minimizes the number of dressing changes, reduces patient discomfort and frees up valuable clinician time.

Exceptional Comfort and Conformability


The soft TheraBond 3D material provides gentle compression, remains dimensionally stable and sets new standards of comfort and conformability, even over joints.

Easy Application

The soft, conformable TheraBond 3D can often be applied by a single clinician with minimal discomfort to the patient.

No Stinging or Staining

TheraBond 3D does not sting or stain upon application.



TheraBond 3D's unmatched advantages can mean:

- ✓ Fewer dressing changes
- ✓ Less medication, pain and anxiety for patients
- ✓ Minimization of clinician time
- ✓ Cost savings for the burn center, while providing better patient care