



HEK293 Wnt TCF Reporter Cell Line-Mutant

Catalog Number:	WRHEK293M
Source:	Human embryonic kidney cell line
Synonyms:	Wnt reporter, TCF reporter, LEF reporter cell line
Background	The WNT gene family consists of structurally related genes that encode secreted signaling proteins, membrane bound receptors, and signaling transduction proteins. These proteins have been implicated in oncogenesis, adipogenesis, etc. and in several other developmental processes, including regulation of cell fate and patterning during embryogenesis. Activity of the Wnt signaling pathway leads to nuclear translocation of β -catenin and the formation of TCF transcription factor complex. The TCF complex interacts with Wnt gene transcriptional response elements and leads to the expression of Wnt-responsive genes.
Product Description	This Wnt reporter cell line is a control cell line, designed to show background luciferase activity for HEK293 Wnt TCF Reporter Cell Line-Active (Cat: WRHEK293A). This human embryonic kidney cell line hosts CMV promoter, tandem repeats of a mutant TCF transcriptional response element, luciferase gene, and GFP gene. GFP expressed constantly can serve as control of cell numbers.
Activity:	This Wnt reporter cell line expresses low luciferase. The luciferase activity does not increase in response to Wnt3a stimulation at 1000 ng/mL (Fig. 1). Endogenous GFP expression from this Wnt reporter cell line is shown in Figure 2.



Mycoplasma Negative. Test reagents: Mycoplasma Detection Kit-Digital Test v2.0 from Biotool

Handling and Storage

The cell line may be shipped in dry ice or RT in either 25 cm² flask or 15 mL tube. If the cell line is shipped in dry ice, after receiving, store cells at -80°C or in Liquid Nitrogen or culture under standard culture conditions.

Luc Assay

Using normal tissue culture-treated plate: Seed 0.5 mL of cells into each well of 24 wells plate at a density of 10×10^4 cells/mL in complete EMEM medium (Corning Ref: 10-010-CV plus 1 mM sodium pyruvate and 1500 mg/L sodium bicarbonate), incubate cell at 5% CO₂, 37°C incubator overnight, replace complete EMEM with 198 µL EMEM without serum, add 2 µL of control buffer or Wnt3a (concentration range: 0.06 to 1 µg/mL), mix well and return plate into 5% CO₂, 37°C incubator and continue to incubate for 6 to 8 hours, suction out medium, lyse cells with 0.2 mL of cell lysis buffer, incubate for 5 to 10 min on rocking shaker at room temperature, transfer 50 µL cell lysate from each well into the wells of a 96 well black plate, read fluorescent first, and then add 50 µL of Luciferase substrate into each well, read Luciferase activity within 30 min. Fluorescent reading can serve as control of cell numbers.

Reference

Molenaar M. XTcf-3 transcription factor mediates beta-catenin-induced axis formation in *Xenopus* embryos. *Cell*. 1996; 86:391-9

Xing-Yao LI. A reporter gene system for screening inhibitors of Wnt signaling pathway. *Nat. Prod. Bioprospect*. 2013; 3: 24–28